

NAT'L INST. OF STAND & TECH



A11106 851425

A UNITED STATES
DEPARTMENT OF
COMMERCE
PUBLICATION

PUBLICATIONS

NBS SPECIAL PUBLICATION 260

1973 Catalog



STANDARD
REFERENCE
MATERIALS

U.S.
DEPARTMENT
OF
COMMERCE

QC 100
U57

al
au
of
ds

no. 260
1973
C.2

NATIONAL BUREAU OF STANDARDS

The National Bureau of Standards¹ was established by an act of Congress March 3, 1901. The Bureau's overall goal is to strengthen and advance the Nation's science and technology and facilitate their effective application for public benefit. To this end, the Bureau conducts research and provides: (1) a basis for the Nation's physical measurement system, (2) scientific and technological services for industry and government, (3) a technical basis for equity in trade, and (4) technical services to promote public safety. The Bureau consists of the Institute for Basic Standards, the Institute for Materials Research, the Institute for Applied Technology, the Center for Computer Sciences and Technology, and the Office for Information Programs.

THE INSTITUTE FOR BASIC STANDARDS provides the central basis within the United States of a complete and consistent system of physical measurement; coordinates that system with measurement systems of other nations; and furnishes essential services leading to accurate and uniform physical measurements throughout the Nation's scientific community, industry, and commerce. The Institute consists of a Center for Radiation Research, an Office of Measurement Services and the following divisions:

Applied Mathematics — Electricity — Mechanics — Heat — Optical Physics —
Linac Radiation² — Nuclear Radiation² — Applied Radiation² — Quantum
Electronics³ — Electromagnetics³ — Time and Frequency³ — Laboratory
Astrophysics³ — Cryogenics³.

THE INSTITUTE FOR MATERIALS RESEARCH conducts materials research leading to improved methods of measurement, standards, and data on the properties of well-characterized materials needed by industry, commerce, educational institutions, and Government; provides advisory and research services to other Government agencies; and develops, produces, and distributes standard reference materials. The Institute consists of the Office of Standard Reference Materials and the following divisions:

Analytical Chemistry—Polymers—Metallurgy—Inorganic Materials—Reactor
Radiation—Physical Chemistry.

THE INSTITUTE FOR APPLIED TECHNOLOGY provides technical services to promote the use of available technology and to facilitate technological innovation in industry and Government; cooperates with public and private organizations leading to the development of technological standards (including mandatory safety standards), codes and methods of test; and provides technical advice and services to Government agencies upon request. The Institute also monitors NBS engineering standards activities and provides liaison between NBS and national and international engineering standards bodies. The Institute consists of a Center for Building Technology and the following divisions and offices:

Engineering and Product Standards—Weights and Measures—Invention and
Innovation—Product Evaluation Technology—Electronic Technology—Techni-
cal Analysis—Measurement Engineering—Building Standards and Code Serv-
ices⁴—Housing Technology⁴—Federal Building Technology⁴—Structures, Mate-
rials and Life Safety⁴—Building Environment⁴—Technical Evaluation and
Application⁴—Fire Technology.

THE INSTITUTE FOR COMPUTER SCIENCES AND TECHNOLOGY conducts research and provides technical services designed to aid Government agencies in improving cost effectiveness in the conduct of their programs through the selection, acquisition, and effective utilization of automatic data processing equipment; and serves as the principal focus within the executive branch for the development of Federal standards for automatic data processing equipment, techniques, and computer languages. The Center consists of the following offices and divisions:

Information Processing Standards—Computer Information—Computer Services
—Systems Development—Information Processing Technology.

THE OFFICE FOR INFORMATION PROGRAMS promotes optimum dissemination and accessibility of scientific information generated within NBS and other agencies of the Federal Government; promotes the development of the National Standard Reference Data System and a system of information analysis centers dealing with the broader aspects of the National Measurement System; provides appropriate services to ensure that the NBS staff has optimum accessibility to the scientific information of the world, and directs the public information activities of the Bureau. The Office consists of the following organizational units:

Office of Standard Reference Data—Office of Technical Information and
Publications—Library—Office of International Relations.

¹ Headquarters and Laboratories at Gaithersburg, Maryland, unless otherwise noted; mailing address Washington, D.C. 20234.

² Part of the Center for Radiation Research.

³ Located at Boulder, Colorado 80302.

⁴ Part of the Center for Building Technology.

NATIONAL BUREAU OF STANDARDS

AUG 8 1973

not acc

QC100

U57

no. 260

1973

e. 2

Standard Reference Materials 1973 Catalog

Office of Standard Reference Materials
Institute for Materials Research
U.S. National Bureau of Standards
Washington, D.C. 20234

t Special publication no. 260.

CAUTION: The values given in the following sections are listed primarily as a guide to purchaser. The values shown are nominal and may differ from those shown on the certificates. Space limitations have required that some values be omitted. For these reasons, the certificates issued with the standards should always be consulted to obtain the certified values.



U.S. DEPARTMENT OF COMMERCE, Frederick B. Dent, Secretary
NATIONAL BUREAU OF STANDARDS, Richard W. Roberts, Director

Issued April 1973

Library of Congress Catalog Card Number: 72-600291

National Bureau of Standards Special Publication 260

(Supersedes NBS Spec. Publ. 260-1970 Edition)

Nat. Bur. Stand. (U.S.), Spec. Publ. 260-1973 catalog, 96 pages (Apr. 1973)

CODEN: XNBSAV

CONTENTS

GENERAL INFORMATION	1	Nickel base	
Renewal and Discontinued Standard		Nickel oxide	
Reference Materials		Selenium base	
New Standard Reference Materials		Tin base	
Catalog Supplements		Titanium base	
Ordering		Zinc base	
Terms		Zirconium base	
Domestic Shipments		Nonferrous Alloys (solid form)	34
Foreign Shipments		Copper base	
GUIDE FOR THE SUBMISSION OF		Lead base	
REQUESTS FOR THE DEVELOP-		Nickel base	
MENT OF NEW OR RENEWAL		Tin base	
STANDARD REFERENCE MATERIALS	12	Titanium base	
OTHER SERVICES OF THE NATIONAL		Titanium base-oxygen and	
BUREAU OF STANDARDS	14	hydrogen only	
Calibration and Test Services		Zinc base	
National Standard Reference Data		Die casting alloys and spelter	
System		Zirconium base	
Standards Information Services		Gases in Metals	38
Technical Information and Publications		High-Purity Metals	40
STANDARDS OF CERTIFIED CHEMICAL		Microprobe Standards	41
COMPOSITION		Primary, Working, and Secondary	
Steels (chip form)	16	Standard Chemicals	42
Plain carbon steel		Microchemical Standards	42
Low alloy steel		Clinical Laboratory Standards	43
High alloy steel		Biological Standards	43
Stainless steel		Metallo-Organic Compounds	44
Tool steel		Environmental Standards	45
Steels (granular form)	18	Analyzed gases	
Steels (solid form)	20	Permeation tubes	
Ingot iron and low-alloy steel		Analyzed liquids	
Special ingot irons and		Analyzed solids	
low-alloy steel		Hydrocarbon Blends	46
Stainless steel		Ores	47
Tool steel		Cements	48
Maraging steel		Fertilizers	48
High-temperature alloys		Minerals, Refractories, Carbides,	
Oxygen standards		and Glasses	48
Cast Irons (chip form)	28	Trace Element Standards	51
Cast Steels, White Cast Irons, Ductile		Nuclear Materials	52
Irons, and Blast Furnace Irons		Special nuclear materials	
(solid form)	28	Plutonium assay	
Steelmaking Alloys	30	Plutonium isotopic	
Nonferrous Alloys (chip form)	31	Uranium assay	
Aluminum base		Uranium isotopic	
Cobalt base		Neutron density standards	
Copper base		Isotopic Reference Standards	53
Lead base		Ion Activity Standards	54
Magnesium base		pH standards	
		pD standards	
		Ion selective electrodes	

STANDARDS OF CERTIFIED PHYSICAL PROPERTIES

Mechanical and Metrology Standards . . .	55
Coating thickness	
Density standards	
Glass viscosity standards	
Polymer standards	
Elasticity standards	
Heat standards	60
Superconductive thermometric fixed point devices	
Freezing point standards	
Defining fixed points	
Determined reference points	
Melting point standards	
Calorimetric standards	
Combustion calorimetry	
Solution calorimetry	
Heat source calorimetry	
Differential thermal analysis	
Vapor pressure standards	
Thermal conductivity standards	
Thermal expansion standards	
Thermocouple materials	
Magnetic Standards	64
Magnetic susceptibility standard	
Optical Standards	65
Spectrophotometric standards	
Color standards for spectrophotometric-tristimulus integrator system	
Glass and liquid filters for spectrophotometry	
Thermal emittance standards	
Refractive index standards	
Reflectance standards	
Radioactivity Standards	68
Alpha-particle standards	
Beta-particle and gamma-ray gas standards	
Beta-particle, gamma-ray, and electron-capture solution standards	
Contemporary standard for carbon-14 dating laboratories	
Environmental standards	
Gamma-ray "point-source" standards	
Radium gamma-ray solution standards	
Radium solution standards for radon analysis	

Metallurgical Standards	71
Mossbauer Standards	71
Permittivity Standards	72

ENGINEERING TYPE STANDARDS

Standard Rubber and Rubber Compounding Materials	73
Reference Magnetic Tape	73
Sizing Standards	74
Glass spheres for particle size	
Turbidimetric and fineness standard (cement)	
Color Standards	74
The ISCC-NBS centroid color charts	
Paint pigment	
Light-sensitive papers	
Light-sensitive plastic chips	
Photographic Standards	75
Surface Flammability Standards	75
Smoke Density Chamber Standards	76
Water Vapor Permeance	76
Internal Tearing Resistance Paper	76
Tape Adhesion Testing Standard	76

RESEARCH MATERIALS	77
High-Purity Materials	
Phosphors	

GENERAL MATERIALS	78
Hydrogen in Steel	
Nickel and Vanadium in Residual Oil	
Absorbent Clay, (ASTM D2007)	

INDEXES

Alphabetical Index	79
Numerical Index to SRM Certificates . . .	85

Standard Reference Materials

Issued by the National Bureau of Standards

This Catalog lists and describes the Standard Reference Materials (SRM's), Research Materials (RM's), and General Materials (GM's) currently distributed by the National Bureau of Standards, as well as many of the materials currently in preparation. SRM's are used to calibrate measurement systems and to provide a central basis for uniformity and accuracy of measurement. The unit and quantity, the type, and the certified characterization are listed for each SRM, as well as directions for ordering. The RM's are not certified, but are issued to meet the needs of scientists engaged in materials research. RM's are issued with a "Report of Investigation," the sole authority of which is the author of the report. The GM's are standardized by some agency other than NBS. NBS acts only as a distribution point and does not participate in the standardization of these materials. Announcements of new and renewal SRM's, RM's, and GM's are made in the semi-annual supplements to this Catalog, SRM Price and Availability List, the NBS Technical News Bulletin, and in scientific and trade journals.

Key words: Analysis; characterization; composition; properties; Standard Reference Materials; Research Materials; General Materials.

General Information

All of the Standard Reference Materials (SRM's), Research Materials (RM's), and General Materials (GM's) listed in this Catalog bear distinguishing names and numbers by which they are permanently identified. Each SRM, RM, or GM bearing a given designation is of identical characterization with every other sample bearing the same designation, within the limits required by the use for which it is intended; or if necessary, it is given a serial number and an individual calibration.

The first SRM's issued by the Bureau were a group of ores, irons, and steels certified for chemical composition, and by custom they came to be called "standard samples." At present, more than 800 SRM's are available, covering a wide range of chemical and physical properties, and the designation, Standard Reference Material, is more appropriate. As the number of SRM's has increased, so has the variety, with such new groups being established as: clinical laboratory standards, nuclear materials, glass viscosity standards, rubber and rubber compounding materials, color standards, and coating thickness standards. These groups are listed under the headings: Standards of Certified Chemical Composition, Standards of Certified Physical Properties, Engineering Type Standards, Research Materials, or General Materials. The groups of materials under these general headings are listed in the Table of Contents. An alphabetical index provides the location of a particular material, or group of similar materials. A numerical index provides the date of the current Certificate issued with these materials.

The detailed listing of materials indicates the nominal certification for which the material is issued, but the Certificate must be consulted for the actual certification. A number of SRM's are issued for which it is not feasible to supply numerical values, or for which such certification would not be useful. These SRM's provide assurance of identity among all samples with the same designation, and permit standardization of test procedures and referral of physical or chemical data on unknown materials to a common basis.

Renewal and Discontinued Standard Reference Materials

The preparation of "renewal" SRM's is intended to be completed by the time the existing supply of each kind of material is exhausted, but this is not always possible. The renewal will not usually be identical to its predecessor, but will be quite similar especially with regard to the characteristics certified, and generally the renewal can be used in place of its predecessor. As an example, when the first 0.1 percent carbon Bessemer steel was prepared in 1909, it was called Standard Sample No. 8. During the following years, a number of renewal batches, 8a, 8b, etc., were prepared; SRM 8j is now available and represents the 10th renewal batch of 0.1 percent carbon Bessemer steel. While each of these batches differ somewhat in detailed analysis from one batch to another, all retain the relatively high level of phosphorus, sulfur, and nitrogen, and low alloy metal content characteristic of this type of material. It is not possible to supply preceding numbers of a renewal series when the stock is exhausted. If little demand exists or an alternate source of supply has become available for a material, production may be discontinued permanently or until sufficient justification is obtained to warrant renewal.

New Standard Reference Materials

New SRM's are issued from time-to-time, and are announced through semi-annual supplements to this Catalog, through the NBS Technical News Bulletin, through news releases to scientific, technical, and trade publications, as well as directly to prospective users.

Catalog Supplements

SRM Price and Availability Lists are usually prepared semi-annually to keep the Catalog current. These supplements provide a complete list of the available SRM's and their prices and provide descriptions of SRM's issued since the latest Catalog was printed.

Ordering

Orders should be addressed to the Office of Standard Reference Materials, Room B311, Chemistry Building, National Bureau of Standards, Washington, D.C. 20234. Telephone (301) 921-2045. Orders should give the amount (number of units), catalog number and name of the standard requested. For example: 1 each, No. 11h, Basic-Open-Hearth Steel, 0.2 percent C. These materials are distributed only in the units listed.

Acceptance of an order does not imply acceptance of any provision set forth in the order contrary to the policy, practice, or regulations of the National Bureau of Standards or the U.S. Government.

Orders received for "out-of-stock" materials are cancelled if only out-of-stock items are ordered. On other orders, shipment is made of available materials and out-of-stock items are cancelled. Back-orders are not accepted for out-of-stock materials; if a renewal lot of material is available, it will be furnished automatically.

Terms

Prices are given in the SRM Price and Availability List. These prices are subject to revision and orders will be billed for prices in effect at the time of shipment. New SRM Price and Availability Lists, when issued, are sent to users who have made purchases during the preceding twelve months, and to persons or organizations who request them. No discounts are given on purchases of Standard Reference Materials.

Remittances of the purchase price need not accompany purchase orders for firms or persons in the North American continent, agents in the United States of foreign firms, or foreign firms with established credit. Payment of invoices is expected within 30 days of receipt of an invoice. Payment on foreign orders may be made by any of the following:

- (a) UNESCO coupons,
- (b) banker's draft against U.S.A. bank,
- (c) bank to bank transfer to a U.S.A. bank,
- (d) letter of credit on a U.S.A. bank, or
- (e) by International Money Order.

Pro-forma invoice service will frequently require 6 to 8 weeks to process, and will be furnished only to those requiring such service, or when credit has not been established.

Domestic Shipments

Shipments of material (except for certain restricted categories, e.g., hydrocarbons, special nuclear materials, compressed gases, rubber, rubber compounding materials, and radioactive standards) intended for the United States, Mexico, and Canada are normally shipped prepaid (providing that the parcel does not exceed the weight limitations as prescribed by Postal Laws and Regulations) unless the purchaser requests a different mode of shipment, in which case the shipment will be sent collect. The Bureau does not prepay such shipping charges. Hydrocarbons, organic sulfur compounds, compressed gases, rubber, rubber compounding materials, radioactive standards, and similar materials are shipped express collect.

Foreign Shipments

Orders for small weight shipments will be shipped by prepaid International Air Parcel Post. Other shipments will be shipped prepaid International Parcel Post, except those shipments exceeding the parcel post weight limitations, which must be handled through an agent (shipping or brokerage firm) located in the U.S.A. as designated by the purchaser. Shipments handled through an agent will be packed for overseas shipment and forwarded via express collect to the U.S.A. firm designated as agent.

NOTE: Orders and inquiries submitted in English will be processed more rapidly than those requiring translations.

Introducción

Todo el Material de Referencia de Normas (Standard Reference Materials--SRM's), Material de Investigación (Research Materials--RM's), y Material General (General Materials--GM's), que figura en este Catálogo lleva nombres y números que lo identifican en forma permanente. Cada SRM, RM o GM designado en forma específica, tiene características idénticas a las de todas las demás muestras que llevan la misma designación, dentro de los límites que requiere el uso al cual está destinada; o, en caso necesario, se le asigna un número de serie y una calibración individual.

Los primeros SRM's emitidos por la Dirección fueron un grupo de minerales, hierros y aceros certificados en cuanto a su composición química, y se tomó la costumbre de llamarlos "muestras de normas" (standard samples). En la actualidad se cuenta con más de 800 SRM's que cubren una amplia gama de propiedades químicas y físicas, y la designación "Material de Referencia de Norma" (Standard Reference Material) es más apropiada. A medida que ha aumentado el número de SRM's, ha aumentado la diversidad de éstos, estableciéndose nuevos grupos, tales como normas para laboratorios clínicos, materiales nucleares, normas de viscosidad de vidrio, normas para goma y compuestos para goma, normas de colores, así como normas de grosor de revestimiento. Estos grupos figuran bajo los encabezamientos: Normas de Composición Química Certificada (Standards of Certified Chemical Composition), Normas de Propiedades Físicas Certificadas (Standards of Certified Physical Properties), Normas de Tipo de Ingeniería (Engineering Type Standards), Materiales de Investigación (Research Materials), o Materiales Generales (General Materials). Los grupos de materiales que corresponden a estos títulos generales aparecen en el Índice de Materias. Un índice alfabético indica la ubicación de determinado material. Un índice numérico ofrece la fecha del Certificado vigente emitido con los materiales en cuestión.

La lista detallada de materiales indica la certificación nominal para la cual se emite el material, pero es necesario consultar el Certificado para encontrar la certificación precisa. Se emite SRM's para los cuales no es factible suministrar valores numéricos o para los cuales tal certificación no sería útil. Estos SRM's ofrecen la seguridad de que son idénticas todas las muestras que llevan la misma designación, lo cual permite normalizar los procedimientos de prueba y referir a una base común los datos físicos o químicos de materiales desconocidos.

Norma de Renovación y Descontinuación Materiales de Referencia

Se tiene el propósito de completar la preparación de SRM's "de renovación" para cuando estén agotadas las existencias de cada clase de material, pero esto no es siempre posible. Por lo general, el material de renovación no será idéntico a su predecesor, pero será bastante similar, en especial en lo que a las características certificadas se refiere, y por lo general la renovación puede ser utilizada en lugar de su predecesor. Por ejemplo, en 1909, cuando se preparó el primer acero Bessemer con 0.1 por ciento de carbono, fue designado como "Standard Sample No. 8." Durante los años subsiguientes, fueron preparados lotes de renovación 8a, 8b, etc. En la actualidad está disponible el SRM 8j, que representa la décima renovación de acero Bessemer con 0.1 por ciento de carbono. Si bien cada uno de estos lotes difiere en algo en cuanto a análisis detallado entre lote y lote, todos ellos tienen el nivel relativamente alto de contenido de fósforo, azufre, nitrógeno y metal de baja aleación que es característico de este material. Una vez agotadas las existencias, no es posible suministrar números precedentes de una serie de renovaciones. Si existe poca demanda o se encuentra disponible una fuente alternativa de suministro de determinado material, es posible que la producción sea descontinuada en forma permanente, o hasta que se obtenga justificación suficiente para una renovación.

Nuevo Material de Referencia de Normas

De cuando en cuando se emite nuevos SRM's, los cuales son anunciados mediante suplementos semestrales de este Catálogo, por el NBS Technical News Bulletin, mediante informativos enviados a publicaciones científicas, técnicas y del ramo, así como directamente a los usuarios en potencia.

Suplementos de Catálogo.

Por lo general, las Listas de Precios y Disponibilidad de SRM's (SRM Price and Availability List) son preparadas semestralmente para mantener al día el Catálogo. Estos suplementos ofrecen una lista completa de los SRM's disponibles, así como sus precios, y dan la descripción de SRM's emitidos desde que fuera impreso el último Catálogo.

Pedidos

Los pedidos deberán ser hechos a la siguiente dirección:

Office of Standard Reference Materials
Room B311, Chemistry Building
National Bureau of Standards
Washington, D.C. 20234

El teléfono es (301) 921-2045. En los pedidos se deberá indicar el número de unidades, número y nombre del SRM. Por ejemplo: 1 muestra, No. 11h, Acero Siemens-Martin Básico, 0.2 por ciento C. Estos materiales son distribuidos solamente en las unidades que figuran en la lista.

Pedidos que se recibe por materiales "agotados" son cancelados y devueltos si se ha pedido solamente artículos agotados. En caso contrario, se remite los materiales disponibles y se cancela los materiales agotados. No se acepta pedidos retroactivos por materiales agotados; si un lote de material de renovación está disponible será suministrado en forma automática.

Condiciones

Los precios aparecen en la Lista de Precios y Disponibilidad de SRM's (SRM Price and Availability List). Estos precios están sujetos a cambio y se facturará los pedidos a los precios que estén en vigencia en el momento en que se hace el embarque. Al ser emitidas las nuevas Listas de Precios y Disponibilidad, son enviadas a los usuarios que han hecho compras en el curso de los doce meses precedentes, y a las personas u organizaciones que las solicitan. No se hace descuentos en las compras de Materiales de Referencia de Normas.

Órdenes de compra hechas por empresas o personas en el continente norteamericano, agentes en los Estados Unidos de empresas extranjeras, o empresas extranjeras de crédito bien establecido, no necesariamente tienen que ir acompañadas del valor de la compra. Se espera que las facturas sean pagadas a los treinta días de recibidas. El pago de pedidos desde el extranjero se puede hacer mediante cualquiera de los siguientes:

- (a) cupones de la UNESCO,
- (b) giro bancario contra un banco estadounidense,
- (c) transferencia de banco a un banco en los Estados Unidos,
- (d) carta de crédito en un banco en los Estados Unidos, o
- (e) Giro Internacional.

Suele requerir 8 semanas procesar facturas pro-forma, y serán suministradas solamente a aquellos que requieren de tal servicio, o cuando no ha sido establecido el crédito.

Remesas dentro del País

Los materiales (a excepción de ciertas categorías restringidas, por ejemplo hidrocarburos, materiales nucleares especiales, gases comprimidos, compuestos orgánicos de azufre y normas radioactivas) destinados a los Estados Unidos, México y el Canadá, por lo general son remesados con porte pagado (siempre y cuando que el paquete no exceda las limitaciones prescritas por las Leyes y el Reglamento Postal), a no ser que el comprador solicite una forma de embarque diferente, en cuyo caso se enviará con flete por pagar. Esta Dirección no paga por tales fletes. Hidrocarburos, compuestos orgánicos de azufre, gases comprimidos, materiales para compuestos de goma, normas radioactivas y material similar, son remitidos por expreso con flete por cobrar.

Embarques al Extranjero

Embarques de poco peso por valor de mas de \$100 dólares con porte pagado serán enviados por encomienda aérea con flete pagado. Los embarques que excedan las limitaciones de encomienda postal deberán ser manejados por intermedio de un agente (empresa de transportes o corredor) ubicado en los Estados Unidos, según indicación del comprador. Tales paquetes serán embalados para embarque marítimo y remesados como expreso por cobrar a la empresa en los Estados Unidos que ha sido designada como agente.

Pedidos que no han sido pagados serán embarcados con flete pagado por Encomienda Postal Internacional, sujeto a limitaciones en cuanto a tamaño, peso y categoría de material. Cualquier otra forma de embarque solicitado por un cliente deberá ser pagada por ese cliente. (Embarques que por cualquier motivo sean excluidos de remesa por Encomienda Postal Internacional deberán ser manejados por intermedio de un agente [empresa de transporte o corredor] ubicado en los Estados Unidos, según indicación del comprador.) Tales paquetes serán embalados para embarque marítimo y remesados como expreso por cobrar a la empresa en los Estados Unidos que ha sido designada como agente.

Nota: Pedidos y consultas presentados en idioma inglés serán tramitados en forma mas rápida que aquellos que requieren traducción.

Avant-Propos

Toutes les mesures étalon des matériaux classés dans la catégorie SRM, (Standard Reference Matériel) celles des matériaux de recherche ou RM (Research Material), et celles des matériaux divers ou GM (General Material) répertoriés dans ce catalogue se sont vus attribuer une désignation et un numéro de référence distincts afin de les identifier à titre permanent. Chacun des matériaux tombant dans les catégories SRM, RM ou GM est identique en ce qui concerne sa caractérisation à tout autre échantillon de même désignation, compte tenu des tolérances établies pour son usage. Si nécessaire, le matériau reçoit un numéro de série et un calibrage individuel.

Les titrages des premiers SRM établis par le Service des Poids et Mesures concernaient certains minerais, du fer et des aciers qui, par la suite, devaient être désignés sous le nom de "échantillon-type" (Standard samples). A l'heure actuelle, il existe plus de 800 SRM aux propriétés chimiques et physiques très variées. Il est donc apparu plus approprié de les designer sous le nom de "Mesures étalon pour matériaux répertoriés" (Standard Reference Materials). Au fur et à mesure que le nombre de SRM augmentait il en était de même de la diversité des nouvelles catégories de normes à établir pour, par exemple : les laboratoires médicaux, les matériaux nucléaires, la viscosité dynamique du verre, le caoutchouc et ses dérivés, les produits colorants et l'épaisseur des enduits et revêtements. Des rubriques spéciales ont été établies à cet effet: Normes de composition chimique (Standards of Certified Chemical Composition); Normes de propriétés physiques (Standards of Certified Physical Properties); Normes d'Engineering (Engineering Type Standards); Matériaux de recherche (Research Materials); ou Matériaux divers (General Materials). Les catégories de matériaux répertoriés sous ses rubriques sont indiquées à la table des matières. Il suffit de consulter l'index présenté par ordre alphabétique pour trouver la page se rapportant à un matériau particulier. Un index numérique fournit la date de délivrance du brevet correspondant aux normes exigées pour tel ou tel matériau.

La liste détaillée des matériaux constituant une homologation à caractère purement nominal, il est nécessaire de consulter le brevet délivré pour chacun des matériaux par le Service des Poids et Mesures. Il existe cependant un certain nombre de SRM pour lesquels il est impossible de fournir des valeurs numériques ou pour lesquels une telle certification serait inutile. Ces SRM sont identiques aux autres échantillons de même désignation. Ils permettent d'une uniformisation des méthodes expérimentales et dans le cas de matériaux inconnus de se baser sur des données de physique ou de chimie communes.

Mises à jour et discontinuation des normes. Matériaux repertoriés.

La mise à jour de la liste concernant les SRM doit en principe être terminée lorsque les stocks de chaque lot de matériaux sont épuisés, cela n'est toutefois pas toujours possible. Habituellement le nouveau lot n'est pas rigoureusement identique au précédent mais en est cependant très proche, ses spécifications étant les mêmes si bien qu'en règle générale il peut être utilisé en remplacement du lot précédent. Par exemple, lorsque le premier lot d'acier Bessemer à 0,1% de teneur en carbone a été préparé en 1909, il s'appelait à l'époque "Standard Sample No. 8". Les lots suivants furent désignés 8a, 8b, etc... Le SRM 8j., actuellement en stock correspond au dixième lot de lingots Bessemer fabriqués à partir d'une teneur en carbone de 0,1%. Bien que chacun de ces lots diffère quelque peu les uns des autres, si on les analyse en détail l'on constate qu'ils comportent tous une teneur relativement élevée en phosphore, soufre et azote ainsi qu'une faible proportion d'alliage, caractéristique à ce type de matériau. Pour les raisons exposées ci-dessus, il n'est donc pas possible de fournir des lingots appartenant à une série précédente épuisée. De surcroît, si la demande pour un type de matériau donné est marginale ou si une autre source d'approvisionnement est disponible pour tel ou tel matériau, la production peut alors en être suspendue à titre permanent ou jusqu'à ce qu'un renouvellement des stocks soit justifié.

Etablissement de normes pour nouveaux matériaux répertoriés.

Les mises à jour relatives aux SRM se font de temps à autre et sont publiées dans les suppléments semestriels de ce catalogue, dans le *NBS Technical News Bulletin* ou dans des articles paraissant dans des revues scientifiques, techniques ou commerciales ou encore directement portées à la connaissance de futures personnes intéressées.

Suppléments du Catalogue.

Les tarifs et listes de matériaux SRM disponibles sont (SRM Price and Availability List), en règle générale, publiés semestriellement afin que le catalogue demeure à jour. Ces suppléments fournissent une liste complète des SRM disponibles et des tarifs en vigueur ainsi que les spécifications des SRM repertoriés depuis la dernière publication du catalogue.

Modalités de commande.

Les lettres et bordereaux de commande doivent être adressés à : Office of Standard Reference Materials, Room B311, Chemistry Building, National Bureau of Standards, Washington, D.C. 20234. Téléphone: (301) 921-2045. Les commandes doivent spécifier la quantité désirée (nombre unitaire), le numéro de référence du catalogue ainsi que les spécifications et titrage. Exemple: Qte. 1, Ref.: no. 11h, Acier sur sole basique à teneur de 0,2% de carbone (en anglais, de préférence).

Les bordereaux de commande reçus pour des matériaux dont les stocks sont épuisés sont annulés et renvoyés à leur expéditeur si tous les articles demandés sont hors stock. Dans ce cas, seuls les matériaux disponibles sont expédiés. Les soldes de commande ne sont pas honorés pour ce qui concerne les matériaux dont les stocks sont épuisés, si un lot est renouvelé, le matériau est alors expédié automatiquement.

Conditions de paiement.

Les tarifs en vigueur relatifs aux SRM figurent à côté de la liste des matériaux disponibles (SRM Price and Availability List). Ces tarifs sont sujets à révision et les factures sont établies en fonction des tarifs applicables au moment de l'expédition de la commande. Les tarifs et listes de matériaux

en stock sont adressés aux clients ayant passé commande dans les douze mois ayant précédé l'établissement de ces tarifs et listes, ainsi qu'aux personnes ou organisations en sollicitant l'envoi. Aucun rabais n'est consenti sur les commandes de SRM.

Pour toute firme ou particulier, établi sur le continent nord américain, représentant aux Etats-Unis des firmes étrangères ou une société étrangère, offrant des garanties bancaires, il n'est pas nécessaire de joindre le règlement au bordereau de commande. Cependant les factures doivent être réglées dans les 30 jours suivant leur réception. Les commandes passées à l'étranger peuvent être réglées de la manière suivante:

- (a) Coupons UNESCO
- (b) Lettre de change à l'ordre d'une banque américaine
- (c) Transfert de banque à banque dans un établissement américain.
- (d) Lettre de crédit à l'ordre d'une banque américaine
- (e) ou par mandat international.

L'établissement de factures pro-forma entraîne fréquemment un délai de 6 à 8 semaines. Elles ne seront fournies que sur demande expresse des intéressés ou en l'absence d'une ouverture de crédit.

Expéditions à destination du continent nord-américain.

Les expéditions de matériaux (à l'exclusion toutefois de certaines catégories faisant l'objet de restrictions tels que : hydrocarbures, matériaux nucléaires spéciaux, gaz comprimés, composés de soufre organique, étalons pour matières radio-actives) à destination des Etats-Unis, du Mexique et du Canada sont en général réglées d'avance (à condition que le colis n'excède pas les limites des poids fixées aux termes des règlements postaux en vigueur) et sauf si le destinataire désire un mode d'envoi différent, auquel cas le colis est expédié contre remboursement, et en expresse.

Expéditions à l'étranger.

Les colis légers dont la valeur excède \$100 et dont le contenu a été réglé d'avance seront expédiés par la poste aérienne avec la mention "marchandise réglée". Les expéditions de matériaux dont le poids excède les limites fixées pour l'envoi de colis postaux doivent être faites par l'intermédiaire d'un agent maritime ou d'un transitaire ayant ses bureaux aux Etats-Unis et désigné par le destinataire. Les colis seront emballés spécialement pour envoi à l'étranger et expédiés en expresse et contre remboursement à la firme américaine agissant en qualité de Transitaire.

L'expédition des commandes non encore réglées se fera par colis postal, régime international, et compte tenu des restrictions imposées, sur le volume, le poids et les catégories de matériaux, par le Service des Postes. Tout autre mode d'expédition sollicité par un client doit être réglé directement par ses soins. (Les expéditions ne pouvant être faites pour une raison quelconque par colis postal, régime international, devront être confiées à un agent maritime ou à un transitaire dont les bureaux sont installés aux Etats-Unis et désigné par le destinataire. Ces colis seront spécialement emballés pour expédition à l'étranger en expresse et contre remboursement à la firme américaine agissant en qualité de transitaire.

N.B. Les commandes et demandes de renseignements rédigées en anglais seront satisfaites plus rapidement que celles exigeant une traduction.

Katalog der Standard-Nachweis-Materialien

Alle Standard-Nachweis-Materialien (Standard Reference Materials oder SRM's), Forschungs-Materialien (Research Materials oder RM's) und Allgemein-Materialien (General Materials oder GM's) die in diesem Katalog aufgeführt sind tragen unterschiedliche Namen und Nummern durch welche sie jederzeit identifizierbar sind. Jedes SRM, RM oder GM das eine gewisse Kennzeichnung

traegt hat dieselben Eigenschaften wie jede andere Probe versehen mit derselben Kennzeichnung innerhalb der Nachweisgrenzen die fuer einen bestimmten Zweck angegeben sind; andernfalls wird einer solchen Probe zweckbedingt eine Seriennummer als auch eine individuelle Eichung zugeteilt.

Die ersten vom Normamt (National Bureau of Standards oder NBS) herausgegebenen SRM's waren eine Reihe von Erzen, verschiedene Eisen und Staehle alle bescheinigt mit Bezug auf die jeweilige chemische Zusammensetzung, so dass im Laufe der Zeit diese Materialien als "Normproben" (Standard Samples) allgemein bekannt wurden. Zur Zeit sind ueber 800 SRM's verfuegbar die ein weites Gebiet chemischer und physikalischer Eigenschaften umfassen, auf Grund dessen die Bezeichnung "Standard-Nachweis-Material" sachgemaesser ist. Ebenso wie die Anzahl der SRM's zugenommen hat so ist auch derer Vielseitigkeit angestiegen: Neue Arten von Normproben wie z.B. klinische Labornormen, Kerntechnische Materialien, Normen zur Bestimmung der Glasviskositaet, Gummi- und Gummiverarbeitungs-Materialien, Farbnormen, als auch Normen zur Oberflaechenschichte-Dichtemessung, sind inzwischen herausgebracht worden. Derartige Gruppen von Standard-Nachweis-Materialien sind hierbei unter den folgenden Ueberschriften angegeben: Proben mit amtlich bescheinigter chemischer Zusammensetzung (Standards of Certified Chemical Composition), Proben mit attestierten physikalischen Eigenschaften (Standards of Certified Physical Properties), Kontroll-Proben mit Pruefungszeugnis fuer Anwendungen in der Technik (Engineering Type Standards), Forschungs-Materialien (Research Materials) und Allgemein-Materialien (General Materials). Die Gruppen von Materialien die diesen Kategorien zugehoeren sind in dem Inhaltsverzeichnis angegeben. Ein alphabetisch angeordneter Index gibt die Stelle jedes gewissen Materials oder auch jeder Gruppe aehnlicher Materialien an. Ein Nummerindex gibt das Datum des gueltigen Beglaubigungs-Zeugnisses an das jedem dieser Materialien zugeschrieben ist.

Die ausfuehrliche Anordnung dieser Materialien gibt den nominellen Richtwert an der dem entsprechenden Material zugesprochen ist; trotzdem muss in jedem einzelnen Falle das originale Beglaubigungs-Zeugnis fuer den rechtmassigen Richtwert in Anbetracht gezogen werden. Eine Anzahl von SRM's werden angeboten fuer die es nicht moeglich ist entweder einen Zahlenwert (Eichwert oder Richtwert) in dem Pruefungszeugnis anzugeben oder fuer die eine derartige Beglaubigung von fraglichem Vorteil ist. SRM's dieser Art gewaehrleisten die Identitaet aller Proben die mit ein und derselben Kennzeichnung versehen sind und demzufolge ermoeglichen die Normung von Versuchsverfahren als auch Bezugnahme auf eine gemeinsame Grundlage von physikalischen und chemischen Daten hinsichtlich unerforschter Materialien.

Erneuerung und Ausscheidung von Standard-Nachweis-Materialien

Die Vorbereitung von "Erneuerungs-" SRM's ist meistens dann abgeschlossen wenn bei einem gewissen Zeitpunkt der vorhandene Betrag eines jeden Materials erschoept ist; das ist jedoch nicht immer moeglich. Im allgemeinen ist das neuere Material dem aelteren Vorlaeufer nicht unbedingt identisch; trotzdem werden sie sich beiderseits sehr aehnlich sein besonders in Bezug auf die Eigenschaften fuer die das jeweilige SRM bescheinigt ist, so dass man im allgemeinen das neuere Material an Stelle des aelteren Materials benutzen kann. Zum Beispiel: Als der erste 0,1 Prozent Kohlenstoff Bessemer Stahl im Jahre 1909 vorbereitet wurde, erhielt diese Analysen-Kontrollprobe die Kennzeichnung "Standard Sample No. 8". Im Laufe der folgenden Jahre wurde eine Anzahl von Erneuerungsproben, gekennzeichnet 8a, 8b, u.s.w., vorbereitet. SRM 8j ist inzwischen erhaeltlich und dieses Material ist daher die zehnte Erneuerungsprobe des 0,1 Prozent Kohlenstoff Bessemer Stahls. Waehrend sich jede dieser Stahlschmelzen in Anbetracht der Richtanalyse von der anderen wenn auch nur gering unterscheidet, so enthaelt jede dieser Schmelzen den relativ hohen Gehalt an Phosphor, Schwefel und Stickstoff und einen niedrigeren Prozentsatz an Legierungsmetallen der kennzeichnend fuer diese Art Material ist. Es ist nicht moeglich vorangegangene Nummern einer gewissen SRM-Erneuerungs-Serie zu liefern sobald der Vorrat eines solchen Materials aufgebraucht ist. Sofern nur geringe Nachfrage besteht, oder eine andere Lieferungsquelle fuer eine bestimmte Analysen-Kontrollprobe verfuegbar ist, so kann die Vorbereitung bei NBS eines solchen Materials entweder ganz eingestellt oder voruebergehend unterbrochen werden bis ausreichender Grund vorhanden ist, eine Erneuerung desselben zu rechtfertigen.

Neue Standard-Nachweis-Materialien

Von Zeit zu Zeit werden neue SRM's herausgegeben die durch halbjaehrige Nachtraege zu diesem Katalog bekannt gemacht werden. Neuerscheinungen dieser Art werden ebenso zur Kenntnis gebracht durch das "NBS Technical News Bulletin" als auch durch Veroeffentlichungen in verschiedenen wissenschaftlichen, technischen Fachschriften und Fachhandelsblaettern sowie Ankuendigungen direkt gerichtet an Kunden und interessierte Kaeufer.

Katalog Nachtraege

Eine SRM Preisliste mit Lieferungsverzeichnis (SRM Price and Availability List) wird gewoehnlich halbjaehrlich veroeffentlicht um den Katalog auf dem laufenden zu halten. Derartige Nachtraege bieten ein vollstaendiges Verzeichnis aller verfuegbarer SRM's und deren Preise und enthalten Beschreibungen derjenigen SRM's die nach dem Druck des letzten Katalogs herausgegeben sind.

Bestellungen

Bestellungen muessen an die

Office of Standard Reference Materials
Chemistry Building, Room B311
National Bureau of Standards
Washington, D.C. 20234

gerichtet werden. Telefon: (301) 921-2045.

Bestellungen muessen die Anzahl (Menge jeder einzelnen Probe) und Bezeichnung (Nummer und Namen) des gewuenschten Standard-Nachweis-Materials angeben. Zum Beispiel: Ein, No. 11h, Siemens Martin Stahl, 0,2 Prozent C. Wir bitten, Bestellungen in Englisch zu erhalten, sofern moeglich. In diesem Falle lautet die obige Bestellung folgendermassen: "1 each, No. 11h, Basic-Open-Hearth Steel, 0.2 percent C." Diese Materialien sind nur in den beschriebenen Mengeinheiten erhaeltlich.

Erhalt einer Bestellung bei NBS schliesst keinerlei Vereinbarung ein mit irgendwelchen Bedingungen denen in besagter Bestellung Ausdruck gegeben ist, insofern als solche Bedingungen im Gegensatz zu den Vorschriften, ueblichen Handhabungen und Regulierungen des National Bureau of Standards oder der U.S. Regierung stehen.

Bestellungen fuer Materialien die aus dem Handel gezogen sind werden als ungueltig erklaert und an den Kunden zurueck geschickt sofern nur Bestellungen fuer nicht mehr erhaeltliche Proben bei NBS in Empfang genommen werden. Bei anderen Bestellungen werden die im Lager vorhandenen Materialien geliefert und die gegebenenfalls nicht mehr vorhandenen Materialien einfach gestrichen. Gleichermassen, Nachbestellungen von Materialien nicht mehr auf Lager koennen nicht angenommen werden; sobald eine Neuausgabe eines angeforderten Materials verfuegbar ist wird dasselbe unverzueglich geliefert.

Zahlungsbedingungen

Preise sind in der SRM Preis- und Lieferungsliste (SRM Price and Availability List) angegeben. Die genannten Preise werden von Zeit zu Zeit einer Beruecksichtigung unterzogen und Bestellungen werden berechnet zu den Preisen in Geltung zur Zeit des Versands. Neu veroeffentlichte SRM Preis- und Lieferungslisten werden an Kunden geschickt die in den vergangenen zwoelf Monaten ein Material angefordert haben. Diese Listen sind ebenso erhaeltlich bei Privatpersonen oder Betrieben die dafuer Antrag stellen. Wir moechten darauf hinweisen, dass wir keinen Rabatt auf Standard-Nachweis-Materialien irgend welcher Art gewaehrleisten koennen.

Zahlung braucht einer Bestellung nicht beizuliegen solange diese von Personen oder Firmen im Nordamerikanischen Kontinent kommt, oder diese von Vertretungen von auslaendischen Firmen in den Vereinigten Staaten, oder von auslaendischen Firmen mit etabliertem Kredit stammt.

Zahlung wird innerhalb von 30 Tagen nach Erhalt der Rechnung erwartet. Zahlungen fuer Bestellungen aus dem Ausland koennen wie folgt verschiedentlich gehandhabt werden:

- (a) durch UNESCO Coupons,
- (b) durch Bankbezug gegen eine Bank in den U.S.A.,
- (c) durch Bankueberweisung an eine Bank in den U.S.A.,
- (d) mittels Kreditbrief an eine Bank in den U.S.A., oder
- (e) durch eine internationale Geldueberweisung.

Pro-forma Abrechnungsdienst nimmt gewoehnlich 6 bis 8 Wochen in Anspruch und wird nur auf Antrag geleistet oder wenn Kredit nicht vorhanden ist.

Inlandversand

Versand von Materialien (mit Ausnahme von gewissen Beschraenkungen, wie z.B. der Versand von Kohlenwasserstoffen, besonderen Materialien der Kerntechnik, verdichteten Gasen, Gummi und Gummi-Verarbeitungsmaterialien, als auch radioaktive Probematerialien) innerhalb der Vereinigten Staaten, Mexiko und Kanada geschieht normalerweise unter Vorauszahlung (solange die Lieferung nicht entsprechende Gewichtsbegrenzungen der Postvorschriften ueberschreitet). Andernfalls, sollte der Kunde eine andere Weise des Versands verlangen, dann wird die Ware so versandt dass sie bei Empfang zahlbar ist. Das Normamt leistet keine Vorauszahlungen derartiger Versandspesen. Kohlenwasserstoffe, organische Schwefelverbindungen, verdichtete Gase, Gummi-Verarbeitungsmaterialien, radioaktive Stoffe und aehnliche Materialien werden per Express geliefert, mit Versandkosten zu Lasten des Empfaengers, zahlbar bei Erhalt.

Ueberseeversand

Lieferungen von nur geringem Gewicht werden als Versandspesen-vorausbezahlte Bestellung per internationale Paket-Luftpost befoerdert. Lieferungen anderer Art werden als Versandspesen-vorausbezahlte Ware durch internationale Paketpost befoerdert. Eine Ausnahme sind Lieferungen die gewisse Paketpost-Gewichtsbestimmungen ueberschreiten; diese muessen durch eine von dem Kunden beauftragte Agentur (Versand-oder Maklerfirma), ansaessig in den U.S.A., gehandhabt werden. Solche von einem Agenten vollzogenen Bestellungen die fuer den Ueberseeversand bestimmt sind werden entsprechend verpackt und via Express an den von dem Kunden angeforderten Vertreter in den U.S.A. weitergeleitet, zahlbar durch den Agenten an Erhalt der Lieferung.

Anmerkung: Bestellungen und Anfragen koennen schneller bearbeitet werden wenn sie auf Englisch erfolgen im Vergleich zu solchen, die eine Uebersetzung erfordern.



Guide for the Submission of Requests for the Development of New or Renewal Standard Reference Materials

August 20, 1964
(June 1, 1970 - Revised)

Introduction

The National Bureau of Standards presently has available more than 800 Standard Reference Materials. It is also working on the development of about 150 new ones and has on hand requests for the preparation of many others. The requests have always far exceeded the Bureau's capacity to produce and certify these materials.

Policy

One of the functions of the National Bureau of Standards is to develop, produce, and distribute Standard Reference Materials that provide a basis for comparison of measurements on materials and aid in the control of production processes in industry. To help carry out this function, the Office of Standard Reference Materials evaluates the requirements of science, industry, and government for carefully characterized reference materials, and directs their production and distribution. Emphasis is given to providing NBS Standard Reference Materials (a) where attainment of needed accuracy of analysis or accuracy of measurement of characteristics is not economically or technically feasible elsewhere, and where such accuracy is generally important to users, (b) where industry-wide standards for commerce are needed from a neutral supplier who is not otherwise available, and (c) where continuing availability of highly characterized material from a common source is important to science, industry, or government.

The National Bureau of Standards recognized and had responded to the need for broadening the present program on reference materials to include all types of well-characterized materials that can be used to calibrate a measurement system or to produce scientific data that can be readily referred to a common base. With this broadening, however, it still remains apparent that the demand for new Standard Reference Materials will continue to greatly exceed the Bureau's capacity for development. Therefore, requests for new Standard Reference Materials that will have limited use and for which the need is not very great will have to be deferred in favor of requests clearly showing a critical need. For the purpose of determining which requests are to receive top priority, the National Bureau of Standards will need, and will rely heavily upon, the information supplied by industry, either through its own representatives or through interested committees, such as those of the American Society for Testing and Materials, the American National Standards Institute, the International Organization for Standardization, etc.

Accordingly, while the Bureau welcomes all requests for the development of new Standard Reference Materials, it will help both the Bureau, and industry as well, if requests are accompanied by such information as will permit an assessment of the urgency and importance of proposed new reference materials.

Information Needed

Those requesting the development of new Standard Reference Materials should supply as much as possible of the following information:

- (1) Short title of Standard Reference Material.
- (2) Purpose for which the new standard material is needed.
- (3) Reasons why the new standard material is needed.
- (4) Special characteristics and/or requirements for the material. Include additional requirements and reasons, if more than one standard material is necessary for standardization in this area.
- (5) Your estimate of the possible present and future (10 year) demand for this new standard in your own operations and elsewhere.
- (6) Whether this standard, or a similar standard, can be produced by, or obtained from a source other than the National Bureau of Standards. If so, give reasons to justify its preparation by NBS.
- (7) Miscellaneous pertinent comments to aid justification for the new Standard Reference Material, such as: (a) an estimate of the range of application, monetary significance, and scientific and/or technological significance including, when feasible, estimates of the impact upon industrial productivity or growth, and (b) supporting letters from industry leaders, trade organizations, interested committees and others.



Other Services of the National Bureau of Standards

The following is a list of some of the services offered by NBS that may be of interest to SRM users. For general information see the entry on Technical Information and Publications.

Calibration and Test Services of the National Bureau of Standards

The measurement services of the National Bureau of Standards include the calibration of standards, test of instruments, and certain interlaboratory testing programs. These services are listed in NBS Special Publication 250, Calibration and Test Services of the National Bureau of Standards. [Available from the Superintendent of Documents, U.S. Government Printing Office, Washington, D.C. 20402, as SD Catalog No. C13.10:250 (1970 Edition) for \$2].

All technical inquiries regarding Calibration and Test Services, except those services related to electrical standards in the radio frequency region (above 30 kHz) or to time and frequency standards, should be directed to:

Office of Measurement Services
Room B154, Metrology Building
National Bureau of Standards
Washington, D.C. 20234
Telephone: (301) 921-2807.

Technical Inquiries related to electrical standards in the radio frequency range (above 30 kHz) or to time and frequency standards, should be directed to:

Office of Measurement Services
National Bureau of Standards
Boulder, Colorado 80302
Telephone: (303) 499-3753.

National Standard Reference Data System

The National Standard Reference Data System (NSRDS) is a nationwide program established to make critically evaluated data in the physical sciences available to the technical community. It publishes compilations of critically evaluated data, critical reviews and bibliographies. A complete listing of the publications of the NSRDS is available from the Office of Standard Reference Data (OSRD). The OSRD responds in a limited way to queries within the scope of the program by providing references, referrals, documentation, or data, as available. The program's monthly newsletter is available on request. Inquiries or requests for further information should be directed to:

Information Services
Office of Standard Reference Data
Room A523, Administration Building
National Bureau of Standards
Washington, D.C. 20234
Telephone: (301) 921-2583.

Standards Information Services Engineering and Product Standards

This service maintains a reference collection of some 20,000 engineering standards issued by more than 380 U.S. technical societies, professional organizations, and trade associations; specifications of the various State purchasing offices; standards and specifications of U.S. civilian government agencies; and the specifications and standards of the major foreign national and international standardizing bodies. The collection is open to the public Monday through Friday from 8:30 a.m. to 5 p.m.

The center publishes general and specialized indexes of standards. Information services consist of searching a Key-Word-In-Context (KWIC) Index to determine whether there are any published standards, specifications, test methods, or recommended practices for a given item or product.

Inquiries are referred to the appropriate source to obtain copies of standards. The Center neither sells nor distributes standards.

Inquiries or requests for additional information should be directed to:

Standards Information Service
Room B151, Technology Building
National Bureau of Standards
Washington, D.C. 20234
Telephone: (301) 921-2587.

Technical Information and Publications

The Office of Technical Information and Publications maintains a general correspondence and inquiry service on the technical activities of the National Bureau of Standards. Inquiries of a general nature and not covered by the services listed above should be directed to:

Office of Technical Information and Publications
Room A607, Administration Building
National Bureau of Standards
Washington, D.C. 20234
Telephone: (301) 921-2318.



Standards of Certified Chemical Composition

Steels (Chip Form)

These SRM's were prepared for the steel industry primarily for use with methods involving sample solutions in checking chemical methods of analysis both for production control and for customer acceptance. These SRM's consist of nominal composition steel alloys selected to provide a wide range of analytical values for the various elements that are of vital concern to the chemist. They are furnished in chips, usually sized between 16- and 40-mesh sieves, prepared from selected portions of commercial ingots.

The Certificate of Analysis, provided with each SRM, gives the chemical composition determined at NBS. Most certificates also include values obtained by other laboratories that cooperated in the certification of the SRM's.

Plain Carbon Steels (150 gram units unless otherwise noted)

Chemical Composition (Nominal Weight Percent)

SRM	Type	C	Mn	P	Grav	S Comb	Si
8j	Bessemer (Simulated), 0.1C	0.081	0.505	0.095	-----	0.077	0.058
10g	Bessemer, 0.2C240	.850	.086	0.109	.109	.020
11h	BOH, 0.2C200	.510	.010	-----	.026	.211
12h	BOH, 0.4C407	.842	.018	-----	.027	.235
13g	BOH, 0.6C61	.85	.006	-----	.030	.35 _s
14e	BOH, 0.8C753	.404	.008	-----	.039	.177
15g	BOH, 0.1C097	.485	.005	-----	.026	.095
16e	BOH, 1.1C	1.09	.381	.021	-----	.029	.20
19g	AOH, 0.2C	0.223	.554	.046	.032	.033	.186
20g	AISI 1045462	.665	.012	-----	.028	.305
51b	Electric Furnace 1.2C	1.21	.573	.013	.014	.014	.246
65d	Basic Electric 0.3C	0.264	.730	.015	.010	.010	.370
152a	BOH, 0.5C (Tin bearing)486	.717	.012	-----	.030	.202
178	Basic oxygen 0.4C395	.824	.012	-----	.014	.163
335	BOH 0.1C (Carbon only) 300g092	-----	-----	-----	-----	-----
337	BOH 1.1C (Carbon only) 300g	1.07	-----	-----	-----	-----	-----

Low Alloy Steels (150 gram units unless otherwise noted)

Chemical Composition (Nominal Weight Percent)

SRM	Type	(Other Forms)	C	Mn	P	Grav	S Comb.	Si	Cu	Ni
30f	Cr-V (SAE 6150)	-----	0.49	0.79	0.010	-----	0.010	0.28	0.76	0.071
32e	Ni-Cr (SAE 3140)	-----	.409	.798	.008	0.022	.021	.278	.127	1.19
33d	Ni-Mo (SAE 4820)	-----	.173	.537	.006	.010	.011	.253	.123	3.58
36b	Cr2-Mo1	-----	.114	.404	.007	-----	.019	.258	.179	0.203
72f	Cr-Mo (SAE X4130)	-----	.301	.545	.014	.024	.024	.256	.062	.055
100b	Manganese (SAE T1344)	-----	.397	1.89	.023	.029	.028	.210	.064	.030
105	High-Sulfur 0.2C (Carbon only)	-----	.193	-----	-----	-----	(.60)	-----	-----	-----
106b	Cr-Mo-Al (Nitalloy G)	-----	.326	0.506	.008	.016	.017	.274	.117	.217
125b	High-Silicon	1134	.028	.278	.029	-----	.008	2.89	.071	0.038
139a	Cr-Ni-Mo (AISI 8640)	-----	.404	.780	.013	.019	.019	0.241	.096	.510
155	Cr0.5-W0.5	-----	.905	1.24	.015	.010	.011	.322	.083	.100
361	AISI 4340	661,1095,1261	.383	0.66	.014	-----	.017	.222	.042	2.00
362	AISI 94B17 (Mod)	662,1096,1262	.160	1.04	.014	-----	.038	.39	.50	0.59
363	Cr-V (Mod)	663,1097,1263	.62	1.50	.029	-----	.009	.74	.10	.30
364	High Carbon (Mod)	664,1098,1264	.87	0.25 _s	.01 ₈	-----	.02 ₉	.06 ₅	.24 ₉	.14 ₄
366	Set: 1 each of 361,362,363,364, and 365, Electrolytic Iron. (SRM 365 is described on page 28.)									

SRM	B	Pb	Sb	Bi	Ag	Se	Te	Ce	La	Nd
361	(0.0005)	(<0.0001)	0.004 ₂	(0.0005)	(0.0004)	(0.004)	(0.0005)	(0.005)	(0.001)	(0.001)
362	(.0025)	(.0006)	.013	(.006)	(.0009)	(.001)	(.001)	(.002)	(.0005)	(.0005)
363	(.0010)	(.0018)	(.0022)	(.0006)	(.0032)	(.0001)	(.0023)	(.002)	(.0012)	-----
364	(.014)	.01 ₉	(.025)	(.002)	(.00005)	(.0003)	(.0002)	(.0005)	(.00007)	(.0002)



Plain Carbon Steels (Continued)

Cu	Ni	Cr	V	Mo	Co	Ti	Sn	Al (total)	N	Other	SRM
0.020	0.113	0.097	0.015	0.038	-----	-----	-----	-----	-----	-----	8j
.008	.005	.008	.007	.002	-----	-----	-----	-----	0.015	-----	10g
-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	11h
.073	.032	.074	.003	.006	-----	-----	-----	(0.038)	.006	-----	12h
-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	13g
.072	.053	.071	.002	.013	-----	-----	-----	.060	-----	-----	14e
-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	15g
-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	16e
.093	.066	.374	.012	.013	0.012	0.027	0.008	.031	-----	Nb 0.026	19g
.034	.034	.036	.002	.008	-----	-----	-----	.040	-----	-----	20g
.071	.053	.455	.002	.014	-----	-----	.008	-----	.011	-----	51b
.051	.060	.049	.002	.025	-----	-----	.004	.059	.013	Al ₂ O ₃ 0.009	65d
.023	.056	.046	.001	.036	-----	-----	.032	-----	-----	-----	152a
.032	.010	.016	.001	.003	-----	-----	-----	-----	-----	-----	178
-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	335
-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	337

Low Alloy Steels (Continued)

Cr	V	Mo	W	Co	Ti	As	Sn	Al Total	Nb	Ta	Zr	N	SRM
0.95	0.18	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	30f
.678	.002	0.023	-----	-----	-----	-----	(0.011)	-----	-----	-----	-----	0.009	32e
.143	.002	.246	-----	-----	-----	-----	-----	-----	-----	-----	-----	(.011)	33d
2.18	.004	.996	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	36b
0.891	.005	.184	-----	-----	-----	-----	-----	-----	-----	-----	-----	.009	72f
.063	.003	.237	-----	-----	-----	-----	-----	-----	-----	-----	-----	.004	100b
-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	105
1.18	.003	.199	-----	-----	-----	-----	-----	1.07	-----	-----	-----	-----	106b
0.019	-----	.008	-----	-----	-----	-----	.003	0.329	-----	-----	-----	-----	125b
.486	.003	.183	-----	-----	-----	-----	-----	-----	-----	-----	-----	.008	139a
.485	.014	.039	0.517	-----	-----	-----	-----	-----	-----	-----	-----	-----	155
.69	.011	.19	(.011)	0.030	0.02	0.01	.01	.02	0.02	(0.021)	0.01	(.0037)	361
.30	.040	.068	(.20)	.30	(.084)	(.079)	(.016)	(.086)	(.28)	(.20)	(.21)	(.0040)	362
1.31	.31	.028	(.047)	.04	(.06)	(.011)	(.094)	(.25)	(.049)	(.04)	(.048)	(.0042)	363
0.06 ₃	.10 ₅	.49	(.10)	.15	(.24)	(.057)	(.005)	(.014)	(.157)	(.11)	(.070)	(.003)	364

Ca	Mg	Zn	Pr	Ge	O	H	Au	Hf	SRM
(0.0001)	(0.0002)	(0.0005)	(0.0005)	(0.006)	(0.001)	(<0.0005)	(<0.00005)	(0.0002)	361
(.0003)	(.0007)	(.001)	(.0003)	(.002)	(.001)	(<.0005)	(<.00006)	(.0040)	362
(<.0001)	(.0005)	(.0004)	(.0005)	(.010)	(.0006)	(<.0005)	(.0006)	(.0042)	363
(.00005)	(.00005)	(.0005)	(.0001)	(.003)	(.0017)	(<.0005)	(.00007)	(.005)	364

High Alloy Steels (150 gram units unless otherwise noted)

Chemical Composition (Nominal Weight Percent)

SRM	Type	C	Mn	P	S		Si	Cu
					Grav	Comb		
126c	High-Nickel Steel (36% Ni)	0.026	0.47	0.004	----	0.006	0.19	0.040
131b	Low Carbon-Silicon (Carbon only) 100 grams0018	-----	-----	-----	-----	-----	-----
344	Cr15-Ni7-Mo2-Al169	.57	.018	-----	.019	.395	.106
345	Cr16-Ni4Cu3048	.224	.018	0.012	.012	.610	3.44
346	Cr22-Ni4-Mn9541	9.15	.018	-----	.063	.239	-----
348	Ni26-Cr15 (A286)044	1.48	.015	-----	.002	.54	0.22

Stainless Steels (150 gram units unless otherwise noted)

Chemical Composition (Nominal Weight Percent)

SRM	Type	(Other Forms)	C	Mn	P	S		Si	Cu
						Grav	Comb		
73c	Cr13 (SAE 420)		0.310	0.330	0.018	-----	0.036	0.181	0.080
121d	Cr17-Ni11-Ti0.3 (AISI 321) ..	1171	.067	1.80	.019	-----	.013	.54	.121
123c	Cr17-Ni11-Nb0.6 (AISI 348) .	1172	.056	1.7	.024	-----	.014	.59	.103
133a	Cr13-Mo0.3-S0.3120	1.03	.026	0.326	.330	.412	.118
160b	Cr19-Ni14-Mo3046	1.64	.020	-----	.018	.50	.172
166c	Low Carbon (AISI 3162) Carbon Only		.0078	-----	-----	-----	-----	-----	-----
339	Cr17-Ni9-Se0.2 (SAE 303Se) .		.052	0.738	.129	-----	.013	.654	.199
343	Cr16-Ni2 (SAE 431)150	-----	-----	-----	-----	-----	-----

Tool Steels (150 gram units unless otherwise noted)

Chemical Composition (Nominal Weight Percent)

SRM	Type	C	Mn	P	S		Si	Cu
					Grav	Comb		
50c	W18-Cr4-V1	0.719	0.342	0.022	0.010	0.009	0.311	0.079
132b	Mo-W-Cr-V86 ₅	.34 ₅	.01 ₃	----	.005	.18	.087
134a	Mo8-W2-Cr4-V1808	.218	.18	.007	.007	.323	.101
153a	Co8-Mo9-W2-Cr4-V2902	.192	.023	.007	.007	.270	.094

Steels (Granular Form – 100 gram units)

These granular-form SRM's are prepared by a pre-alloyed powder metallurgical process, which generally includes argon atomization and hydrogen annealing. The material normally is sized between 25 and 200 mesh sieves to ensure satisfactory homogeneity.

The Certificate of Analysis, provided with each of these SRM's gives the chemical composition as determined at NBS and values obtained by other laboratories that cooperated in the certification of the SRM's.

Chemical Composition (Nominal Weight Percent)

SRM	Type	Wt/Unit (grams)	C	Mn	P	S	Si	Cu	Ni
163	Low Alloy, 1.0 Cr	100	0.933	0.897	0.007	0.027	0.488	0.087	0.081
101f	Stainless, (AISI 304L)	100	.014	.087	.008	.008	.876	.030	9.96

SRM	Cr	V	Mo	W	Co	N	As	Sb	Ga
163	0.982	-----	0.029	-----	-----	0.007	-----	-----	-----
101f	18.49	0.034	.007	(0.0002)	0.088	-----	(0.003)	(0.0009)	(0.004)

High Alloy Steels (Continued)


Ni	Cr	V	Mo	Co	Ti	Al Total)	Nb	Ta	B	Fe	SRM
36.05	0.06 ₄	0.001	0.011	0.008	-----	-----	-----	-----	-----	-----	126c
7.28	14.95	0.040	2.40	-----	0.076	1.16	-----	0.002	-----	-----	131b
4.24	16.04	.041	0.122	.089	-----	-----	0.231	-----	-----	-----	344
3.94	21.61	.058	-----	-----	-----	-----	-----	N .441	-----	-----	345
25.8	14.54	.25	1.3	-----	2.24	0.23	-----	-----	0.0031	53.3	346
											348

Stainless Steels¹(Continued)

Ni	Cr	V	Mo	Co	Ti	Nb	Ta	Pb	Se	N	SRM
0.246	12.82	0.030	0.091	-----	-----	-----	-----	-----	-----	0.037	73c
11.17	17.4 ₃	-----	.165	0.10	0.342	-----	-----	-----	-----	-----	121d
11.3 ₄	17.4 ₀	-----	.22	.12	-----	0.65	<0.001	-----	-----	-----	123c
0.241	12.89	.026	.294	-----	-----	-----	-----	-----	-----	.032	133a
12.2 ₆	18.4 ₅	.047	2.38	.10 ₁	-----	-----	-----	0.001	-----	.03 ₉	160b
-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	166c
8.89	17.42	.058	0.248	.096	-----	-----	-----	-----	0.247	-----	339
2.14	15.76	.036	-----	-----	-----	-----	-----	-----	-----	.074	343

Tool Steels (Continued)

Ni	Cr	V	Mo	W	Co	Sn	As	N	SRM
0.069	4.13	1.16	0.082	18.44	-----	0.018	0.022	0.012	50c
.23	4.38	1.84	4.9 ₃	6.2 ₈	0.028	-----	-----	-----	132b
.088	3.67	1.25	8.35	2.00	-----	-----	-----	-----	134a
.168	3.72	2.06	8.85	1.76	8.47	-----	-----	.024	153a




Steels (Solid Form)

Several groups of SRM's have been prepared to meet the basic needs of the steel industry for analytical control primarily by optical emission and x-ray spectroscopic methods of analysis. Both nominal composition and analytical range SRM's are provided for ingot iron, low-alloy steel, stainless steel, tool steel, and specialty steel.

These SRM's are furnished in various forms. The 400 series is intended for optical emission spectroscopic methods of analysis utilizing the "point-to-point" technique. The 600 series is intended for microchemical methods of analysis such as electron probe microanalysis, spark source mass spectrometric analysis, and laser probe analysis. The 800, 1100, and 1200 series are intended for "point-to-plane" optical emission spectroscopic methods of analysis. The D800 series, and the 1100 and 1200 series also are intended for x-ray spectroscopic methods of analysis.

Because of the special homogeneity requirements, most of these materials have been prepared by using the most modern techniques of melting, casting, fabrication, and heat treatment to ensure adequate uniformity of composition. A Certificate of Analysis, which gives the chemical composition as determined at the NBS, is furnished for each SRM; many certificates also include values obtained by outside laboratories which cooperated in the certification of the SRM's. (Values in parentheses are not certified as they are based on the results from a single laboratory. Values in brackets are not certified but are nominal values obtained from heat analyses. These values are given for additional information on the chemical composition.)



Ingot Iron and Low-Alloy Steels

The preparation of these original spectroscopic SRM's began in about 1944 when the cores remaining after lathe cutting the materials for chip form standards were tested for homogeneity. Those found satisfactory were fabricated to the final shapes and sizes. To meet the urgent need in the mid-1950's for calibration standards for x-ray spectroscopic methods of analysis, portions of the material from five of these SRM's were converted to the applicable disk form. Although entirely satisfactory for conventional spectroscopic methods of analysis, these SRM's generally do not meet the stringent requirements for homogeneity necessary for use with the newer microchemical methods of analysis. These standards will be discontinued when the supply is exhausted.

SRM			Type	Chemical Composition (Nominal Weight Percent)				
7/32 in D × 4 in Long	1/2 in D × 2 in Long	1 1/4 in D × 1/4 in Disk		Mn	Si	Cu	Ni	Cr
----	803a	D803a	Acid Open Hearth, 0.6C	1.04	0.34	0.096	0.190	0.101
404a	804a	----	Basic Electric	0.88	.44	.050	.040	.025
405a	805a	D805a	Medium Manganese	1.90	.27	.032	.065	.037
407a	807a	D807a	Chromium-Vanadium	0.76	.29	.132	.169	.92
408a	808a	----	Chromium-Nickel76	.28	.10	1.20	.655
409b	809b	D809b	Nickel46	.27	.104	3.29	.072
----	810a	----	Cr2-Mo1	----	.36	.11	0.24	2.39
413	----	----	Acid Open Hearth, 0.4C67	.22	.25	.18	0.055
414	----	----	Cr-Mo (SAE 4140)67	.26	.11	.080	.99
417a	817a	----	Basic Open Hearth, 0.4C78	----	.13	.062	.050
418	----	----	Cr-Mo (SAE X4130)52	.28	----	.11	.96
418a	----	----	Cr-Mo (SAE X4130)52	.27	.040	.125	1.02
420a	820a	D820a	Ingot Iron017	----	.027	.0092	0.0032
----	821	----	Cr-W, 0.9C	1.24	----	.080	.10	.49
427	827	----	Cr-Mo (SAE 4150) (B only)	----	----	----	----	----

SRM			Chemical Composition (Nominal Weight Percent)						
7/32 in D × 4 in Long	1/2 in D × 22 in Long	1 1/4 in D × 1/4 in Disk	V	Mo	W	Co	Sn	Al Total	B
----	803a	D803a	0.005	0.033	----	----	-----	-----	-----
404a	804a	----	.002	.007	----	----	-----	-----	-----
405a	805a	D805a	-----	.005	----	----	-----	0.056	-----
407a	807a	D807a	.146	-----	----	----	-----	-----	-----
408a	808a	----	.002	.065	----	----	-----	-----	-----
409a	809a	D809a	.002	.009	----	0.025	0.012	-----	-----
----	810a	----	-----	.91	----	----	-----	-----	-----
413	----	----	.007	.006	----	----	-----	-----	-----
414	----	----	.003	.32	----	----	.014	.020	-----
417a	817a	----	-----	.013	----	----	.036	-----	-----
418	----	----	-----	.22	----	----	-----	-----	-----
418a	----	----	-----	.21	----	----	-----	-----	-----
420a	820a	D820a	-----	.0013	----	.006	.0017	.003	-----
----	821	----	.012	.040	0.52	----	-----	-----	-----
427	827	----	-----	-----	----	----	-----	-----	0.0027

Special Ingot Irons and Low-Alloy Steels

The planning of the 1100 series SRM's began in late 1952 to meet critical requirements of calibration in the iron and steel industry. Steel for these SRM's was prepared by the most modern melting, casting, and fabrication techniques to provide large quantities of material of the highest possible homogeneity. The materials were fully characterized and included investigations by means of electron probe microanalysis and quantitative metallographic techniques. It was concluded that, for example, SRM's 461 and 463 are sufficiently homogeneous that any present microanalytical technique can be carried out with little chance of inaccuracy caused by inhomogeneity. Details of the metallographic and homogeneity characterization are given in NBS Miscellaneous Publication 260-3 and 260-10, respectively (see inside back cover for ordering instructions).

These SRM's were first issued in 1957 and they have been in great demand ever since. Several years ago it became apparent that supplies of some of the 1100 series would be exhausted. To remedy this situation, the new 1200 series (SRM's 1261-1265) was prepared.

The 1200 Series consists of four low alloy steels and an electrolytic iron containing a graded series of 40 elements. Material from the same melts are available in three other forms: chip form, 361-365, for chemical methods of analysis, (pages 16 and 28); rods, 661-665, 3.2 mm (1/8 in) in diameter and 51 mm (2 in) long for microchemical methods of analysis such as electron probe microanalysis, spark source mass spectrometric analysis, and laser probe analysis (see below); and rods (1095-1099), 6.4 mm (1/4 in) in diameter and 102 mm (4 in) long for determining gases in metals by vacuum fusion and neutron activation methods of analysis (page 27).

SPECIAL INGOT IRONS AND LOW ALLOY STEELS

Chemical Composition (Nominal Weight Percent)

Sizes: 400 Series — 7/32 in D × 4 in long
 600 Series — 3.2 mm D × 51 mm long
 1100 Series — 1 1/4 in D × 3/4 in thick
 1200 Series — 31 mm D × 19 mm thick

Set s: 666
 667
 668
 1266

Set of 2: 661 and 665
 Set of 2: 662 and 663
 Set of 5: 661, 662, 663, 664, and 665
 Set of 5: 1261, 1262, 1263, 1264, and 1265

SRM	Type	(Other Forms)	C	Mn	P	S	Si	Cu	Ni	Cr
1134	High Silicon	125b	0.026	0.277	0.028	0.009	2.89	0.070	0.038	0.019
1135	High Silicon	(179)	.027	.094	.006	.026	3.19	.056	.050	.022
461	Low Alloy A15	.36	.053	(.02)	0.047	.34	1.73	.13
462	Low Alloy B40	.94	.045	(.02)	.28	.20	0.70	.74
463	Low Alloy C19	1.15	.031	(.02)	.41	.47	.39	.26
464	Low Alloy D54	1.32	.017	(.02)	.48	.094	.13 _s	.078
465	1165 Ingot Iron E037	0.032	.008	(.01)	.029	.019	.026	.004
466	1166 Ingot Iron F065	.11 _s	.012	(.01)	.025	.033	.051	.011
467	1167 Low Alloy G11	.27 _s	.033	(.01)	.26	.067	.088	.036
468	Low Alloy H26	.47 _s	.023	(.02)	.075	.26	1.03	.54
661	1261 AISI 4340	361,1095	.38 ₂	.66	.015	.017	.223	.042	1.99	.69
662	1262 AISI 94B17 (Mod) ...	362,1096	.16 ₀	1.04	.042	.038	.39	.50	0.59	.30
663	1263 Cr-V (Mod)	363,1097	.62	1.50	.02 _s	.008	.74	.09 _s	.32	1.31
664	1264 High Carbon (Mod) ...	364,(1098)	.87 ₀	0.25 _s	.01 _s	.028	.067	.24 _s	.14 _s	0.06 _s
665	1265 Electrolytic Iron	365,1099	.0067	.0057	.002 _s	.0059	.008 ₀	.0058	.041	.007 ₂

SRM	B	Pb	Ag	Ge	O	N	H
1134	-----	-----	-----	-----	-----	-----	-----
1135	-----	-----	-----	-----	-----	-----	-----
461	0.000 ₂	(0.003)	(0.001 _s)	(0.001 _s)	(0.02 _s)	(0.00 ₆)	-----
462	.000 _s	.006	(<.0002)	(.003 ₀)	(.006)	(.00 ₈)	-----
463	.0012	.012	(<.0002)	(.002 _s)	(.007)	(.00 ₆)	-----
464	.005	.020	(.003 _s)	(.001 _s)	(.006)	(.00 ₇)	-----
465	.000 ₁	(<.0005)	(.0002 _s)	(.003 _s)	(.003)	(.00 ₅)	-----
466	(.000 ₂)	(.001 _s)	(.0004 _s)	(.003 ₀)	(.005)	(.00 ₆)	-----
467	(.000 ₂)	.000 ₆	(.004 ₀)	(.003 ₀)	(.004)	(.00 ₄)	-----
468	.009	(<.0005)	(<.0005)	(.001 ₀)	(.004)	(.00 ₆)	-----
661	.0005	.00002 _s	.0004	[.006]	(.0009) [†]	(.0037) [†]	[<0.0005] [†]
662	.0025	.0004 _s	(.0010)	[.002]	(.0011) [†]	(.0041) [†]	[<.0005] [†]
663	.0009 ₁	.0022	(.0038)	[.010]	(.007) [†]	(.0041) [†]	[<.0005] [†]
664	.011	.024	(.00002)	[.003]	[.0017] [†]	[.003] [†]	[<.0005] [†]
665	.00013	.00001 _s	(~.000002)	(~.0014)	(~.0063) [†]	(~.0011) [†]	(~.0001) [†]

SRM		Sb	Bi	Ca	Mg	Se	Te
661	1261	0.004 ₂	0.0004	(<0.0001)	(0.0001)	0.004	0.0006
662	1262	.012	(.002)	(.0002)	(.0006)	[.001]	(.0005)
663	1263	.001 ₆	(.0008)	(<0.0001)	(.0005)	[.0001]	(.0022)
664	1264	(.035)	(.0009)	(<0.0001)	(.0001)	[.0003]	[.0002]
665	1265	-(<.00005)	-(<.00001)	-(<.00001)	-(<.00002)	-(<.00001)	-(<.00001)

SRM		Zn	Au	Ce	Hf	La	Nd	Pr	Fe
661	1261	(0.0001)	(<.0.00005)	0.001 ₃	[0.0002]	0.0004	0.0003	(0.00014)	(95.6)
662	1262	0.0005	(<.0.00005)	(.0011)	[.006]	.0004	(.0005)	(.00012)	(95.3)
663	1263	(.0004)	.0005	(.0016)	[.0015]	.0006	(.0007)	(.00018)	(94.4)
664	1264	.001	.0001	(.00025)	[.005]	.00007	(.00012)	(.00003)	(96.7)
665	1265	(.0001)	-(<.000002)	-(<.000005)	-(<.00002)	-(<.000005)	-(<.000005)	-(<.000005)	(99.9)

(Continued)

V	Mo	W	Co	Ti	As	Sn	Al Total	Nb	Ta	Zr	
-----	0.008	-----	-----	-----	-----	0.003	0.329	-----	-----	-----	1134
<0.01	.014	-----	-----	-----	-----	.004	.0028	-----	-----	-----	1135
.024	.30	0.012	0.26	(0.01)	0.028	.022	.005	0.011	0.002	(<0.005)	461
.058	.080	.053	.11	.037	.046	.066	.02 ₃	.096	.036	.063	462
.10	.12	.10 ₅	.01 ₃	.010	.10	.013	.02 ₇	.19 ₅	.15	.20	463
.29 ₅	.029	.022	.02 ₈	.004	.018	.043	.005	.037	.069	.010	464
.002	.005	(.001)	.008	.20	.010	.001	.19	(.001)	.001	(.002)	465 1165
.007	.011	(.006)	.04 ₆	.057	.014	.005	.01 ₅	.005	.002	(<.005)	466 1166
.041	.021	.20	.07 ₄	.26	.14	.10	.16	.29	.23	.094	467 1167
.17	.20	.077	.16	.011	.008	.009	.04 ₂	.006	.005	(<.005)	468
.011	.19	.01 ₅	.030	.020	.017	.01 ₁	.02 ₁	.022	.020	.009	661 1261
.04 ₁	.06 ₈	.21	.30	.084	.09 ₂	.016	.09 ₅	.29	.20	.19	662 1262
.31	.030	.04 ₅	.048	.050	.010	(.095)	.24	.049	(.053)	.049	663 1263
.10 ₅	.49	.10	.15	.24	.05 ₂	[.005]	(.008)	.15 ₇	.11	.068	664 1264
.0006	.0050	(~.00004)	.007 ₀	.0006	(.0002)	(~.0002)	(.0007)	(.00001)	-(<.00005)	-(<.00001)	665 1265

NOTE: Values in parentheses not certified, based on a single analytical method.

Values in brackets not certified, approximate values from the heat analyses.

† From Gasometric Certificates: SRM's 1095 through 1099.

- Not detected, value given is conservative "Upper Limit" of detection by specific methods of analysis.

Stainless Steels

Three groups of stainless steel SRM's designed primarily for calibration in spectroscopic methods of analysis are available.

Group I, SRM 442 through 444, consists of three SRM's of the 18 Cr - 8 Ni type stainless steel available only in rod form for use with the "point-to-point" technique in emission spectroscopy.

Group II is comprised of six standards, each available in three different physical forms; the 400, the 800, and the D800 series.

Both Group I and Group II standards have been extensively tested for homogeneity and found satisfactory for application in conventional spectroscopic methods of analysis. Neither group, however, has been tested for microanalytical methods and their use in these applications is not recommended.

Group III consists of six stainless steels available only in disk form and two in both chip and disk form for the "point-to-plane" technique of emission spectroscopy and for x-ray spectroscopy. They were prepared by melting, casting, and fabrication techniques known to produce material of high homogeneity.

(Values in parentheses are not certified, but are given for additional information only).

GROUP I

SRM		Chemical Composition (Nominal Weight Percent)									
7/32 in D× 4 in Long	Name	Mn	Si	Cu	Ni	Cr	V	Mo	W	Co	
442	Cr16-Ni10	2.88	(0.09)	0.11	9.9	16.1	0.032	0.12	(0.08)	0.13	
443	Cr18.5-Ni9.5	3.38	(.15)	.14	9.4	18.5	.064	.12	(.09)	.12	
444	Cr20.5-Ni10	4.62	(.65)	.24	10.1	20.5	.12	.23	(.17)	.22	


GROUP II

SRM				Chemical Composition (Nominal Weight Percent)						
7/32 in D × 4 in Long	1/2 in D × 2 in Long	1 1/4 in D × 1/4 in Disks		Name	Mn	Si	Cu	Ni	Cr	V
445	845	D845		Cr13-Mo0.9 (Mod. AISI 410)	0.77	0.52	0.065	0.28	13.31	(0.05)
446	846	D846		Cr18-Ni9 (Mod. AISI 321)53	1.19	.19	9.11	18.35	(.03)
447	847	D847		Cr24-Ni13 (Mod. AISI 309)23	0.37	.19	13.26	23.72	(.03)
448	----	D848		Cr9-Mo0.3 (Mod. AISI 403)	2.13	1.25	.16	.52	9.09	(.02)
449	849	D849		Cr5.5-Ni6.5	1.63	0.68	.21	6.62	5.48	(.01)
450	850	D850		Cr3-Ni25	----	.12	.36	24.8	2.99	(.006)

GROUP III

Chemical Composition (Nominal Weight Percent)

SRM	Type	(Other Forms)	C	Mn	P	S	Si	Cu	Ni	Cr
1152	Cr18-Ni10	----	0.163	1.19	0.017	0.017	0.654	0.497	10.21	18.49
1154	Cr19-Ni10	----	.094	1.74	.038	.033	1.09	.560	10.25	19.58
1155	Cr18-Ni12-Mo2 (AISI 316)	----	.046	1.63	.020	.018	0.50	.169	12.18	18.45
1185	Cr17-Ni13-Mo2 (AISI 316)	----	.11	1.22	.019	.016	.40	.067	13.18	17.09
1171	Cr17-Ni11-Ti0.3	121d	.067	1.8 ₀	.018	.01 ₃	.54	.121	11.2	17.4
1172	Cr17-Ni11-Nb0.6	123c	.056	1.7 ₆	.025	.01 ₄	.59	.10 ₅	11.3 ₅	17.4 ₀



GROUP I (Continued)

Chemical Composition (Nominal Weight Percent)								SRM
Ti	Sn	Nb	Ta	B	Pb	Zr	Zn	7/32 in D × 4 in Long
0.002	0.0035	0.032	(0.0006)	0.0005	0.0017	(0.004)	(.003)	442
.003	.006	.056	(.0008)	.0012	.0025	-----	(.005)	443
.019	.014	.20	(.004)	.0033	.0037	(.011)	(.004)	444

GROUP II (Continued)

Chemical Composition (Nominal Weight Percent)						SRM		
Mo	W	Ti	Sn	Nb	Ta	7/32 in D × 4 in Long	1/2 in D × 2 in Long	1 1/4 in D × 1/4 in Disks
0.92	(0.42)	(0.03)	----	0.11	(0.002)	445	845	D845
.43	(.04)	(.34)	(0.02)	.60	(.030)	446	846	D846
.059	(.06)	(.02)	----	.03	(.002)	447	847	D847
.33	(.14)	(.23)	(.05)	.49	(.026)	448	----	D848
.15	(.19)	(.11)	(.07)	.31	(.021)	449	849	D849
----	(.21)	(.05)	(.09)	.05	(.002)	450	850	D850

GROUP III (Continued)

V	Mo	Co	Ti	As	Sn	Al	Nb	Ta	B	Pb	Zr	SRM
0.044	0.366	(0.095)	(0.12)	(0.01)	(0.004)	(0.003)	(0.20)	(0.085)	(0.005)	(0.001)	(0.03)	1152
.061	.463	(.12)	(.48)	(.03)	(.023)	(.035)	(.26)	(.045)	(.0006)	(.012)	(.022)	1154
.047	2.38	.101	-----	-----	-----	-----	-----	-----	-----	.001	-----	1155
-----	2.01	-----	<.001	-----	-----	-----	<.001	<.001	-----	-----	-----	1185
-----	0.16 _s	.10	.34	-----	-----	-----	-----	-----	-----	-----	-----	1171
-----	.22	.12	-----	-----	-----	-----	0.65	<.001	-----	-----	-----	1172

Tool Steels

A group of six high-speed tool steel SRM's is available in three different physical forms. A wide concentration range is covered by combining the concentration ranges of three American Iron and Steel Institute (AISI) designations with three other tool steels of tailored composition.

SRM			Type	Chemical Compositions (Nominal Weight Percent)							
7/32 in D× 4 in Long	1/2 in D× 2 in Long	1 1/4 in D× 1/4 in Disk		Mn	Si	Cu	Cr	V	Mo	W	Co
436		D836	Special (Cr6-Mo3-W10)	0.21	0.32	0.075	6.02	0.63	2.80	9.7	----
437	837	D837	Special (Cr8-Mo2-W3-Co3)48	.53	-----	7.79	3.04	1.50	2.8	2.9
438	838	D838	Mo High Speed (AISI-SAE-M30)20	.17	.17	4.66	1.17	8.26	1.7	4.9
439	839	D839	Mo High Speed (AISI-SAE-M36)18	.21	.12	2.72	1.50	4.61	5.7	7.8
440	840	D840	Special W High Speed (Cr2-W13-Col 12)15	.14	.059	2.12	2.11	0.070	13.0	11.8
441	841	D841	W High Speed (AISI-SAE-TI)27	.16	.072	4.20	1.13	.84	18.5	----

Maraging Steel

This alloy derives its name from the formation of martensite on age hardening. They attain remarkable metallurgical properties by a simple heat treatment. Extensive use of these alloys is expected, particularly in submarines, missiles, and aircraft. This Maraging Steel, SRM No. 1156, of the 19 percent nickel type, is designed primarily for calibration in optical emission and x-ray spectroscopic methods of analysis.

SRM		Chemical Compositions (Nominal Weight Percent)							
1 1/4 in D× 3/4 in Disk	Type	C	Mn	P	S	Si	Cu		
1156	Maraging, (Ni 19)	0.023	0.21	0.011	0.012	0.184	0.025		
	Ni	Cr	Mo	Co	Ti	Al	Zr	B	Ca
1156	19.0	0.20	3.1	7.3	0.21	0.047	0.004	0.003	<0.001

High-Temperature Alloys (Solid Form)

High-temperature alloy SRM's were prepared to meet the critical needs of industry, particularly the aerospace industry, and government agencies. These SRM's are useful in instrument calibration, primarily for x-ray and optical emission spectroscopic methods of analysis.

Chemical Composition (Nominal Weight Percent)								
SRM	Type	C	Mn	P	S	Si	Cu	
1206-2	Rene-41	0.21 ₇	0.030	(0.004)	0.006	0.21 ₆	0.040	
1207-1	Waspaloy(1)043	.34	.005	.009	.47 ₂	.026	
1207-2	Waspaloy(2)083	.29 ₅	.005	.009	.61 ₅	.033	
1208-1	Inco 718(1)046	.38 ₅	.003	.01 ₁	.43 ₄	.14 ₇	
1208-2	Inco 718(2)022	.23 ₀	.003	.007	.08 ₃	.077	

1209 Set, 1 each: 1206-2, 1207-1, 1207-2, 1208-1, and 1208-2


SRM	Ni	Cr	Mo	Co	Ti	Al	Nb	Ta	Fe
1206-2	53.3	19.7	10.3	11.5 ₅	2.9	1.7 ₄	----	-----	0.46
1207-1	56.1	18.88	4.50	13.0 ₅	3.09	1.26	----	-----	2.22
1207-2	55.7	19.4 ₄	4.34	13.5 ₀	2.54	1.3 ₈	----	-----	2.09
1208-1	51.9	17.5	3.2 ₄	0.82	0.46	(0.15)	5.3 ₈	(0.012)	19.2
1208-2	51.5	17.4	3.13	0.76	(0.8 ₅)	(0.8 ₅)	4.9 ₈	(0.012)	19.8

Oxygen Standards

These SRM's are issued for the determination of oxygen and nitrogen by vacuum fusion, inert gas fusion, and neutron activation methods. SRM's 1095-1099 were prepared from the same melt as the 1200 series (1261-1265), page 22.

(Values in parentheses are not certified, but are given for information only.)

SRM	Type	Unit	Oxygen (ppm)	Nitrogen (ppm)
1090	Ingot iron	Rods 1/4 in D × 4 in Long	491	(60)
1091	Stainless Steel (AISI 431)	Rods 5/16 in D × 4 in Long	131	(945)
1092	Vacuum-melted steel	Rods 1/4 in D × 4 in Long	28	(4)
1093	Valve steel	Rods 1/4 in D × 4 in Long	60	(4807)
1094	Maraging steel	Rods 1/4 in D × 4 in Long	4.5	(71)
1095	AISI 4340 Steel	Rods 6.4 mm in D × 102 mm Long	9	(37)
1096	AISI 94B17 (Mod) Steel	Rods 6.4 mm in D × 102 mm Long	10.7	40.4
1097	Cr-V (Mod) Steel	Rods 6.4 mm in D × 102 mm Long	6.6	(41)
1098	High Carbon (Mod) Steel	Rods 6.4 mm in D × 102 mm Long	(17)	(30)
1099	Electrolytic Iron	Rods 6.4 mm in D × 102 mm Long	61	(13)
1089	Set of 5: 1095, 1096, 1097, 1098, and 1099			



Cast Irons (Chip Form — 150 gram units, unless otherwise noted)

This group of cast iron SRM's is similar to the chip-form steels and was prepared for use in checking chemical methods in the cast iron industry. These SRM's are in the form of chips, usually sized between 16- and 25-mesh sieves. They are prepared by lathe cutting of chips with a multiple-tooth cutting tool from thin-wall cylindrical castings especially made for the purpose. Supplied with each SRM is a Certificate of Analysis listing the chemical composition determined at NBS and other laboratories that cooperated in the certification of the SRM's. For SRM 365, Electrolytic Iron, the Certificate provides information on these additional elements: W, Nb, Ag, Zn, Ge, O, H, Ta, Nd, Zr, Sb, Bi, Ca, Mg, Se, Te, Ce, La, Pr, Au, Hf, and Fe.

(Values in parentheses are not certified, but are given for information only.)

SRM	Type	Wt/Unit (grams)	Chemical Composition (Nominal Weight Percent)							
			C		Mn	P	S		Si	Cu
			Total	Graphitic			Grav.	Comb.		
3b	White	110	2.44	----	0.353	0.086	-----	0.088	1.04	0.050
4j	Cast	150	2.99	2.38	.79	.17	-----	.062	1.31	.24
5L	Cast	150	2.59	1.99	.68	.280	-----	.123	1.83	1.01
6g	Cast	150	2.84	2.00	1.06	.56	-----	.123	1.06	0.50
7g	Cast (High Phosphorus) ...	150	2.69	2.59	0.612	.794	0.061	.060	2.41	.128
55e	Ingot	150	0.0112	----	.035	.003	.012	.011	0.001	.065
82b	Cast (Ni-Cr)	150	2.85	2.37	.745	.025	-----	.007	2.10	.038
107b	Cast (Ni-Cr-Mo)	150	2.75	1.87	.510	.058	.067	.067	1.35	.235
115a	Cast (Cu-Ni-Cr)	150	2.62	1.96	1.00	.086	.064	.065	2.13	5.52
122e	Cast (Car Wheel)	150	3.51	2.78	0.528	.349	-----	.074	0.510	0.033
341	Ductile	150	1.81	1.23	.92	.024	.007	.007	2.44	.152
342	Nodular	150	2.45	2.14	.369	.020	.014	.014	2.85	.14
342a	Nodular	150	1.86	1.38	.275	.018	-----	.006	2.73	.14
365	Electrolytic Iron	150	0.0070	----	.0057	.003	-----	.006	0.007	.0058


Cast Steels, White Cast Irons, Ductile Irons and Blast Furnace Irons (Solid Form)

These chill-cast SRM's were prepared for use in analytical control of cast steels and cast irons by rapid instrumental methods. Although employed in x-ray spectroscopic analysis, they are particularly useful for calibrating vacuum optical emission spectrometers because they permit the determination of carbon, phosphorus, and sulfur in addition to the metallic elements.

These SRM's are chill-cast sections. Details of the preparation and intended use of the SRM's are given in the NBS Miscellaneous Publication 260-1. (See inside back cover for ordering instructions.)

(Values in parentheses are not certified, but are given for information only.)

SRM	Type	Chemical Composition (Nominal Weight Percent)									
		C	Mn	P	S	Si	Cu	Ni	Cr	V	Mo
1 1/4 in × 1 1/4 in × 1/2 in											
1147	White (4i)	3.60	0.78	.160	.059	1.31	0.23	0.070	0.093	.032	0.078
1148	White (5L)	2.89	.66	.300	(.11)	1.82	.99	.091	.146	.036	.022
1149	White (6g)	3.28	1.05	.564	.127	1.04	.49	.138	.363	.055	.036
1140	Ductile (No. 1)	3.18	0.725	.0070	.010	1.92	.10	.028	.030	.030	.090
1141	Ductile (No. 2)	3.64	.480	.072	.020	1.11	.21	.54	.145	.0090	.05
1142	Ductile (No. 3)	2.94	.18	.20	.015	3.33	1.02	1.65	.053	.006	.022
1138	Cast Steel (No. 1)	0.120	.43	.053	.053	0.34	0.09	0.10	.12	.020	.05
1139	Cast Steel (No. 2)792	.98	.011	.013	.85	.40	.93	1.96	.24	.51
1143	Blast Furnace (No. 1) ...	3.91	.414	.158	.028	1.68	.144	.115	0.145	.008	(.005)
1144	Blast Furnace (No. 2) ...	4.27	1.33	.112	.021	0.276	.090	.021	.019	.004	.007



(Continued)

Chemical Composition (Nominal Weight Percent)

Ni	Cr	V	Mo	Co	Ti	As	Sn	Al (Total)	Mg	N	SRM
0.010	0.052	0.006	0.002	-----	-----	-----	-----	-----	-----	-----	3b
.068	.09	.03	.080	-----	0.05	0.03	-----	-----	-----	-----	4j
.086	.15	.036	.020	-----	.05	<.005	-----	-----	-----	0.006	5L
.136	.37	.06	.035	-----	.06	.04	-----	-----	-----	.006	6g
.120	.048	.010	.012	-----	.044	.014	-----	-----	-----	.004	7g
.038	.006	<.001	.011	0.007	-----	.007	0.007	0.002	-----	.004	55e
1.22	.333	.027	.002	-----	.027	-----	-----	-----	-----	-----	82b
2.12	.560	.008	.750	-----	.016	-----	-----	-----	-----	(.008)	107b
14.49	1.98	.014	.050	-----	.020	-----	-----	-----	-----	-----	115a
0.080	(0.038)	(.032)	(.001)	-----	(.026)	(.018)	-----	-----	-----	(.009)	122e
20.32	1.98	.012	.010	-----	.018	-----	-----	-----	0.068	-----	341
0.023	0.032	.005	.009	-----	.019	-----	-----	-----	.053	-----	342
.06	.034	-----	-----	-----	.020	-----	-----	-----	.069	-----	342a
.041	.0072	.0006	.0050	.0070	.0006	.0002	(.0002)	(.0007)	N 0.001	Pb 0.00002	365



(Continued)

Chemical Composition (Nominal Weight Percent)

														SRM
Co	Ti	As	Sb	Sn	Al	Te	Zr	B	Bi	Ce	Y	Pb	Mg	1 1/4 in × 1 1/4 in × 1/2 in
0.009	0.011	0.024	0.17	0.23	(0.001)	0.072	(0.02)	0.040	(0.008)	-----	-----	(0.01)	-----	1174a
.11	.35	.19	.022	.025	(.03)	.009	(.03)	.005	(.017)	-----	-----	.006	-----	1175a
-----	.049	.022	-----	-----	-----	.016	-----	-----	-----	-----	-----	-----	-----	1147
-----	.050	(.022)	-----	-----	-----	.015	-----	-----	-----	-----	-----	-----	-----	1148
-----	.062	.036	-----	-----	-----	.013	-----	-----	-----	-----	-----	-----	-----	1149
-----	.10	(.07)	-----	-----	(.01)	-----	-----	-----	-----	(0.09)	(<0.002)	-----	0.019	1140
-----	.013	(.04)	-----	-----	(.005)	-----	-----	-----	-----	(.05)	.040	-----	.044	1141
-----	.008	(.015)	-----	-----	(.09)	-----	-----	-----	-----	(.015)	.01	-----	.10	1142
-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	1138
-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	1139
-----	.17	(.004)	-----	-----	-----	.020	-----	-----	-----	-----	-----	-----	-----	1143
-----	.44	(.004)	-----	-----	-----	.020	-----	-----	-----	-----	-----	-----	-----	1144


Steelmaking Alloys

These SRM's provide standards of known chemical composition primarily for checking chemical methods of analysis for the major constituents and for selected minor elements covered by ASTM specifications. They are furnished as fine powders, sized to about 100 mesh or finer.

These SRM's are finding increasing application in x-ray and optical emission spectroscopic methods of analysis when procedures are used in which the samples to be analyzed are in the same form or can be converted to the same form; that is, to pellets, solutions, or powders.

SRM	Type	Wt/Unit (grams)	Chemical Composition (Nominal Weight Percent)						
			C	Mn	P	S	Si	Cu	Ni
57	Refined Silicon	60	0.087	0.034	0.008	0.005	96.80	0.02	0.002
58a	Ferrosilicon (75%)	IN PREP							
59a	Ferrosilicon (50%)	50	.04	.76	.016	-----	48.2	.05	.03
195	Ferrosilicon (75%)	IN PREP							
64b	Ferrochromium (HC)	100	4.30	.208	.012	.062	1.42	----	----
196	Ferrochromium(LC)	100	0.035	.28	-----	-----	.38	----	----
71	Calcium Molybdate	60	-----	-----	-----	-----	-----	-----	-----
90	Ferrophosphorus	75	-----	-----	26.2	-----	-----	-----	-----
340	Ferroniobium	100	.060	1.71	.035	-----	4.39	-----	-----

SRM	Cr	V	Mo	Ti	Al	Nb	Zr	Ca	Mg	Fe	B	N
57	0.025	----	----	0.10	0.67	-----	0.025	0.73	0.01	0.65	----	-----
58a	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----
59a	.08	-----	-----	-----	.35	-----	-----	.04	-----	50.0	0.06	-----
195	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----
64b	68.03	0.15	-----	-----	-----	-----	-----	-----	-----	-----	-----	0.033
196	70.87	.12	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----
71	-----	-----	35.3	.06	-----	-----	-----	-----	-----	1.92	-----	-----
90	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----
340	-----	-----	-----	.89	-----	57.51	Ta 3.73	-----	-----	-----	-----	-----



Nonferrous Alloys (Chip Form)

These SRM's provide materials of known composition for checking the performance of chemical methods of analysis both for production control and for customer acceptance. The aluminum-, magnesium-, and zinc-base alloys are furnished as approximately 14- to 40-mesh chips prepared by cutting thin wall castings or wrought bar stock. Certificates of Analysis provided with these standards give the composition as determined at NBS, and most give values obtained by industrial and other outside laboratories cooperating in certification of the standards.

Aluminum-Base Alloys

SRM	Type	Wt/Unit (grams)	Chemical Composition (Nominal Weight Percent)												
			Mn	Si	Cu	Ni	Cr	V	Ti	Sn	Ga	Fe	Pb	Mg	Zn
85b	Wrought	75	0.61	0.18	3.99	0.084	0.211	0.006	0.022	----	0.019	0.24	0.021	1.49	0.030
86c	Casting	75	.041	.68	7.92	.030	.029	-----	.035	----	-----	.90	.031	0.002	1.50
87a	Al-Si	75	.26	6.24	0.30	.57	.11	<.01	.18	0.05	.02	.61	.10	.37	0.16

Cobalt-Base Alloys

SRM	Type			Wt/Unit (grams)	Chemical Composition (Nominal Weight Percent)						
					Co	Ni	Cr	Mo	W	Nb	
168	Co41-Mo4-Nb3-Ta1-W4			150	41.20	20.25	20.33	3.95	3.95	2.95	

SRM	Ta	Fe	Mn	C	P	S	Si	Cu	V	Ti
168	0.95	3.43	1.50	0.37	0.008	0.005	0.80	0.035	0.03	0.06

Copper-Base Alloys

SRM	Type	Wt/Unit (grams)	Chemical Composition (Nominal Weight Percent)						
			Mn	P	S	Si	Cu	Ni	
37e	Brass, Sheet	150	----	-----	-----	----	69.61	0.53	
124d	Bronze, Ounce Metal ..	150	----	.02	.093	----	83.60	.99	
157a	Nickel Silver	135	0.174	.009	-----	----	58.61	11.82	
158a	Bronze, Silicon	150	1.11	.026	-----	3.03	90.93	0.001	
184	Bronze, Leadcd Tin ...	150	-----	.009	-----	----	88.96	.50	

SRM	Co	As	Sn	Fe	Al	Pb	Sb	Ag	Zn
37e	-----	----	1.00	0.004	----	1.00	----	----	27.85
124e	-----	0.02	4.56	.18	----	5.20	0.17	0.02	5.06
157a	0.022	----	0.021	.174	----	0.034	----	----	29.09
158a	-----	----	.96	1.23	0.46	.097	----	----	2.08
184	-----	----	6.38	0.005	----	1.44	----	----	2.69

Lead-Base Alloys

SRM	Type	(Other Forms)	Chemical Composition (Nominal Weight Percent)							
			Cu	Ni	As	Sn	Sb	Bi	Ag	Fe
53e	Bearing Metal (84Pb-10Sb-6Sn)	1132	0.054	0.003	0.057	5.84	10.26	0.052	----	<0.001
127b	Solder (40Sn-60Pb)	1131	.011	.012	.01	39.3	0.43	.06	0.01	-----

Magnesium-Base Alloys

SRM	Type	Wt/Unit (grams)	Chemical Composition (Nominal Weight Percent)							
			Mn	Si	Cu	Ni	Al	Pb	Fe	Zn
171	Alloy	100	0.45	0.0118	0.0112	0.0009	2.98	0.0033	0.0018	1.05

Nickel-Base Alloys

SRM	Type	Wt/Unit (grams)	Chemical Composition (Nominal Weight Percent)							
			C	Mn	P	S	Si	Cu	Ni	Cr
162a	Monel-type (Ni64-Cu31) . . .	150	0.079	1.60	-----	0.007	0.93	30.61	63.95	0.042
349	Ni57-Cr20	150	.08	0.43	0.002	-----	.29	0.006	57.15	19.50

SRM	V	Mo	W	Co	Ti	Al	B	Ca	Fe	Nb	Ta	Zr
162a	-----	-----	-----	0.076	0.005	0.50	-----	0.013	2.19	-----	-----	-----
349	0.081	4.04	<0.01	13.95	3.05	1.23	0.0046	-----	0.13	<0.01	<0.01	0.081

Nickel Oxide

The nickel oxide SRM's are available primarily for application in the electronics industry to the analysis of cathode grade nickel. The "Standard Method for Spectrochemical Analysis of Thermionic Nickel Alloys by the Powder-D-C Arc Technique," ASTM Designation E129 is based on calibration with these standards. The values given are for the percentage of the element in nickel oxide.

SRM	Type	Wt/Unit (grams)	Chemical Composition (Nominal Weight Percent)								
			Mn	Si	Cu	Cr	Co	Ti	Al	Fe	Mg
671	Oxide 1	25	0.13	0.047	0.20	0.025	0.31	0.024	0.009	0.39	0.030
672	Oxide 2	25	.095	.11	.018	.003	.55	.009	.004	.079	.020
673	Oxide 3	25	.0037	.006	.002	.0003	.016	.003	.001	.029	.003

Selenium Base

This SRM is intended to bridge the gap between commercial materials available in bulk and selenium available in primary or purer grades. It should prove useful to the small research laboratory, or to the individual engaged in purification, as a characterized starting material. It should be useful also as a homogeneous material in analytical procedures when a high-purity primary grade is neither necessary nor available.

SRM	Type	Wt/Unit (grams)	Chemical Composition (Nominal Parts Per Million)									
			Mn	S	Cu	Ni	Cr	V	Mo	Co	As	Sn
726	Selenium, Intermediate Purity	450	<0.3	12 ± 3	<1	<0.5	<1	N.D.	<0.3	N.D.	<2	<1

SRM	Al	B	Pb	Bi	Ag	Ca	Mg	Te	Fe	Cl	Tl	Be	Cd	In
726	<1	<1	<1	N.D.	<1	<1	<1	0.3 ± 0.1	1	<0.5	<0.5	N.D.	N.D.	N.D.

N.D. = Not detected at limits of detection of <0.5 ppm.

Tin-Base Alloys

SRM	Type	Chemical Composition (Nominal Weight Percent)								
		Pb	Sn	Sb	Bi	Cu	Fe	As	Ag	Ni
54d	Bearing metal	0.62	88.57	7.04	0.044	3.62	0.027	0.088	0.0032	0.0027

Titanium-Base Alloys

SRM	Type	Wt/Unit (grams)	Chemical Composition (Nominal Weight Percent)									
			C	Mn	Si	Cu	V	Mo	Sn	Al	Fe	N
173a	6Al-4V	100	0.025	-----	0.037	0.002	4.06	0.005	----	6.47	0.15	0.018
174	4Al-4Mn	100	-----	4.57	.015	-----	----	-----	----	4.27	.175	.012
176	5Al-2.5Sn	100	.015	0.0008	-----	.003	----	.0003	2.47	5.16	.070	.010

Zinc-Base

SRM	Type	Wt/Unit (grams)	Chemical Composition (Nominal Weight Percent)										
			Mn	Cu	Ni	Sn	Al	Cd	Fe	Pb	Ag	Mg	Ti
94c	Die Casting Alloy	150	0.014	1.01	0.006	0.006	4.13	0.002	0.018	0.006	-----	0.042	-----
728	Zinc	450	-----	0.00057	-----	(.000002)	-----	.00012	.00027	.00111	0.00011	-----	-----

Zirconium-Base

SRM	Type	Wt/Unit (grams)	Chemical Composition (Nominal Parts Per Million)										
			C	Mn	Si	Cu	Ni	Cr	Ti	Sn (Wt %)	Fe	N	U
360a	Zircaloy-2	100	136	3	51	140	554	1060	27	1.42	1441	43	0.15


Nonferrous Alloys (Solid Form)

These SRM's are designed to fill the basic needs of the nonferrous primary and secondary metals industries for analytical control, primarily with optical emission and x-ray spectroscopic methods. Both nominal chemical composition and analytical range SRM's have been prepared for many of the commercially important nonferrous alloy systems.

Copper-Base Alloys

Eight groups of copper-base alloy SRM's have been prepared to provide for analytical control by rapid instrumental methods in the copper industry. These SRM's are intended primarily for calibration of optical emission and x-ray spectroscopic equipment. These SRM's have been prepared in two forms: chill-cast form (with "C" prefix) for the producer (1 1/4 X 1 1/4 in blocks), and wrought form for the consumer (disks, 1 1/4 in in diameter and 3/4 in thick). Both forms have nearly identical chemical compositions. Each of the eight principal copper-base alloys are covered by three SRM's comprised of a "nominal-composition" together with a low- and a high-composition standard. To make the cartridge-brass SRM's more widely applicable, a number of trace elements were purposely added to and certified for these SRM's. (The low-composition cartridge-brass A, C1100, is no longer available.) The beryllium copper SRM's are representative of the nominal chemical composition of three Copper and Brass Research Association (CABRA) alloy designations. (Values in parentheses are not certified, but are given for information only.)

			Chemical Composition (Nominal Weight Percent)								
SRM		Type	Cu	Zn	Pb	Fe	Sn	Ni	Al	Sb	As
1101	C1101	Cartridge Brass B	69.50	30.30	.05	.037	.016	.013	.0006	.012	.009
1102		Cartridge Brass C	72.85	27.10	.020	.011	.006	.005	.0007	.005	.004
1103	C1103	Free-Cutting Brass A	59.23	35.7	3.73	.26	.88	.16	-----	-----	-----
1104	C1104	Free-Cutting Brass B	61.33	35.3	2.77	.088	.43	.070	-----	-----	-----
1105	C1105	Free-Cutting Brass C	63.7	34.0	2.0	.044	.21	.043	-----	-----	-----
1106	C1106	Naval Brass A	59.08	40.08	0.032	.004	.74	.025	-----	-----	-----
1107	C1107	Naval Brass B	61.21	37.34	.18	.037	1.04	.098	-----	-----	-----
1108	C1108	Naval Brass C	64.95	34.42	.063	.050	0.39	.033	-----	-----	-----
1109	C1109	Red Brass A	82.2	17.4	.075	.053	.10	.10	-----	-----	-----
1110	C1110	Red Brass B	84.59	15.20	.033	.033	.051	.053	-----	-----	-----
1111	C1111	Red Brass C	87.14	12.81	.013	.010	.019	.022	-----	-----	-----
1112	C1112	Gilding Metal A	93.38	6.30	.057	.070	.12	.100	-----	-----	-----
1113	C1113	Gilding Metal B	95.03	4.80	.026	.043	.064	.057	-----	-----	-----
1114	C1114	Gilding Metal C	96.45	3.47	.012	.017	.027	.021	-----	-----	-----
1115	C1115	Commercial Bronze A . . .	87.96	11.73	.013	.13	.10	.074	-----	-----	-----
1116	C1116	Commercial Bronze B . . .	90.37	9.44	.042	.046	.044	.048	-----	-----	-----
1117	C1117	Commercial Bronze C . . .	93.01	6.87	.069	.014	.021	.020	-----	-----	-----
1118	C1118	Aluminum Brass A	75.1	21.9	.025	.065	-----	-----	2.80	.010	.007
1119	C1119	Aluminum Brass B77.1	20.5	.050	.030	-----	-----	2.14	.050	.040
1120	C1120	Aluminum Brass C80.1	18.1	.105	.015	-----	-----	1.46	.100	.090
1121	C1121	Beryllium Copper CA-17297.49	(0.01)	(.002)	.085	.01	.012	0.07	-----	-----
1122	C1122	Beryllium Copper CA-17097.45	(.01)	(.003)	.16	(.01)	(.01)	.17	-----	-----
1123	C1123	Beryllium Copper CA-17597.10	(.01	(.001)	.04	(.01)	(.01)	.02	-----	-----



(Continued)

Chemical Composition (Nominal Weight Percent)

Be	Bi	Cd	Mn	P	Si	Ag	Te	Co	Cr	SRM	
.00055	.0004	.0055	.0055	.0020	(.005)	.003	.0015	-----	-----	1101	C1101
.00003	.0005	.0045	.0045	.0048	(.002)	.0010	.0003	-----	-----	1102	
-----	-----	-----	-----	.003	-----	-----	-----	-----	-----	1103	C1103
-----	-----	-----	-----	.005	-----	-----	-----	-----	-----	1104	C1104
-----	-----	-----	-----	.003	-----	-----	-----	-----	-----	1105	C1105
-----	-----	-----	.005	-----	-----	-----	-----	-----	-----	1106	C1106
-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	1107	C1107
-----	-----	-----	.025	-----	-----	-----	-----	-----	-----	1108	C1108
-----	-----	-----	-----	.006	-----	-----	-----	-----	-----	1109	C1109
-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	1110	C1110
-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	1111	C1111
-----	-----	-----	-----	.009	-----	-----	-----	-----	-----	1112	C1112
-----	-----	-----	-----	.008	-----	-----	-----	-----	-----	1113	C1113
-----	-----	-----	-----	.009	-----	-----	-----	-----	-----	1114	C1114
-----	-----	-----	-----	.005	-----	-----	-----	-----	-----	1115	C1115
-----	-----	-----	-----	.008	-----	-----	-----	-----	-----	1116	C1116
-----	-----	-----	-----	.002	-----	-----	-----	-----	-----	1117	C1117
-----	-----	-----	-----	.13	.0021	-----	-----	-----	-----	1118	C1118
-----	-----	-----	-----	.070	.0015	-----	-----	-----	-----	1119	C1119
-----	-----	-----	-----	.018	.0011	-----	-----	-----	-----	1120	C1120
1.90	-----	-----	(.004)	(.005)	.11	(.005)	-----	0.295	(0.002)	1121	C1121
1.75	-----	-----	(.004)	(.004)	.17	(.005)	-----	.220	(.002)	1122	C1122
0.46	-----	-----	(.002)	(.002)	.03	(.009)	-----	2.35	(.002)	1123	C1123

Lead-Base Alloy

These SRM's are designed primarily for the calibration of optical emission and x-ray spectro-metric methods of analysis. SRM 1131, Solder, and SRM 1132, Bearing Metal, are available in chip form as: 53e (1132) and 127b (1131), see Lead-base Alloys, page 32.

SRM	Type	Other Forms	Size	Chemical Composition (Nominal Weight Percent)							
				Cu	Ni	As	Sn	Sb	Bi	Ag	Fe
1131	Solder Pb60-Sn40	127b	31.4 mm D × 19 mm thick	0.011	0.012	0.01	39.3	0.43	0.06	0.01	-----
1132	Bearing Metal	53e	31.4 mm D × 19 mm thick	.054	.003	.057	5.84	10.2	.052	---	<0.001

Nickel-Base Alloys

These SRM's are designed primarily for calibration in optical emission and x-ray spectroscopic methods of analysis. They are issued in disk form.

SRM			Chemical Composition (Nominal Weight Percent)										
	Type	Unit Size	C	Mn	P	S	Si	Cu	Ni	Cr	Mo	Co	Fe
1159	Ni48, balance Fe	1 1/4 in D x 3/4 in thick	0.007	0.305	0.003	0.003	0.32	0.038	48.2	0.06	0.010	0.022	51.0
1160	Ni80, Mo4, balance Fe . . .	1 1/4 in D x 3/4 in thick	.019	.550	.003	.001	.37	.021	80.3	.05	4.35	.054	14.3

Tin-Base Alloys

These tin metal SRM's have been prepared primarily for the tin-plate industry; they are useful for the calibration of optical emission spectroscopic equipment by the "point-to-point" technique. They are furnished in the form of rods, 1/4 in in diameter and 4 in long.

SRM	Type	Cu	Ni	Co	As	Pb	Chemical Composition (Nominal Weight Percent)				
							Sb	Bi	Ag	Zn	Cd
431	Tin A	0.19	0.038	0.021	0.16	0.19	0.19	0.020	0.015	0.041	0.020
432	Tin B097	.020	.011	.075	.094	.095	.0098	.0095	.020	.0095
433	Tin C055	.0095	.0045	.047	.055	.050	.0052	.0055	.0095	.0053
434	Tin D019	.0044	.0020	.019	.022	.019	.0020	.0018	.0046	.0020
435	Tin E0077	.0024	.0011	.0090	.015	.010	.0011	.0010	.0020	.0011

Titanium-Base Alloys

These SRM's are issued primarily for the aerospace industries for analytical control and equipment calibration purposes. These SRM's are in disk form intended as calibration materials for optical emission and x-ray spectroscopic methods of analysis of similar materials.

SRM	Type	Unit Size	Chemical Composition (Nominal Weight Percent)					
			Mn	Cr	Fe	Mo	Al	V
641	8Mn (A)	1 1/4 in D × 3/4 in Disks	6.68	----	----	----	----	----
642	8Mn (B)	1 1/4 in D × 3/4 in Disks	9.08	----	----	----	----	----
643	8Mn (C)	1 1/4 in D × 3/4 in Disks	11.68	----	----	----	----	----
644	2Cr-2Fe-2Mo (A) . . .	1 1/4 in D × 3/4 in Disks	-----	1.03	1.36	3.61	----	----
645	2Cr-2Fe-2Mo (B) . . .	1 1/4 in D × 3/4 in Disks	-----	1.96	2.07	2.38	----	----
646	2Cr-2Fe-2Mo (C) . . .	1 1/4 in D × 3/4 in Disks	-----	3.43	2.14	1.11	----	----
654a	6Al-4V (B)	1 1/4 in D × 1/4 in Disks	(<0.1)	(0.20)	(0.20)	(<0.05)	6.3 ₄	3.9 ₅

Titanium-Base Alloys – Oxygen and Hydrogen Only

SRM's intended for determination of hydrogen and oxygen in titanium-base alloys are available in sheet and rod form. These were designed primarily for calibration of vacuum fusion or inert gas fusion equipment.

A group of iron-base alloys certified for oxygen also are available. See Oxygen Standards, page 27.

SRM	Type	Unit Size	Wt/Unit (grams)	Oxygen (ppm)	Hydrogen (Wt %)
352	Unalloyed titanium for hydrogen	1/4 in square × 0.05 in thick	20	----	0.0032
353	Unalloyed titanium for hydrogen	1/4 in square × 0.05 in thick	20	----	.0098
354	Unalloyed titanium for hydrogen	1/4 in square × 0.05 in thick	20	----	.0215
355	Unalloyed	Rod-1/2 in D × 2 in long	-----	3031	-----
356	Alloy, 6Al-4V	Rod-.425 in D × 1 3/4 in long	-----	1332	-----

Zinc-Base Alloys

Zinc-base alloy SRM's are available ranging from very high-purity zinc to commercial materials such as spelter and die-casting alloy compositions. They are supplied as bar segments (disks) intended for calibrating and checking optical emission and x-ray spectroscopic techniques. The certificate of analysis supplied with each gives the chemical composition determined at NBS and values determined by other laboratories that have cooperated in the certification of the SRM's. For high-purity Zinc, see High-Purity Metals, page 40.

Die Casting Alloys and Spelter

Zinc-base die casting alloys and a spelter SRM are available. They were designed for calibration of optical emission spectroscopic techniques primarily for analysis of such alloys as ASTM Designations AG 40A and AC 41A. The SRM's were prepared by a continuous chill-casting process into square bars which then were cut into segments. The certified portion of each segment is that part included between 3/16 inch and 11/16 inch from each side. The center core, 3/16 inch square, and the outer rim, 3/16 inch from the outer surface, are parts which may differ in chemical composition for some elements from the certified portion, and should not be used.

(Values in parentheses are not certified, but are given for information only.)

SRM	Type	Unit Size	Chemical Composition (Nominal Weight Percent)					
			Cu	Al	Mg	Fe	Pb	Cd
625	Zinc-base A-ASTM AG 40A ..	1 3/4 in square × 3/4 in thick	0.034	3.06	0.070	0.036	0.0014	0.0007
626	Zinc-base B-ASTM AG 40A ..	1 3/4 in square × 3/4 in thick	.056	3.56	.020	.103	.0022	.0016
627	Zinc-base C-ASTM AG 40A ..	1 3/4 in square × 3/4 in thick	.132	3.88	.030	.023	.0082	.0051
628	Zinc-base D-ASTM AC 41A ..	1 3/4 in square × 3/4 in thick	.611	4.59	.0094	.066	.0045	.0040
629	Zinc-base E-ASTM AC 41A ..	1 3/4 in square × 3/4 in thick	1.50	5.15	.094	.017	.0135	.0155
630	Zinc-base F-ASTM AC 41A ..	1 3/4 in square × 3/4 in thick	0.976	4.30	.030	.023	.0083	.0048
631	Zinc spelter (modified)	1 3/4 in square × 3/4 in thick	.0013	0.50	(<0.001)	.005	(.001)	.0002

SRM	Chemical Composition (Nominal Weight Percent)									
	Sn	Cr	Mn	Ni	Si	In	Ga	Ca	Ag	Ge
625	0.0006	0.0128	0.031	0.0184	0.017	-----	-----	-----	-----	-----
626	.0012	.0395	.048	.047	.042	-----	-----	-----	-----	-----
627	.0042	.0038	.014	.0029	.021	-----	-----	-----	-----	-----
628	.0017	.0087	.0091	.030	.009	-----	-----	-----	-----	-----
629	.012	.0008	.0017	.0075	.078	-----	-----	-----	-----	-----
630	.0040	.0031	.0106	.0027	.022	-----	-----	-----	-----	-----
631	.0001	.0001	.0015	(<0.0005)	<0.002	(0.0023)	(0.002)	<0.001	(<0.0005)	(0.0002)

Zirconium-Base Alloys

A zirconium metal SRM certified at the parts-per-million level, is available in the form of a wrought disk for checking and calibrating of optical emission and x-ray spectroscopic instruments used in the analytical control of trace level constituents of zirconium for use in nuclear power applications.

(Values in parentheses are not certified, but are given for information only.)

SRM			Chemical Composition (Nominal Parts Per Million)										
	Type	Unit Size	Mn	Si	Cu	Ni	Cr	W	Ti	Al	B	U	Fe
1210	Zirconium A	1 1/4 in D×3/4 in thick	(5)	(30)	10	8	95	(4)	26	(60)	(≤0.25)	1.8	2500

Gases in Metals

Values in parentheses are not certified but are given for
additional information on composition.

Certified for Nitrogen

Steels, (Chip Form)

SRM	Type	Nitrogen (Wt %)
10g	Bessemer	0.015
12h	Basic Open Hearth, 0.4C006
32e	Ni-Cr (SAE 3140)009
33d	Ni-Mo (SAE 4820)	(.011)
50c	W18-Cr4-V1012
51b	Electric Furnace, 1.2C011
65d	Basic Electric Furnace, 0.3C . .	.013
72f	Cr-Mo (SAE X4130)009
73c	Stainless (Cr13) (SAE 420) . .	.037
100b	Manganese (SAE T1340)004
133a	Stainless (Cr13-Mo0.3-S0.3) . .	.032
139a	Cr-Ni-Mo (AISI 8640)008
153a	Co8-Mo9-W2-Cr4-V2 (Tool) . .	.024
160b	Stainless (AISI 316)039
343	Stainless (SAE 431)074
346	Valve (Cr22-Ni4-Mn9)441

Cast Irons (Chip Form)

SRM	Type	Nitrogen (Wt %)
5L	Cast iron	0.006
6g	Cast iron006
7g	Cast iron (high phosphorus) . .	.004
55e	Ingot iron004
107b	Cast iron (Ni-Cr-Mo)	(.008)

Titanium Base (Chip Form)


SRM	Type	Nitrogen (Wt %)
173a	6Al-4V	0.018
174	4Al-4Mn012
176	5Al-2.5Sn010

Steel, (Granular Form)

SRM	Type	Nitrogen (Wt %)
163	Low alloy, Cr1	0.007

Zirconium Base Alloys, (Chip Form)

SRM	Type	Nitrogen (Wt %)
360a	Zircaloy-2	0.0043



Certified for Oxygen and Nitrogen

Steels (solid form)

SRM	Type	Oxygen (ppm)	Nitrogen (ppm)
1090	Ingot iron	491	(60)
1091	Stainless steel (AISI 431)	131	(945)
1092	Vacuum melted steel	28	(3.6)
1093	Valve steel	60	(4807)
1094	Maraging steel	4.5	(71)
1095	AISI 4340 Steel	9	(37)
1096	AISI 94B17 (Mod) Steel	10.7	40.4
1097	Cr-V (Mod) Steel	6.6	(41)
1098	High Carbon (Mod) Steel	(17)	30
1099	Electrolytic Iron	61	(13)

Certified for Hydrogen or Oxygen

Titanium Base (solid form)

SRM	Type	Hydrogen (Wt %)	Oxygen (ppm)
352	Unalloyed	0.0032	----
353	Unalloyed0098	----
354	Unalloyed0215	----
355	Unalloyed	-----	3031
356	6Al-4V Alloy	-----	1332

High-Purity Metals

Very high-purity metal SRM's are being made available to fill the needs of analysts determining impurity elements in high-purity metal materials. They are intended to serve as bench marks in calibration of methods and equipment; also, they are expected to be valuable in the development of new or improved methods and techniques for extending the sensitivity of detection in the determination of trace constituents in various materials by chemical, optical emission, solid mass spectroscopy, activation, and resistivity methods of analysis.

The certificate of analysis supplied with each high-purity SRM gives the state-of-the-art information on chemical composition in the cooperating laboratories for the various trace determinations reported.

High-purity gold is available in both wire and rod form. The wire form (W), is intended for applications such as in spark source mass spectroscopic techniques. The low levels of impurities make it important for evaluating instrument and system blanks. The rod form (R), is intended for application in other methods of characterization.


Platinum is available in wire form as a high-purity material and as a doped composition material.

Zinc is available in a high-purity and in a less pure version. Both were prepared from the same starting material. The high-purity material is the result of further purification by vacuum distillation, zone refining, and degasification. The zinc is supplied in the form of semicircular bar segments.

SRM	Type	Unit Size	Chemical Compositions (Nominal Parts Per Million by Weight)				
			Cu	Ni	Sn	Pb	Zr
685W*	High-Purity Gold (Wire)	1.4mm D × 102mm long	0.1	----	----	----	----
685R*	High-Purity Gold (Rod)	5.9mm D × 25mm long	.1	----	----	----	----
680L1	High-Purity Platinum (Wire)	0.51mm D × 102mm long	.1	<1	----	<1	<0.1
680L2	High-Purity Platinum (Wire)	0.51mm D × 1.0m long					
681L1	Doped-Platinum (Wire)	0.51mm D × 102mm long	5.1	0.5	----	12	11
681L2	Doped-Platinum (Wire)	0.51mm D × 1.0m long					
682*	High-Purity Zinc	Semi circular segments 57mm D × 19mm long	0.042	----	(0.02)	----	----
683*	Zinc Metal	Semi circular segments 57mm D × 19mm long	5.9	----	(.02)	11.1	----

SRM	Chemical Compositions (Nominal Parts Per Million by Weight)										
	Ag	Mg	In	Fe	O	Pd	Au	Rh	Ir	Cd	Tl
685W*	[0.1]	---	0.007	0.3	[2]	---	---	---	-----	----	----
685R*	[.1]	---	.007	.2	[<2]	---	---	---	-----	----	----
680L1	.1	<1	-----	.7	4	0.2	<1	<0.2	<0.01	----	----
680L2											
681L1	2.0	12	-----	5	7	6	9	9	11	----	----
681L2											
682*	(0.02)	---	-----	(0.1)	----	---	---	---	-----	(0.1)	----
683*	1.3	---	-----	2.2	----	---	---	---	-----	1.1	(0.2)

*Certificate gives upper limits for other elements found to be present.



Microprobe Standards

These SRM's provide a highly homogeneous material at about the micrometer of spatial resolution. They are intended primarily for use in calibration of quantitative electron microprobe analytical techniques.

Fe-Cr-Ni Alloy

The Fe-Cr-Ni alloy, SRM 479, is a wafer (4.6 mm in diameter and 1 mm thick) and is characterized for chemical homogeneity of iron, chromium, and nickel at the micrometer level of spatial resolution. It is satisfactory for use as a homogeneous material for electron probe microanalysis.

Tungsten - 20% Molybdenum

The tungsten-20% molybdenum alloy, SRM 480, consists of a core of tungsten-20% molybdenum wire embedded in pure molybdenum onto which pure tungsten has been deposited by electroplating to provide a composite. The Certificate supplied with this SRM gives the values for tungsten and molybdenum as determined by analysis and a summary of homogeneity testing results from approximately 1500 determinations for tungsten and molybdenum by electron probe microanalysis. Additional details on homogeneity characterization are given in NBS Spec. Publ. 260-16. (See inside back cover for ordering instructions.)

Gold-Silver

Six color-coded wires comprise this set, SRM 481. The wires consist of a high-purity gold and a high-purity silver wire and four wires with nominal chemical composition differences in steps of 20%. The Certificate of Analysis supplied with each SRM gives the results of a cooperative program of analysis and a summary of the extensive homogeneity testing performed in certifying the wires.

Although designed for quantitative elemental microprobe analysis, the wires should be equally useful for other microtechniques.

Gold-Copper

This set of color-coded wires, SRM 482, is similar to the gold-silver set. In both sets special precautions were taken to achieve homogeneity on a microscopic scale.

Iron - 3% Silicon

The iron-3% silicon microprobe, SRM 483, is a platelet (3 mm X 3 mm X 0.28 mm), and is characterized for chemical homogeneity of iron and silicon at the micrometer level of spatial resolution. It is satisfactory for use as a homogeneous material for electron probe microanalysis.

SRM	Type	Unit Size	Chemical Composition (Nominal Weight Percent)								
			Au	Cu	Ag	W	Mo	Si	Fe (by difference)	Cr	Ni
479	Fe-Cr-Ni Alloy	Wafer-4.6 mm D X 1 mm thick -----							71.0	18.3	10.7
480	Tungsten-20 Mo Alloy . .	Wafer-1 mm D X 1 mm thick -----				78.5	21.5				
481	Au100A	Wire-0.5mm D X 50 mm long	100.00								
	Au80-Ag20B	Wire-0.5mm D X 50 mm long	80.05		19.96						
	Au60-Ag40C	Wire-0.5mm D X 50 mm long	60.05		39.92						
	Au40-Ag60D	Wire-0.5mm D X 50 mm long	40.00		59.90						
	Au20-Ag80E	Wire-0.5mm D X 50 mm long	22.43		77.58						
	Ag100F	Wire-0.5mm D X 50 mm long			100.00						
482	Au100A	Wire-0.5mm D X 50 mm long	100.00								
	Au80-Cu20B	Wire-0.5mm D X 50 mm long	80.15	19.83							
	Au60-Cu40C	Wire-0.5mm D X 50 mm long	60.36	39.64							
	Au40-Cu60D	Wire-0.5mm D X 50 mm long	40.10	59.92							
	Au20-Cu80E	Wire-0.5mm D X 50 mm long	20.12	79.85							
	Cu100F	Wire-0.5mm D X 50 mm long		100.00							
483	Iron-3% Silicon	Wafer-3mm X 3mm X 0.28m -----						3.22	96.7-96.8		

Primary, Working, and Secondary Standard Chemicals

These SRM's are high-purity chemicals defined as primary, working, and secondary standards in accordance with recommendations of the Analytical Chemistry Section of the International Union of Pure and Applied Chemistry [Ref. Analyst 90, 251 (1965)]. These definitions are as follows:

Primary Standard:

a commercially available substance of purity 100 ± 0.02 percent (Purity 99.98+ percent).

Working Standard:

a commercially available substance of purity 100 ± 0.05 percent (Purity 99.95+ percent).

Secondary Standard:

a substance of lower purity which can be standardized against a primary grade standard.

SRM	Type	Wt/Unit (grams)	Certified Use	Purity
				Stoichiometric
17	Sucrose	60	Polarimetric Value	a
40h	Sodium Oxalate	60	Reductometric Value	99.95
41a	Dextrose (D-glucose)	70	Reductometric Value	b
83c	Arsenic trioxide	75	Reductometric Value	99.99
84h	Acid potassium phthalate	60	Acidimetric Value	99.99
136c	Potassium dichromate	60	Oxidimetric Value	99.98
350	Benzoic Acid	30	Acidimetric Value	99.98
723	Tris(hydroxymethyl)aminomethane	50	Basimetric Value	99.97
944	Plutonium Sulfate Tetrahydrate	0.5	Assay	100
949d	Plutonium Metal	0.5	Assay	99.99
950a	Uranium oxide (U_3O_8)	25	Uranium Oxide Standard Value	99.94
951	Boric Acid	100	Acidimetric and Boron Isotopic Value	100.00
960	Uranium Metal	26	Assay	99.975
984	Rubidium Chloride	1	Assay and Isotopic	99.90
987	Strontium Carbonate	1	Assay and Isotopic	99.98
988	Strontium-84 Spike	0.010	Assay and Isotopic	99.9
999	Potassium Chloride	60	Assay Standard for Potassium Chloride	99.98 99.99

^aSucrose - Moisture <0.01 percent, Reducing Substances <0.02 percent, Ash 0.003 percent.

^bDextrose - Moisture <0.2 percent, Ash <0.01 percent.

Microchemical Standards

These SRM's are furnished as fine crystals of suitable homogeneity for use as standards for conventional microchemical methods of analysis employing samples of approximately 5 mg. See also Microprobe Standards, page 41 .

SRM	Type	Wt/Unit (grams)	Elements Certified
140b	Benzoic acid	2	C,H
141b	Acetanilide	2	N,C,H
142	Anisic acid	2	Methoxyl(CH_3O-)
143b	Cystine	2	S,C,H,N
147	Triphenyl phosphate	2	P
148	Nicotinic acid	2	N,C,H
2141	Urea	2	N
2142	o-Bromobenzoic Acid	2	Br
2143	p-Fluorobenzoic Acid	2	F(IN PREP)
2144	m-Chlorobenzoic Acid	2	Cl(IN PREP)

Clinical Laboratory Standards

These SRM's are intended for use in the calibration of apparatus and checking methods of analysis used in clinical and pathological laboratories, and to assist manufacturers of clinical products in meeting the chemical and physical specifications required for clinical chemicals. (For details on SRM's 930 and 931, see Spectrophotometric Filters, page 65.)

SRM	Type	Purity %	Wt/Unit
911	Cholesterol	99.4	0.5g
912	Urea	99.7	25g
913	Uric Acid	99.7	10g
914	Creatinine	99.8	10g
915	Calcium Carbonate*	99.9	20g
916	Bilirubin	99.0	100mg
917	D-Glucose	99.9	25g
918	Potassium Chloride	99.9	30g
919	Sodium Chloride	99.9	30g
920	D-Mannitol	99.8	50g
921	Cortisol	IN PREP	
922	Tris(hydroxymethyl)aminomethane	99.9	25g
923	Tris(hydroxymethyl)aminomethane HCl ...	99.7	35g
924	Lithium Carbonate	100.5	30g
925	VMA (4-hydroxy-3-methoxymandelic acid)	IN PREP	
930a	Glass Filters for Spectrophotometry		Set of 3
931	Liquid Filters for Spectrophotometry		3 Sets of 4

*SRM 915, Calcium Carbonate, was used to develop the first referee method of analysis in clinical chemistry. This work is described in NBS Special Publication 260-36, A Referee Method for the Determination of Calcium in Serum. (See inside of back cover for ordering instructions.)

Biological Standards

These SRM's are intended for use in the calibration of apparatus and methods used in the analysis of biological materials for major, minor, and for trace elements.

(Values in parentheses are not certified, but are given for information only.)

Content in $\mu\text{g/g}$ (or where noted, wt %)

SRM	Type	Wt/Unit (grams)	As	Bi	B	Br	Cd	Ca	Cl
1571	Orchard Leaves	75	14	(0.1)	33	(10)	0.11	2.09%	(700)
1573	Tomato Leaves	IN PREP							
1577	Liver, Bovine	50	(0.055)	---	--	---	.27	(123)	(2600)

SRM	Co	Cu	Cr	F	Fe	Pb	Li	Mg	Mn	Hg	Mo	Ni
1571	(0.2)	12	(2.3)	(4)	300	45	(14)	0.62%	91	0.155	-----	1.3
1573												
1577	(.18)	193	-----	---	270	0.34	----	(605)	10.3	.016	(3.2)	---

SRM	N	P	K	Rb	Se	Ag	Na	Sr	S	Tl	U	Zn
1571	2.76%	0.21%	1.47%	12	0.08	-----	82	(37)	(2300)	-----	0.029	25
1573												
1577	10.6%	-----	0.97%	18.3	1.1	(0.06)	0.243%	(0.14)	-----	(0.05)	(.0008)	130

Metallo-organic Compounds

These SRM's are intended for the preparation of solutions in oils of known and reproducible concentrations of metals. Because "matrix" effects occur, it is desirable to prepare the standard solutions in oil identical or similar to the oil being studied. Possession of an adequate collection of these metallo-organic SRM's permits the preparation of any desired blend of known concentrations of metal in the appropriate lubricating oil. They are used primarily for the calibration of spectrochemical equipment used in the determination of metals in lubricating oil. This technique is used extensively in the defense program, the transportation industry, and other industries where the consequences of failure of a moving metal part may range from inconvenient to catastrophic. Details of the selection, preparation, and analysis of the compounds can be found in NBS Monograph 54. (See inside back cover for ordering instructions.)

The Certificate supplied with each SRM gives the percentage of the element of interest and directions for preparing a solution of known concentration in lubricating oil.

SRM	Constituent Certified		Wt/Unit (grams)	Type
	Element	(wt. percent)		
1075a	Al	8.1	5	Aluminum 2-ethylhexanoate
1051b	Ba	28.7	5	Barium cyclohexanebutyrate
1063a	B	2.4	5	Menthyl borate
1053a	Cd	24.8	5	Cadmium cyclohexanebutyrate
1074a	Ca	12.5	5	Calcium 2-ethylhexanoate
1078b	Cr	9.6	5	Tris(1-phenyl-1,3-butanediono)chromium(III)
1055b	Co	14.8	5	Cobalt cyclohexanebutyrate
1080	Cu	16.5	5	Bis(1-phenyl-1,3-butanediono)copper(II)
1079b	Fe	10.3	5	Tris(1-phenyl-1,3-butanediono)iron(III)
1059b	Pb	36.7	5	Lead cyclohexanebutyrate
1060a	Li	4.1	5	Lithium cyclohexanebutyrate
1061c	Mg	6.5	5	Magnesium cyclohexanebutyrate
1062a	Mn	13.8	5	Manganous cyclohexanebutyrate
1064	Hg	36.2	5	Mercuric cyclohexanebutyrate
1065b	Ni	13.9	5	Nickel cyclohexanebutyrate
1071a	P	9.5	5	Triphenyl phosphate
1066a	Si	14.1	5	Octaphenylcyclotetrasiloxane
1076	K	10.1	5	Potassium erucate
1077a	Ag	42.6	5	Silver 2-ethylhexanoate
1069b	Na	12.0	5	Sodium cyclohexanebutyrate
1070a	Sr	20.7	5	Strontium cyclohexanebutyrate
1057b	Sn	23.0	5	Dibutyltin bis(2-ethylhexanoate)
1052b	V	13.0	5	Bis(1-phenyl-1,3-butanediono)oxovanadium(IV)
1073b	Zn	16.7	5	Zinc cyclohexanebutyrate



ENVIRONMENTAL STANDARDS

Analyzed Gases

These SRM's are intended for the calibration of apparatus used for the measurement of various components in gas mixtures, and in some cases for particular atmospheric pollutants. Each SRM is accurately certified and is primarily intended to monitor and correct for long-term drifts in instruments used.

SRM	Type	Vol/Unit (liters at STP)	Certified Constituents
1601	Carbon dioxide in Nitrogen	68	CO ₂ , 308 ppm
1602	Carbon dioxide in Nitrogen	68	CO ₂ , 346 ppm
1603	Carbon dioxide in Nitrogen	68	CO ₂ , 384 ppm
1604a	Oxygen in Nitrogen	68	O ₂ , 1.5 ppm
1605	Oxygen in Nitrogen	68	O ₂ , 10 ppm
1606	Oxygen in Nitrogen	68	O ₂ , 112 ppm
1607	Oxygen in Nitrogen	68	O ₂ , 212 ppm
1608	Oxygen in Nitrogen	68	O ₂ , 978 ppm
1609	Oxygen in Nitrogen	68	O ₂ , 20.95 mole percent
1610	Hydrocarbon in Air	68	Methane, 0.103 mole percent
1611	Hydrocarbon in Air	68	Methane, 0.0107 mole percent
1613	Hydrocarbon in Air	68	Methane, 0.000102 mole percent

Note: SRM's 1665 through 1669, Propane in Air, and SRM's 1673 through 1675, Carbon Dioxide in Nitrogen, are available for calibrating equipment used to monitor automotive emission gases.

Permeation Tubes

These SRM's are intended for calibrating air-pollution monitoring apparatus, and may be used to verify air-pollution analytical methods and procedures. Each tube is individually calibrated and its permeation rate is certified over the temperature range of 20 to 30 °C.

The following table is provided for information only. The values given in the table do not represent certified values for any individual SRM. The concentrations of SO₂ in ppm are based on an approximate permeation rate of 0.28 µg per cm per minute at 25 °C, for flow rates of 1, 5, and 10 liters per minute.

SRM	Type	Tube Length (cm)	Permeation Rate (µg per min)	Typical Concentrations (ppm) Flow Rates (liters per minute)		
				1	5	10
1625	Sulfur Dioxide Permeation Tube	10	2.8	1.07	0.214	0.107
1626	Sulfur Dioxide Permeation Tube	5	1.4	0.535	.107	.0535
1627	Sulfur Dioxide Permeation Tube	2	0.56	.214	.0428	.0214

Analyzed Liquids

These SRM's are intended for use in the analysis of liquids for elements that, when liberated, could become environmental pollutants.

SRM	Type	Element Certified	Wt %	Vol/Unit (ml)
1621	Sulfur in Residual Fuel Oil	S	1.05	100
1622	Sulfur in Residual Fuel Oil	S	2.14	100
1623	Sulfur in Residual Fuel Oil	S	0.268	100
1624	Sulfur in Distillate Fuel Oil	S	.211	100

Analyzed Solids

These SRM's are intended for use in the analysis of materials for elements of interest in health or environmental problems. (See also Clinical SRM's, page 43.)

SRM	Type	Element Certified	Content	Wt/Unit
1579	Powdered Lead Based Paint	Pb	11.87 wt%	35
1630	Mercury in Coal	Hg	0.13 ppm	50
1631	Sulfur in Coal	S		set(3)

Hydrocarbon Blends

Eight standard hydrocarbon blends are available for calibration of mass spectrometers and gas chromatographic procedures used in the analysis of gasolines, naphthas and blending stocks. The even numbered SRM's, 592, 594, 596, and 598, are representative of typical virgin naphthas and the odd numbered SRM's 593, 595, 597, and 599, are representative of typical catalytically cracked naphthas in the C₇ and C₈ paraffin and cycloparaffin series.

Each SRM is supplied in a unit of ten sealed ampoules. Each ampoule contains 0.03 ml of the blend. Each ampoule is intended to provide material for only one calibration analysis so that possible fractionation of components will be avoided.

For individual components present in the mixtures in the amount of 10% or less (by volume), the limits of error in composition are not greater than ± 0.01 percent and for components present in more than 10 percent, the limits of error are not greater than ± 0.10 percent.

SRM	592	593	594	595	596	597	598	599
Blend No.	1	2	3	4	5	6	7	8
Unit (Ampoules)	10	10	10	10	10	10	10	10
Hydrocarbon	Volume Percent (Nominal)							
n-Heptane	45	17	----	----	----	----	----	----
2-Methylhexane	23	25	----	----	----	----	----	----
3-Methylhexane	16	30	----	----	----	----	----	----
2,2-Dimethylpentane	4	----	----	----	----	----	----	----
2,3-Dimethylpentane	6	20	----	----	----	----	----	----
2,4-Dimethylpentane	5	8	----	----	----	----	----	----
3,3-Dimethylpentane	1	----	----	----	----	----	----	----
n-Octane	----	----	39	12	----	----	----	----
2-Methylheptane	----	----	19	25	----	----	----	----
3-Methylheptane	----	----	16	23	----	----	----	----
4-Methylheptane	----	----	8	8	----	----	----	----
3-Ethylhexane	----	----	3	3	----	----	----	----
2,3-Dimethylhexane	----	----	4	9	----	----	----	----
2,4-Dimethylhexane	----	----	5	5	----	----	----	----
2,5-Dimethylhexane	----	----	6	9	----	----	----	----
3,4-Dimethylhexane	----	----	----	6	----	----	----	----
Methylcyclohexane	----	----	----	----	57	32	----	----
Ethylcyclopentane	----	----	----	----	9	14	----	----
1,1-Dimethylcyclopentane	----	----	----	----	4	3	----	----
1,trans-2-Dimethylcyclopentane	----	----	----	----	14	30	----	----
1,trans-3-Dimethylcyclopentane	----	----	----	----	16	21	----	----
Ethylcyclohexane	----	----	----	----	----	----	20	17
1,trans-2-Dimethylcyclohexane	----	----	----	----	----	----	18	7
1,cis-3-Dimethylcyclohexane	----	----	----	----	----	----	25	19
1,trans-4-Dimethylcyclohexane	----	----	----	----	----	----	11	14
1-Methyl-cis-2-ethylcyclopentane	----	----	----	----	----	----	7	20
1,1,3-Trimethylcyclopentane	----	----	----	----	----	----	5	4
1,trans-2-cis-3-Trimethylcyclopentane	----	----	----	----	----	----	9	6
1,trans-2-cis-4-Trimethylcyclopentane	----	----	----	----	----	----	5	13

Ores

These SRM's are intended for use in checking the accuracy of assay methods. They are certified for their content of elements of economic interest, and occasionally, have additional data given for information only. These SRM's are supplied in the form of fine powders, usually passing a 100-mesh or finer sieve.

Chemical Composition (Nominal Weight Percent)

SRM	Type	Wt/Unit (Gram)	CaF	Fe	Mn	Li ₂ O	SiO ₂	P ₂ O ₅	P	Available Oxygen
25c	Manganese	100	-----	-----	57.85	---	2.36	0.22	-----	16.7
27e	Iron (Sibley)	100	-----	66.58	-----	---	3.65	-----	0.042	---
79a	Fluorspar	120	97.39	-----	-----	---	-----	-----	-----	---
113a	Zinc, concentrate	IN PREP	-----	-----	-----	-----	-----	-----	-----	---
180	Fluorspar, high-grade	120	98.8	-----	-----	---	-----	-----	-----	---
181	Lithium (Spodumene)	45	-----	-----	-----	6.4	-----	-----	-----	---
182	Lithium (Petalite)	45	-----	-----	-----	4.3	-----	-----	-----	---
183	Lithium (Lepidolite)	45	-----	-----	-----	4.1	-----	-----	-----	---

Chemical Composition
(Nominal Weight Percent as the Oxide)

SRM	Type	Wt/Unit (gram)	Al ₂ O ₃	CaO	P ₂ O ₅	SiO ₂	Fe ₂ O ₃	F	CO ₂	TiO ₂	Na ₂ O	MgO
69a	Bauxite	50	55.0	0.29	0.08	6.01	5.8	----	----	2.78	<0.01	0.02
120b	Phosphate Rock (Florida)	90	1.06	49.40	34.57	4.68	1.10	3.84	2.79	0.15	< .35	.28
	ZrO ₂	K ₂ O	K ₂ O	Cr ₂ O ₃	SO ₃	MnO	V ₂ O ₅	BaO	CdO	Loss on Ignition		
69a	0.18	<0.01	-----	0.05	0.04	0.02	0.03	0.01	-----	29.55		
120b	----	< .12	0.090	-----	-----	.28	----	----	0.002	-----		

Concentrations in weight percent
(except as noted)

SRM	TYPE	Wt/Unit (grams)	Total Cu	Re	Mo	Au	Ag
330	Copper, millheads.	100	0.84	(0.3)*	-----	-----	-----
331	Copper, milltails.	100	.091	(.05)*	-----	-----	-----
332	Copper, concentrate	50	28.45	(10.2)*	-----	-----	-----
333	Molybdenum, concentrate ...	35	1.038	0.037	(+)	-----	-----

* parts per million

+ Nominal value 93% MoS₂



Cements

These SRM's are furnished for x-ray spectroscopic analysis and for chemical analysis of cements and related materials. Because these SRM's are hygroscopic, each unit consists of three sealed vials each containing approximately 5 g of material. The supply of the 1011, 1013-1016 Cements will soon be exhausted and replaced with seven new Cements, SRM's 633-639. (Values in parentheses are not certified, but are given for information only.)

Cements

SRM	Type	Wt/Unit (grams)	Chemical Composition (Nominal Weight Percent as the Oxide)				
			SiO ₂	Al ₂ O ₃	Fe ₂ O ₃	TiO ₂	P ₂ O ₅
633	Portland B (red)	}	NOW AVAILABLE				
634	Portland C (gold)						
635	Portland D (blue)						
636	Portland F (yellow)						
637	Portland G (pink)						
638	Portland I (green)	}	NOW AVAILABLE				
639	Portland J (clear)						
1011	Portland	15	21.03	5.38	2.07	0.25	0.33
1013	Portland	15	24.17	3.30	3.07	.20	.20
1014	Portland	15	19.49	6.38	2.50	.25	.32
1015	Portland	15	20.65	5.04	3.27	.26	.05
1016	Portland	15	21.05	4.97	3.71	.34	.13

Fertilizer Standards

These SRM's are intended for use in the fertilizer industry as working standards for the determination of the certified constituents.

SRM	Type	Wt/Unit (g)	Certified Composition (Wt percent)		
			N	P	K
193	Potassium Nitrate	90	13.85	-----	38.66
194	Ammonium Dihydrogen Phosphate	90	12.15	29.92	-----

Minerals, Refractories, Carbides, and Glasses

These SRM's are supplied in the form of powders, usually 100 mesh or finer. They are intended to provide materials for checking the accuracy of methods used in the analysis of similar materials, primarily in the glass, ceramics, and steel industries.

Minerals

SRM	Type	Wt/Unit (grams)	Chemical Composition (Nominal Weight Percent as the Oxide)					
			SiO ₂	Fe ₂ O ₃	Al ₂ O ₃	TiO ₂	MnO	CaO
1b	Limestone, argillaceous	50	4.92	0.75	1.12	0.046	0.20	50.9
88a	Limestone, dolomitic	50	1.20	.28	0.19	.02	.03	30.1
70a	Feldspar, potash	40	67.1	.075	17.9	.01	----	0.11
99a	Feldspar, soda	40	65.2	.065	20.5	.007	----	2.14
97a	Clay, flint	60	43.7	.45	38.8	1.90	----	0.11
98a	Clay, plastic	60	48.9	1.34	33.2	1.61	----	.31
81a	Glass sand	IN PREP						
165a	Glass sand (low iron)	IN PREP						
154b	Titanium dioxide	90	-----	-----	-----	99.74	-----	----

Cements (Continued)

CaO (+SrO)	SrO	MgO	SO ₃	Mn ₂ O ₃	Na ₂ O	K ₂ O	Li ₂ O	Rb ₂ O	Loss on Ignition	SRM
										633
										634
										635
										636
										637
										638
										639
66.60	0.11	1.12	1.75	0.03	0.08	0.26	(0.002)	(0.001)	1.13	1011
64.34	.08	1.39	1.80	.05	.20	.32	(.001)	(.004)	0.99	1013
63.36	.26	2.80	2.70	.07	.24	.99	(.005)	(.007)	.81	1014
61.48	.11	4.25	2.28	.06	.16	.87	(.004)	(.005)	1.70	1015
65.26	.25	0.42	2.27	.04	.55	.04	(.012)	(<.001)	1.20	1016



Minerals (Continued)

SrO	MgO	Cr ₂ O ₃	Na ₂ O	K ₂ O	Li ₂ O	ZrO ₂	BaO	Rb ₂ O	P ₂ O ₅	CO ₂	Loss on Ignition	SRM
0.14	0.36	----	0.04	0.25	----	----	----	----	0.08	40.4	41.1	1b
.01	21.3	----	.01	.12	----	----	----	----	.01	46.6	46.7	88a
----	----	----	2.55	11.8	----	----	0.02	0.06	----	----	0.40	70a
----	0.02	----	6.2	5.2	----	----	.26	----	.02	----	0.26	99a
.18	.15	0.03	0.037	0.50	0.11	0.063	.078	----	.36	----	13.32	97a
.039	.42	.03	.082	1.04	.070	.042	.03	----	.11	----	12.44	98a
												81a
												165a
----	----	----	----	----	----	----	----	----	----	----	----	154b

Refractories

Chemical Composition
(Nominal Weight Percent as the Oxide)

SRM	Type	Wt/Unit (grams)	SiO ₂	Al ₂ O ₃	Total as Fe ₂ O ₃	FeO	TiO ₂
76a	Burned Refractory (Al ₂ O ₃ -40%)	IN PREP					
77a	Burned Refractory (Al ₂ O ₃ -60%)	IN PREP					
78a	Burned Refractory (Al ₂ O ₃ -70%)	IN PREP					
103a	Chrome refractory	60	4.6	29.96	----	12.43	0.22
198	Silica refractory	45	----	0.16	0.66	----	.02
199	Silica refractory	45	----	.48	.74	----	.06
104	Burned magnesite	60	2.54	.84	7.07	----	.03

SRM	ZrO ₂	MnO	P ₂ O ₅	Cr ₂ O ₃	CaO	MgO	Li ₂ O	Na ₂ O	K ₂ O	Loss on Ignition
76a										
77a										
78a										
103a	0.01	0.11	0.01	32.06	0.69	18.54	-----	-----	-----	----
198	<.01	.008	.022	-----	2.71	0.07	0.001	0.012	0.017	0.21
199	.01	.007	.015	-----	2.41	.13	.002	.015	.094	.17
104	----	.43	.057	0.026	3.35	85.67	.001	.015	.015	----

Carbides

Chemical Composition
(Nominal Weight Percent)

SRM	Type	Wt/Unit (grams)	Carbon		Silicon		Fe	Al	Ti	Zr	Ca	Mg
			Total	Free	Total	SiC						
112	Silicon carbide	85	29.10	0.09	69.11	96.85	0.45	0.23	0.025	0.027	0.03	0.02

Glasses

Chemical Composition (Nominal Weight Percent)

SRM	Type	Wt/Unit (grams)	SiO ₂	PbO	Al ₂ O ₃	Fe ₂ O ₃	ZnO	MnO	TiO ₂	ZrO ₂
89	Lead-Barium	45 g	65.35	17.50	0.18	0.049	----	0.088	0.01	0.005
91	Opal	45 g	67.53	0.097	6.01	.081	0.08	.008	.019	.0095
92	Low-Boron	45 g	-----	-----	-----	-----	-----	-----	-----	-----
93a	High-Boron	Wafers 32 mm D x 6 mm	80.8	-----	2.3	.029	-----	-----	.012	.03
620	Soda-Lime, Flat	3 platelets 35 x 35 x 3 mm	72.1	-----	1.8	.04	-----	-----	.02	-----
621	Soda-Lime, Container	IN PREP								

SRM	CaO	BaO	MgO	K ₂ O	Na ₂ O	B ₂ O ₃	P ₂ O ₅	As ₂ O ₅	As ₂ O ₃	SO ₃	Cl	F	Loss on Ignition
89	0.21	1.40	0.03	8.40	5.70	-----	0.23	0.36	0.03	0.03	0.05	----	0.32
91	10.48	----	.008	3.25	8.48	-----	.022	.102	.091	----	.014	5.72	----
92	-----	-----	-----	-----	-----	0.70	-----	-----	-----	-----	-----	-----	----
93a	<0.02	----	<.01	0.01	4.0	12.6	-----	-----	-----	-----	-----	-----	----
620	7.1	----	3.7	.4	14.4	----	-----	-----	.06	.3	-----	-----	----
621													

Trace Element Standards

These SRM's are intended for use in calibrating instruments and in checking analytical techniques and procedures employed in the determination of trace elements in various inorganic matrices. [Note: For trace elements in biological matrices, see Biological Standards, page 43.]

SRM	Type - Matrix	Size	Unit of Issue
606	Trace Elements in Calcium Carbonate		10 gram
607	Trace Elements in Feldspar		5 gram
608	Trace Elements in Glass, Set	Wafers 3 mm Diameter	Set: 2 each 614 and 616
609	Trace Elements in Glass, Set	Wafers 1 mm Diameter	Set: 2 each 615 and 617
610	Trace Elements in Glass, 500 ppm	Wafers 3 mm Diameter	6 Wafers
611	Trace Elements in Glass, 500 ppm	Wafers 1 mm Diameter	6 Wafers
612	Trace Elements in Glass, 50 ppm	Wafers 3 mm Diameter	6 Wafers
613	Trace Elements in Glass, 50 ppm	Wafers 1 mm Diameter	6 Wafers
614	Trace Elements in Glass, 1 ppm	Wafers 3 mm Diameter	6 Wafers
615	Trace Elements in Glass, 1 ppm	Wafers 1 mm Diameter	6 Wafers
616	Trace Elements in Glass, 0.02 ppm	Wafers 3 mm Diameter	6 Wafers
617	Trace Elements in Glass, 0.02 ppm	Wafers 1 mm Diameter	6 Wafers
618	Trace Elements in Glass, Set	Wafers 3 mm Diameter	Set: 1 each 610, 612, 614 and 616
619	Trace Elements in Glass, Set	Wafers 1 mm Diameter	Set: 1 each 611, 613, 615 and 617

NOTE: Glass - Nominal Composition: 72% SiO₂, 12% CaO, 14% Na₂O, and 2% Al₂O₃.

Element	606	607	610-611 500 ppm	612-613 50 ppm	614-615 1 ppm	616-617 0.02 ppm
Antimony	-----	-----	-----	-----	(1.06)	(0.078)
Barium	-----	-----	-----	(41)	-----	-----
Boron	-----	-----	(35.1)	(32)	(1.30)	(0.20)
Cadmium	-----	-----	-----	-----	(0.55)	-----
Cerium	-----	-----	-----	(39)	-----	-----
Chromium	-----	-----	-----	-----	(0.99)	-----
Cobalt	-----	-----	(390)	(35.5)	0.71	-----
Copper	-----	-----	(444)	(37.7)	1.34	(0.65)
Dysprosium	-----	-----	-----	(35)	-----	-----
Erbium	-----	-----	-----	(39)	-----	-----
Eurpoium	-----	-----	-----	(36)	(0.99)	-----
Gadolinium	-----	-----	-----	(39)	-----	-----
Gallium	-----	-----	-----	-----	(1.3)	(0.23)
Gold	-----	-----	(25)	(5)	(0.5)	(0.18)
Indium	-----	-----	-----	-----	(0.75)	(0.26)
Iron	-----	-----	458	51	13.5	(11)
Lanthanum	-----	-----	-----	(36)	(0.83)	(0.034)
Lead	(0.374)	-----	426	38.57	2.32	1.85
Manganese	-----	-----	485	(39.6)	(1.41)	(0.65)
Molybdenum	-----	-----	(111)	-----	-----	-----
Neodymium	-----	-----	-----	(36)	-----	-----
Nickel	-----	-----	458.7	38.8	(0.95)	-----
Potassium	-----	-----	(461)	(64)	30	29
Rhenium	-----	-----	(49.1)	-----	-----	-----
Rubidium	-----	523.90	425.7	31.4	0.855	(0.0998)
Samarium	-----	-----	-----	(39)	-----	-----
Scandium	-----	-----	-----	-----	(0.59)	(0.026)
Silver	-----	-----	(254)	22.0	0.42	-----
Strontium	-----	65.485	515.5	78.4	45.8	41.72
Thallium	-----	-----	(61.8)	(15.7)	(0.269)	(0.0082)
Thorium	-----	-----	457.2	37.79	0.748	0.0252
Titanium	-----	-----	(437)	(50.1)	(3.1)	(2.5)
Uranium	-----	-----	461.5	37.38	0.823	0.0721
Ytterbium	-----	-----	-----	(42)	-----	-----
Zinc	-----	-----	(433)	-----	(2.43)	-----

In addition to the 35 elements listed above, the Glass SRM's contain the following 26 elements: As, Be, Bi, Cs, Cl, F, Ge, Hf, Hg, Li, Lu, Mg, Nb, P, Pr, Se, S, Ta, Te, Tb, Tm, Sn, W, V, Y, and Zr.

Nuclear Materials

Special Nuclear Materials

These SRM's consist of four groups: Plutonium Assay Standards, Plutonium Isotopic Standards, Uranium Assay Standards, and Uranium Isotopic Standards.

These SRM's are available to AEC contractors, AEC or State Licensees, and foreign governments that have entered an Agreement for Cooperation with the U.S. Government concerning the Civil Uses of Atomic Energy. The purchase request for these SRM's must be made on special forms obtainable from the Office of Standard Reference Materials, Room B311, Chemistry Building, National Bureau of Standards, Washington, D.C. 20234.

Plutonium Assay Standards

SRM	Type	Certified for	Wt/Units (grams)	Purity (%)
944	Plutonium sulfate tetrahydrate	Plutonium Content	0.5	47.50*
945	Plutonium metal, standard matrix	Impurities	5	(99.9)
949d	Plutonium metal assay	Plutonium Content	0.5†	99.99

*Stoichiometric

†Nominal weight, each SRM is issued with an individual weight.

(Values in parentheses are not certified, but are given for information only.)

Plutonium Isotopic Standards

SRM	Type	Wt/Units (grams)	Atom Percent				
			²³⁸ Pu	²³⁹ Pu	²⁴⁰ Pu	²⁴¹ Pu	²⁴² Pu
946	Plutonium Sulfate Tetrahydrate	0.25	0.247	83.128	12.069	3.991	0.565
947	Plutonium Sulfate Tetrahydrate25	.296	75.696	18.288	4.540	1.180
948	Plutonium Sulfate Tetrahydrate25	.011	91.574	7.914	0.468	0.0330

Uranium Assay Standards

SRM	Type	Certified For	Wt/Unit (grams)	Purity (%)
950a	Uranium Oxide	Uranium Oxide	25	99.94 (U ₃ O ₈)
960	Uranium Metal	Uranium	26	99.975 (U)

Uranium Isotopic Standards

SRM	Uranium Oxide (U ₃ O ₈)	Wt (grams)	Atom Percent			
			²³⁴ U	²³⁵ U	²³⁶ U	²³⁸ U
U-0002	Depleted	1.0	0.00016	0.01755	<0.00001	99.9823
U-005	Depleted	1.0	.00218	.4895	.0046	99.504
U-010	Enriched	1.0	.00541	1.0037	.00681	98.984
U-015	Enriched	1.0	.00850	1.5323	.0164	98.443
U-020	Enriched	1.0	.0125	2.038	.0165	97.933
U-030	Enriched	1.0	.0190	3.046	.0204	96.915
U-050	Enriched	1.0	.0279	5.010	.0480	94.915
U-100	Enriched	1.0	.0676	10.190	.0379	89.704
U-150	Enriched	1.0	.0993	15.307	.0660	84.528
U-200	Enriched	1.0	.1246	20.013	.2116	79.651
U-350	Enriched	1.0	.2498	35.190	.1673	64.393
U-500	Enriched	1.0	.5181	49.696	.0755	49.711
U-750	Enriched	1.0	.5923	75.357	.2499	23.801
U-800	Enriched	1.0	.6563	80.279	.2445	18.820
U-850	Enriched	1.0	.6437	85.137	.3704	13.848
U-900	Enriched	1.0	.7777	90.196	.3327	8.693
U-930	Enriched	1.0	1.0812	93.336	.2027	5.380
U-970	Enriched	1.0	1.6653	97.663	.1491	0.5229

Neutron Density Standards

These SRM's are provided as reference sources of a cobalt-in-aluminum alloy to serve as a neutron density monitor wire SRM. Accurate determination of thermal neutron densities is essential in irradiation tests to obtain a basis for comparison of densities among reactors, in applying data in the design of reactors, in understanding the mechanisms of radiation damage, and for use in neutron activation analysis. The wire is 0.5 mm in diameter and is available in four lengths.

SRM	Type	Unit	Cobalt Content (Weight percent)
953	Neutron density monitor wire (Co in Al)	1 meter	0.116
953L1	Neutron density monitor wire (Co in Al)	5 meters	.116
953L2	Neutron density monitor wire (Co in Al)	10 meters	.116
953L3	Neutron density monitor wire (Co in Al)	25 meters	.116

Isotopic Reference Standards


SRM's for isotopic ratio are natural-ratio materials, unless otherwise noted, and are furnished with a certificate of isotopic composition.

The isotopic composition of these SRM's has been determined by mass spectrometry, by comparison with mixtures prepared from high-purity separated isotopes. They are useful for those looking for small variations in the isotopic composition of the elements, and for the evaluation of mass discrimination effects encountered in the operation of mass spectrometers.

A description of the methods used in the characterization of the normal and enriched boric acid, SRM's 951 and 952, is given in NBS Special Publication 260-17. (See inside of back cover for ordering instructions.)

SRM	Isotopic Reference Standards	Element Certified	Wt/Unit (grams)
951	Boric Acid	Boron	1.0
952	Boric Acid, 95% Enriched ^{10}B	Boron	100
975	Sodium Chloride	Chlorine	0.25
976	Copper Metal	Copper	.25
977	Sodium Bromide	Bromine	.25
978	Silver Nitrate	Silver	.25
979	Chromium Nitrate	Chromium	.25
980	Magnesium Metal	Magnesium	.25
*981	Lead Metal, Natural	Lead	1.0
*982	Lead Metal, Equal Atom (206/208)	Lead	1.0
*983	Lead Metal, Radiogenic (92%-206)	Lead	1.0
984	Rubidium Chloride, assay and isotopic	Rubidium	0.25
987	Strontium Carbonate, assay and isotopic	Strontium	.25
988	Strontium-84 Spiked, assay and isotopic	Strontium	10

*Sold as a set only of three 981, 982, and 983



Ion Activity Standards

These SRM's are intended for use in the preparation of solutions for the calibration of specific-ion electrodes. This includes the pH and pD measuring systems.

pH Standards

These SRM's are furnished as crystals for the preparation of solutions of known hydrogen ion concentration for calibrating and checking the performance of commercially available pH materials and instruments. They are furnished with certificates giving directions for preparation of the solutions and tables of pH values at various temperatures.

SRM's 186Ic and 186IIc, 191 and 192, and 922 and 923, are certified for use in admixture only. At an equimolar (0.025 molal) mixture of SRM's 186Ic and 186IIc, a pH(S) of 6.863 at 25 °C is obtained. Directions also are furnished for the preparation of a physiological reference solution from 186Ic and 186IIc having a pH(S) of 7.415 at 25 °C.

SRM	Type	pH(S) (at 25 °C)	Wt/Unit (grams)
185e	Potassium Acid Phthalate	4.004	60
186Ic	Potassium dihydrogen phosphate	{ 6.863 }	30
186IIc	Disodium hydrogen phosphate		30
187b	Borax	9.183	30
188	Potassium hydrogen tartrate	3.557	60
189	Potassium tetroxalate	1.679	65
191	Sodium bicarbonate	10.01	30
192	Sodium carbonate		30
922	Tris(hydroxymethyl)aminomethane	7.699	25
923	Tris(hydroxymethyl)aminomethanehydrochloride		35

pD Standards

These SRM's are furnished as crystals for preparation of solutions of known deuterium-ion concentration for the calibration and correction of pH indicating equipment to indicate pD data. SRM's 2186I and 2186II, and 2191 and 2192, are certified for use in admixtures only.

SRM	Type	pD(S) Values	Wt/Unit (grams)
2186-I	Potassium dihydrogen phosphate	7.43	30
2186-II	Disodium hydrogen phosphate		30
2191	Sodium bicarbonate	10.74	30
2192	Sodium carbonate		30

Ion - Selective Electrodes

These SRM's are certified for the calibration of ion-selective electrodes and have conventional ionic activities based on the Stokes-Robinson hydration theory for ionic strengths greater than 0.1 mole per liter.

SRM	Type	Certified Property	Wt/Unit (grams)
2201	Sodium Chloride	pNa, pCl pK, pCl	125
2202	Potassium Chloride		160

Standards of Certified Physical Properties

Mechanical and Metrology Standards

These SRM's are intended to relate measurement units made in industrial, university, and government laboratories to the mechanical and metrological units related through a National Measurement System¹ to the base units of mass, length, and time.

Coating Thickness Standards

These SRM's have a specimen size of 3 X 3 cm and are for calibrating coating thickness gages of the magnetic type for the measurement of thickness of nonmagnetic coatings on steel, nickel on steel, or nickel on nonmagnetic substrates. The steel substrates have the magnetic properties of AISI 1010 steel and the nickel coatings have the magnetic properties of an annealed, Watts nickel electrodeposit free of cobalt and iron.

The magnetic type thickness gages are often used to measure the thickness of paint and other organic coatings on steel, as well as zinc (galvanized) and other nonmagnetic metallic coatings. SRM's in the 1301 to 1320 series (sets 1351, 1361-64) are used to calibrate these gages. The number of different thicknesses required for these calibrations depends on the type of gage and the coating thicknesses to be measured.

The magnetic type thickness gages can be used to estimate magnetic properties of austenitic stainless steel weld metal. Because the magnetic properties of the weld metal are closely related to the ferrite content of the weld, these instruments are used to estimate the ferrite content. For these measurements, the coating thickness SRM's 1312-1319 (as sets 1368-70) are used to calibrate the instrument. The ferrite contents having magnetic properties similar to those of the various coating thickness SRM's have been established by other laboratories. For sets 1351 to 1369, the specimens are mounted on one card. Set 1370 is mounted on two cards, but packed in one box.

SRM's with gold and tin coatings on various substrates have a specimen size of 15 X 15 mm and are for calibrating coating thickness gages of the beta-backscatter type and for calibrating x-ray fluorescence methods for the measurement of the weight per unit area of gold or tin coatings. For gold and tin sets, the specimens are mounted on separate cards, but are packed in one box.

The gold coating standards were measured by beta-ray backscatter and x-ray fluorescence techniques relative to NBS gold coating materials for which the average weights per unit area were determined by weight and area measurements. They are suitable for the direct calibration of equipment used to measure weight per unit area of gold coating of equivalent purity. From the density and weight per unit area, the instruments can be calibrated in terms of the thickness of the standard.

For the tin coating standards, x-ray fluorescence techniques were used to measure the thickness of the tin coating relative to NBS tin coating material for which the average weights per unit area were determined by weight and area measurements.

ASTM Methods of Measuring Coating Thickness

Instrumental methods of measuring coating thickness are set forth in the following ASTM Methods of Test. [NOTE: Metric units. ASTM plating specifications use μm and mil, but do not go to thicknesses greater than about 75 μm : 1 mil = 0.001 inch = 25.4 μm .]

- B499 Method of Measurement of Coating Thicknesses by the Magnetic Method; Nonmagnetic Coatings on Magnetic Basis Metals.
- B529 Measurement of Coating Thicknesses by the Eddy-Current Test Method: Nonconductive Coatings on Nonmagnetic Basis Metals.
- B530 Measurement of Coating Thicknesses by the Magnetic Method: Electrodeposited Nickel Coatings on Magnetic and Nonmagnetic Substrates.
- B504 Measuring the Thickness of Metallic Coatings by the Coulometric Method.
- B244 Measuring Thickness of Anodic Coatings on Aluminum with Eddy-Current Instruments.
- E376 Recommended Practice for Measuring Coating Thickness by Magnetic-Field or Eddy-Current (Electromagnetic) Test Methods.
- D1400 Measurement of Dry Film Thickness of Nonmetallic Coatings of Paint, Varnish, Lacquer, and Related Products Applied on a Nonmagnetic Metal Base.
- D1186 Measurement of Dry Film Thickness of Nonmagnetic Organic Coatings Applied on a Magnetic Base.

¹ "Concept of a National Measurement System," Science 158, 67-71 (1967).

Nonmagnetic Coatings on Steel, or Nickel on Nonmagnetic Substrate

Nonmagnetic Coating on Magnetic Substrate (Copper and Chromium on Steel)

SRM	Nominal Coating Thickness		
	(inch)	(metric)	(mil)
1301	0.00010	2.54 μm	0.10
1302	.00025	6.25 μm	.25
1303	.00050	12 μm	.50
1304	.00075	19 μm	.75
1305	.0010	25 μm	1.0
1306	.0015	38 μm	1.5
1307	.0020	51 μm	2.0
1308	.0025	62 μm	2.5
1309	.0027	68 μm	2.7
1310	.0032	81 μm	3.2
1311	.0055	0.14 mm	5.5
1312	.0080	.23 mm	8.0
1313	.010	.25 mm	10.0
1314	.015	.38 mm	15.0
1315	.020	.51 mm	20.0
1316	.025	.62 mm	25.0
1317	.03	.76 mm	30.0
1318	.04	1.0 mm	40.0
1319	.06	1.5 mm	60.0
1320	.08	2.0 mm	80.0
1351	Set of 2: 1307 and 1311		
1361	Set of 4: 1302, 1303, 1305, and 1307		
1362	Set of 4: 1306, 1310, 1311, and 1312		
1363	Set of 4: 1313, 1314, 1315, and 1316		
1364	Set of 4: 1317, 1318, 1319, and 1320		
1368	Set of 4: 1312, 1313, 1314, and 1315		
1369	Set of 4: 1316, 1317, 1318, and 1319		
1370	Set of 8: 1312, 1313, 1314, 1315, 1316, 1317, 1318, and 1319		

Magnetic Coating on Magnetic Substrate (Nickel on Steel)

SRM	Nominal Coating Thickness		
	(inch)	(micrometer)	(mil)
1331	0.00012	3.0	0.12
1332	.00035	8.9	.35
1333	.00055	14	.55
1334	.00075	19	.75
1335	.0010	25	1.0
1336	.0013	33	1.3
1337	.0016	40	1.6
1338	.0020	51	2.0
1339	.0025	62	2.5
1352	Set of 2: 1332 and 1334		
1353	Set of 2: 1335 and 1339		
1365	Set of 4: 1331, 1332, 1333, and 1334		
1366	Set of 4: 1335, 1336, 1337, and 1338		

Magnetic Coating on Non-Magnetic Substrate (Nickel and Chromium on Brass)

SRM	Nominal Coating Thickness		
	(inch)	(micrometer)	(mil)
1341	0.00012	3.0	0.12
1342	.00035	8.9	.35
1343	.00065	16.5	.65
1344	.0010	25	1.0
1345	.0015	38	1.5
1346	.0020	51	2.0
1367	Set of 4: 1341, 1342, 1343, and 1344		

Gold and Tin Coating Thickness Standards

Gold Coating on Glass Sealing Alloy – ASTM Designation F15; Fe-53, Ni-29, and Co-17.

SRM	Nominal Coating Weight (mg/cm ²)	Nominal Coating Thickness	
		(micrometers)	(microinches)
1371	1.5	0.8	30
1372	3.0	1.5	60
1373	6.0	3.0	120
1374	14.0	7.1	280
1381	Set of 2: 1371 and 1372		
1382	Set of 2: 1372 and 1373		
1383	Set of 2: 1373 and 1374		
1398	Set of 4: 1371, 1372, 1373, and 1374		

Gold Coating on Nickel

SRM	Nominal Coating Weight (mg/cm ²)	Nominal Coating Thickness	
		(micrometers)	(microinches)
1375	1.5	0.8	30
1376	3.0	1.5	60
1377	6.0	3.0	120
1378	17.0	8.9	350
1384	Set of 2: 1375 and 1376		
1385	Set of 2: 1376 and 1377		
1386	Set of 2: 1377 and 1378		
1399	Set of 4: 1375, 1376, 1377, and 1378		

Gold Coating on Copper-Clad, Glass-Epoxy Laminate

SRM	Nominal Coating Weight (mg/cm ²)	Nominal Coating Thickness	
		(micrometers)	(microinches)
2301	1.5	0.8	30
2302	3.0	1.5	60
2303	6.0	3.0	120
2304	14.0	7.1	280
2305	Set of 2: 2301 and 2302		
2306	Set of 2: 2302 and 2303		
2307	Set of 2: 2303 and 2304		
2308	Set of 4: 2301, 2302, 2303, and 2304		

Gold Coating on Copper

SRM	Nominal Coating Weight (mg/cm ²)	Nominal Coating Thickness	
		(micrometers)	(microinches)
2311	1.5	0.8	30
2312	3.0	1.5	60
2313	6.0	3.0	120
2314	14.0	7.1	280
2315	Set of 2: 2311 and 2312		
2316	Set of 2: 2312 and 2313		
2317	Set of 2: 2313 and 2314		
2318	Set of 4: 2311, 2312, 2313, and 2314		

Tin Coating on Steel

SRM	Nominal Coating Weight (mg/cm ²)	Nominal Coating Thickness	
		(micrometers)	(microinches)
2331	1.1	1.5	60
2332	2.0	2.8	110
2333	3.0	4.0	160
2334	5.0	7.0	275
2335	12.0	16.5	650
2336	14.0	19.0	750
2338	Set of 2: 2332 and 2335		
2339	Set of 4: 2331, 2333, 2334, and 2336		
2340	Set of 6: 2331, 2332, 2333, 2334, 2335, and 2336		

Density Standards

These SRM's are certified with respect to values of density, for air-saturated material at 1 atm, at 20, 25, and 30 °C, to ± 0.00002 g/ml. These SRM's may be used to calibrate picnometers and density balances. A Certificate is supplied with each of these SRM's. SRM 217b-8S is contained in a special ampoule with an internal breakoff tip, the others are sealed "in vacuum" in plain glass ampoules.

SRM	Type	Approx.d ²⁰	Amount, ml
217b-5	2,2,4-Trimethylpentane	0.6918	5
217b-8S	2,2,4-Trimethylpentane6918	8
217b-25	2,2,4-Trimethylpentane6918	25
217b-50	2,2,4-Trimethylpentane6918	50

Glass Viscosity Standards

SRM's 710, 711, and 717 are furnished as rectangular-shaped bars, and are certified for viscosity between values of 10^2 and 10^{12} poises. They are furnished to check the performance of high-temperature viscosity equipment (rotating cylinders) and low-temperature viscosity equipment (fiber elongation). In addition, values are furnished for the softening point, annealing point, and strain point by ASTM Designations (C338-61 and C336-61). Certificates of Data from participating laboratories are furnished for each of the three glasses.

SRM's 712, 713, 714, 715, and 716 are furnished in cane, gobs, or patties as listed, and are certified only for softening point, annealing point, and strain point. Certificates of Data from three laboratories are furnished for each of these glasses.

SRM	Type	Unit of Issue
710	Soda-lime silica glass-type 523/586	2 lb
711	Lead-silica glass-type 617/366	3 lb
712	Mixed alkali lead silicate glass, 1/4 in patties (6 pcs.)	0.5 lb
713	Dense barium crown 620/603 glass, 1 3/8 in diam \times 5/8 in thick gobs (4 pcs.)5 lb
714	Alkaline earth alumina silicate glass, 1/4 in diam cane (16 pcs-6 in long)5 lb
715	Alkali-free aluminosilicate glass, 1/4 in diam cane (13 pcs-6 in long)	200 g
716	Neutral (borosilicate) glass, 1/2 in diam cane (6 pcs-6 in long)	250 g
717	Borosilicate glass, 4.2 cm \times 4.2 cm \times 12.5 cm bar	500 g

SRM	Viscosity (Poises at Indicated Temperature (°C))											Softening Point °C	Annealing Point °C	Strain Point °C
	10 ²	10 ³	10 ⁴	10 ⁵	10 ⁶	10 ⁷	10 ⁸	10 ⁹	10 ¹⁰	10 ¹¹	10 ¹²			
710	1434.3	1181.7	1019.0	905.3	821.5	757.1	706.1	664.7	630.4	601.5	576.9	724	546	504
711	1327.1	1072.8	909.0	794.7	710.4	645.6	594.3	552.7	518.2	489.2	464.5	602	432	392
712	----	----	----	----	----	----	----	----	----	----	----	528	386	352
713	----	----	----	----	----	----	----	----	----	----	----	738	631	599
714	----	----	----	----	----	----	----	----	----	----	----	908	710	662
715	----	----	----	----	----	----	----	----	----	----	----	961	764	714
716	----	----	----	----	----	----	----	----	----	----	----	794	574	530
717	1545.1	1248.8	1059.4	927.9	831.2	757.1	698.6	651.1	611.9	579.0	550.9	720	516	471

Polymer Standards

Four Polymer SRM's are available: Two polystyrenes, one with a narrow molecular weight distribution (SRM 705), the other with a broad distribution (SRM 706), and two polyethylenes, linear (SRM 1475) and branched (SRM 1476).

These materials are certified for the properties indicated in the table, such as weight and number average molecular weight, molecular weight distribution, limiting viscosity numbers (intrinsic viscosities) in several solvents, density, and melt flow.

These SRM's have wide application not only in the calibration of instruments used in polymer characterization, such as light scattering photometers, osmometers, gel permeation chromatographs, but also wherever a well characterized polymer material is needed, as for example in studies of dilute solution behavior, rheology, and polymer crystal physics.

The certificate for SRM 1475 is accompanied by a series of papers, reprinted from the Journal of Research of the National Bureau of Standards, which describe how the measurements were obtained.

SRM	Type	Wt/Unit (grams)
705	Polystyrene, narrow molecular weight distribution, $M_w \approx 179,000$, $M_w/M_n \approx 1.07$	5
706	Polystyrene, broad molecular weight distribution, $M_w \approx 258,000$, $M_w/M_n \approx 2.1$	18
1475	Polyethylene, linear, $M_w \approx 52,000$, $M_w/M_n = 2.9$	50
1476	Polyethylene, branched.	50

The following table lists the properties (and method) certified for these SRM's.

Property (and method)		705	706	1475	1476
Molecular Weight					
Weight Average	(Light Scattering)	X	X	X	---
	(Sedimentation Equilibrium)	X	X	---	---
	(Gel Permeation Chromatography-GPC)	---	---	X	---
Number Average	(Osmometry)	X	---	---	---
	(GPC)	---	---	X	---
Molecular Weight Distribution	(GPC)	---	---	X	---
Limiting Viscosity Number	(Capillary Viscometer)				
Benzene 25 °C		X	X	---	---
Benzene 35 °C		X	---	---	---
Cyclohexane 35 °C		X	X	---	---
1-Chloronaphthalene 130 °C		---	---	X	X
1,2,4-trichlorobenzene 130 °C		---	---	X	X
Decahydronaphthalene 130 °C		---	---	X	X
Melt Flow	(ASTM)	---	---	X	X
Density	(ASTM)	---	---	X	X

Elasticity Standards

This SRM is polycrystalline alumina prepared from a single block of material by isostatically cold pressing and then sintering alumina powder containing 0.1 percent magnesium oxide. It is intended for the calibration of apparatus used in the measurement of resonance frequencies from which elastic moduli are calculated. Each bar has been individually measured and calibrated, and all surfaces were machined flat and parallel.

SRM	Type	Size
718	Polycrystalline Alumina	12.7 × 1.27 × 0.32 cm

Heat Standards

These SRM's are intended to relate heat and temperature measurements made in industrial, university, and government laboratories with the International Practical Temperature Scale-1968¹.

Superconductive Thermometric Fixed Point Devices

The SRM is a device composed of small cylinders of high purity lead, indium, aluminum, zinc, and cadmium mounted in a threaded copper stud and enclosed by a mutual inductance coil set. It provides convenient temperature calibrations in the range 0.5 - 7 K with precision of ± 0.001 K. It should prove particularly valuable to users of ³He-⁴He dilution refrigerators, in which direct calibrations on the liquid helium vapor pressure-temperature scales are difficult, and to those who wish to determine the temperature reproducibility of physical phenomena or of cryogenic equipment.

SRM	Type	Element	Nominal Temperature (K)
767	Superconductive Thermometric Fixed Point Device	Lead	7.2
		Indium	3.4
		Aluminum	1.2
		Zinc	0.8
		Cadmium	.5

Freezing Point Standards

Defining Fixed Points – International Practical Temperature Scale

These SRM's are of such purity that they are suitable for defining fixed points for the International Practical Temperature Scale of 1968.¹

SRM	Type	Temperature C	Wt/Unit (grams)
740	Zinc	419.58	350
741	Tin	231.968 1	350

Determined Reference Points

These SRM's are intended for use in calibration of thermometers, thermocouples, and other temperature measuring devices. The temperatures certified are in accord with the International Practical Temperature Scale of 1968.

SRM	Type	Temperature C	Wt/Unit (grams)
42g	Tin	231.967	350
43h	Zinc	419.58*	350
44f	Aluminum	660.3	200
45d	Copper	1084.8	450
49e	Lead	327.493	600

* SRM 43h is less pure than SRM 740 and has a freezing point 0.001 C lower.

¹ "International Practical Temperature Scale of 1968," Metrologia, 5 35-44 (1969).

Melting Point Standard

This SRM is calcined alpha alumina the purity of which (99.9+ percent) makes it a suitable pyrometric standard for melting point on the International Practical Temperature Scale of 1968.

SRM	Type	Temperature °C	Wt/Unit (grams)
742	Alumina	2053	10

Calorimetric Standards

These SRM's are intended to relate the gain or loss of energy and work experienced during a chemical reaction or by change of temperature to the units of energy and work as defined by the National Measurement System. This system uses the units prescribed by the International System of Units (SI). The unit for energy and work under this system is the joule, which is related to the historically defined calorie by the equation: 4.184 joule = 1 calorie.

Combustion Calorimetric Standards

These SRM's are issued primarily to check the performance of calorimetric methods for the determination of the heat of combustion. SRM 217b-8S is contained in a special ampoule with an internal break-off tip, the others are sealed "in vacuum" in plain glass ampoules.

SRM	Type	Unit Amount
39i	Benzoic acid, 26.434 absolute kilojoules/gram	30 g
217b-5	2,2,4-Trimethylpentane, 47.713 absolute kilojoules/gram	5 ml
217b-8S	2,2,4-Trimethylpentane	8 ml
217b-25	2,2,4-Trimethylpentane	25 ml
217b-50	2,2,4-Trimethylpentane	50 ml

Solution Calorimetric Standards

These SRM's are issued primarily to check the performance of calorimetric methods used for the determination of heats of solution and heats of reactions in solution.

SRM	Type	Wt/Unit (grams)
724	tris(hydroxymethyl)aminomethane	50
1654	α -Quartz for HF acid solution calorimetry	25

Heat Source Calorimetric Standards

SRM	Type	Wt/Unit (grams)
1651	Zirconium-barium chromate heat source powder (ca 350 cal/g)	50
1652	Zirconium-barium chromate heat source powder (ca 390 cal/g)	50
1653	Zirconium-barium chromate heat source powder (ca 425 cal/g)	50

Enthalpy and Heat Capacity Standards

This SRM has been certified for enthalpy and heat capacity of 99.95+ percent α -alumina over a temperature range from 273.15 to 2250 K.

SRM	Type	Wt/Unit (grams)
720	Sapphire, synthetic (Al_2O_3)	15 g

Differential Thermal Analysis Standards

The first two of these SRM's, 755 and 756, were issued upon the recommendation of the International Confederation of Thermal Analysis (ICTA) as standards for checking the performance of differential thermal analysis equipment. SRM's 758, 759, and 760 have been issued by NBS in cooperation with ICTA as standards for calibrating differential thermal analysis and related thermo-analytical equipment under operating conditions. SRM's 758, 759, and 760 comprise a total of eight inorganic substances and two metals.

SRM	Type	Temperature/Range (°C)	Unit of Issue
755	Quartz SiO_2	575	2 g
756	Potassium Nitrate	130	5 g
758	DTA Temperature Standard	125-435	Set of 5 (See below)
759	DTA Temperature Standard	295-675	Set of 5 (See below)
760	DTA Temperature Standard	570-940	Set of 5 (See below)

758 (125-435 °C)	759 (295-675 °C)	760 (570-940 °C)	Peak Temp. °C	Wt(g)
KNO_3	-----	-----	135	10
In (Metal)	-----	-----	159	3
Sn (Metal)	-----	-----	237	3
KClO_4	KClO_4	-----	309	10
Ag_2SO_4	Ag_2SO_4	-----	433	3
-----	SiO_2	SiO_2	574	3
-----	K_2SO_4	K_2SO_4	588	10
-----	K_2CrO_4	K_2CrO_4	673	10
-----	---	BaCO_3	819	10
-----	---	SrCO_3	938	10

Vapor Pressure Standards

These SRM's are intended for use in the testing and calibration of vapor pressure measurement apparatus and techniques. The materials ultimately will include gold, cadmium, platinum, silver, and tungsten, and will cover a temperature range of 600 to 3,000 K.

SRM	Type	Pressure Range (atmosphere)	Temperature Range (K)	Unit Size
745	Gold	10^{-3} to 10^{-8}	1300 - 2100	Wire 1.44 mm x 152 mm
746	Cadmium	10^{-4} to 10^{-11}	350 - 594	Rod 6.4 mm x 64 mm
747	Platinum	-----	-----	IN PREP
748	Silver	10^{-3} to 10^{-12}	800 - 1600	Rod 6.4 mm x 64 mm
749	Tungsten	-----	-----	IN PREP

Thermal Conductivity Standards

These SRM's cover the high, medium, and low conductivity ranges. They will be useful for intercomparing thermal conductivity apparatus, debugging new apparatus, and calibrating comparative apparatus.

SRM	Type	Temperature Range (K)	Diameter (mm)	Length (mm)
734-S	Electrolytic Iron	6-280	6.4	305
734-L1	Electrolytic Iron	6-280	31.8	152
734-L2	Electrolytic Iron	6-280	31.8	305
735-S	Stainless Steel	5-280	6.5	300
735-M1	Stainless Steel	5-280	12.5	150
735-M2	Stainless Steel	5-280	12.5	300
735-L1	Stainless Steel	5-280	35	50
735-L2	Stainless Steel	5-280	35	100

Thermal Expansion Standards


These SRM's cover the temperature range from 20 to 1900 K having coefficients of thermal expansion over the range of 0.5 to $25 \times 10^{-6}/K$.

SRM	Type	Temperature Range (K)	Diameter (mm)	Length (mm)
731-L1	Borosilicate Glass	80-680	6.4	51
731-L2	Borosilicate Glass	80-680	6.4	102
731-L3	Borosilicate Glass	80-680	6.4	152
732-L1	Sapphire	IN PREP		
732-L2	Sapphire	IN PREP		
732-L3	Sapphire	IN PREP		
736-L1	Copper	20-800	6.4	51
736-L2	Copper	20-800	6.4	102
736-L3	Copper	20-800	6.4	152
737-L1	Tungsten	IN PREP		
737-L2	Tungsten	IN PREP		
737-L3	Tungsten	IN PREP		
739-L1	Fused Silica	80-1000	6.4	51
739-L2	Fused Silica	80-1000	6.4	102
739-L3	Fused Silica	80-1000	6.4	152

Thermocouple Materials

This SRM is intended to serve as a convenient mechanism for the comparison of manufactured wire to standard reference thermocouple tables.

SRM	Type	Form
733	Silver-28 Atomic Percent Gold	Wire: 32AWG(0.2019 mm) diameter 3 meters long




Magnetic Standards

Magnetic Susceptibility Standards

These SRM's are intended for use in the calibration of instruments used to measure magnetic susceptibility. The susceptibility values shown are nominal for 297 K.

SRM	Type	Gram Susceptibility $\chi, \text{cm}^3 \cdot \text{g}^{-1}$	Volume Susceptibility k	Form/Unit	
763-1	Aluminum *	0.605×10^{-6}	1.63×10^{-6}	Cylinder	3 mm diameter \times 3 mm
763-2	Aluminum			Wire	0.5 mm diameter \times 250 mm
763-3	Aluminum			Rod	6 mm diameter \times 175 mm
764-1	Platinum	0.990×10^{-6}	21.2×10^{-6}	Cylinder	3 mm diameter \times 3 mm
764-2	Platinum			Wire	0.5 mm diameter \times 50 mm
765-1	Palladium	5.25×10^{-6}	63.1×10^{-6}	Cylinder	3 mm diameter \times 3 mm
765-2	Palladium			Wire	0.5 mm diameter \times 50 mm
765-3	Palladium			Sponge	1 gram
766-1	Manganese Fluoride	123.5×10^{-6}	484×10^{-6}	Cube	3 \times 3 \times 3 mm

*At 77.7K, Gram Susceptibility for SRM 763: $\chi = 0.695 \times 10^{-6} \text{ cm}^3 \cdot \text{g}^{-1}$



Optical Standards

Spectrophotometric Standards

Color Standards for Spectrophotometer-Tristimulus Integrator Systems

This SRM is a set of 5 transparent colored glass filters to check the performance of spectrophotometer-tristimulus integrator systems, the automatic recording and computing devices used in routine color measurements. Each glass filter is 2-inches (5 cm) square (approximately 3.0 mm thick) with polished faces. A chart of tristimulus values for CIE sources A, B, and C, representing incandescent-lamp light, noon sunlight, and average daylight; and a detailed report on the changes in tristimulus values caused by errors in the 100-percent and zero adjustments of the photometric scale, wavelength errors, slit-width errors, errors due to stray energy, and inertia errors of the recording mechanisms are furnished with each set of glasses. Through the use of these filters the user of a spectrophotometer-integrator combination will be able to determine when the instrument goes out of adjustment. From the pattern of the discrepancies between measured and reported tristimulus values, he will also be able to obtain some clue as to the type of maladjustment.

This SRM is available only as a set of five filters.

SRM	Type	Unit Size
2101	Orange-red glass	Supplied only as a set - one each of 5 filters
2102	Signal yellow glass	
2103	Sextant green glass	
2104	Cobalt blue glass	
2105	Selective neutral glass	

Glass and Liquid Filters for Spectrophotometry

These SRM's are intended primarily for use in checking the accuracy of the photometric scale of spectrophotometers and to provide a means of interlaboratory comparisons of spectrophotometric data.

The Glass Filters, SRM 930a, consists of three filters having transmittances of approximately 10, 20, and 30 percent. Each filter is individually calibrated and certified for absorbance and transmittance over a spectral wavelength range from 440 to 635 nanometers.

The Liquid Filters, SRM 931, are absorbance standards for use in ultraviolet and visible spectrophotometry. This SRM consists of three sets of four vials containing approximately 10 ml of solution. Each set contains a "Blank" and three concentrations of the liquid filter. The net absorbance of these filters are certified at 302, 395, 512, and 678 nanometers.

SRM	Type	Unit of Issue
930a	Glass filters for Spectrophotometry	Set of 3 filters
931	Liquid Absorbance Standards for Ultraviolet and Visible Spectrophotometry	3 Sets of 4 vials



Thermal Emittance Standards

SRM's of normal spectral emittance are available in three materials, platinum-13 percent rhodium alloy having low emittance, sandblasted and oxidized Kanthal (an iron-chromium-aluminum alloy) having intermediate emittance, and sandblasted and oxidized Inconel (a nickel-chromium-iron alloy) having high emittance. SRM's of all three materials have been calibrated for normal spectral emittance at 800 and 1100 K; the Kanthal and Inconel standards at 1300 K and the platinum-13 percent rhodium at 1400 and 1600 K. Normal spectral emittance data is supplied at 156 wavelengths in the one to fifteen micron range for all the combinations listed above. In addition, data for the platinum-13 percent rhodium SRM's is supplied in the fifteen to thirty-five micron range at 1100 K.

SRM	Type	Unit Size
1402	Emittance standards	1/2 in disks Pt-13% Rh
1403	Emittance standards	7/8 in disks Pt-13% Rh
1404	Emittance standards	1 in disks Pt-13% Rh
1405	Emittance standards	1 1/8 in disks Pt-13% Rh
1406	Emittance standards	1 1/4 in disks Pt-13% Rh
1407	Emittance standards	2 in x 2 in squares Pt-13% Rh
1408	Emittance standards	1 in x 10 in strips Pt-13% Rh
1409	Emittance standards	3/4 in x 10 in strips Pt-13% Rh
1420	Emittance standards	1/2 in disks Kanthal
1421	Emittance standards	7/8 in disks Kanthal
1422	Emittance standards	1 in disks Kanthal
1423	Emittance standards	1 1/8 in disks Kanthal
1424	Emittance standards	1 1/4 in disks Kanthal
1425	Emittance standards	2 in x 2 in squares Kanthal
1427	Emittance standards	3/4 in x 10 in strips Kanthal
1428	Emittance standards	1/4 in x 8 in strips Kanthal
1440	Emittance standards	1/2 in disks Inconel
1441	Emittance standards	7/8 in disks Inconel
1442	Emittance standards	1 in disks Inconel
1443	Emittance standards	1 1/8 in disks Inconel
1444	Emittance standards	1 1/4 in disks Inconel
1445	Emittance standards	2 in x 2 in squares Inconel

Refractive Index Standards

These SRM's are certified for refractive index for each of seven wavelengths: helium 668 and 502, hydrogen 656(C) and 486(F), mercury 546(e) and 436(g), and sodium 589(D₁, D₂) at 20, 25, and 30 °C to ± 0.00002 . A certificate is supplied with each of these samples. 217b-8S is contained in a special ampoule with an internal breakoff tip, the others are sealed "in vacuum" in plain glass ampoules.

SRM	Type	Approx. n_D^{20}	Quantity (ml)
217b-5	2,2,4-Trimethylpentane	1.3915	5
217b-8S	2,2,4-Trimethylpentane	1.3915	8
217b-25	2,2,4-Trimethylpentane	1.3915	25
217b-50	2,2,4-Trimethylpentane	1.3915	50

These standards are also certified for density.

Reflectance Standards

These SRM's are intended primarily for calibration of (1) reflectometers used in the evaluation of the appearance properties of polished metals and metal plated objects, and (2) reflectometers and other equipment used in the evaluation of thermal radiation properties of materials. These properties are of particular importance in the automotive and aerospace industries, although they also have many other applications.

The SRM's are mirrors produced by vacuum deposition of gold on glass and aluminum on glass, and are calibrated for near-normal (9°) specular reflectance. They are certified in terms of absolute reflectance over the wavelength range from 0.25 to 30 micrometers.

Specular Spectral Reflectance Standards

SRM	Type	Blank Size (cm)	Coated Area Size (cm)
2001	Aluminum on Glass	7.6 × 10.2 × 1.9	5.2 × 7.6
2002	Aluminum on Glass	3.8 × 3.8 × 1.3	2.5 × 2.5
2003	Aluminum on Glass	Disk 2.9 dia. × 1.0	entire surface
2005	Gold on Glass	7.6 × 10.2 × 1.9	5.1 × 7.6
2006	Gold on Glass	3.8 × 3.8 × 1.3	2.5 × 2.5
2007	Gold on Glass	Disk 2.9 dia. × 1.0	entire surface
2008	Gold on Glass	Disk 2.4 dia. × 0.6	entire surface



Radioactivity Standards

Information concerning the SRM appears on it or its container. A Certificate containing pertinent information on the SRM is sent under separate cover; a photocopy of the certificate is sent with the SRM. Copies of these Certificates and information concerning the applications of these SRM's are available on request to the NBS Office of Standard Reference Materials. These Materials (except the carbon-14 contemporary dating standard) are shipped only by express or air freight (shipping charges collect). The prices of SRM's may change as current stocks are depleted and are replaced. Purchasers will be billed at the prices in effect at the time of shipment.

Certain Radioactivity SRM's are not economical to maintain in stock because of short half lives and low demand. When sufficient demand for them exists, based on letters of inquiry, these SRM's are renewed and those who have expressed an interest in them are notified of their availability.

If you need radioactivity standards other than those listed below, or additional technical information on the listed Radioactivity SRM's, please contact the Radioactivity Section, Room C114, Radiation Physics Building, National Bureau of Standards, Washington, D.C. 20234. (Telephone 301-921-2668). NOTE: Chemically stable solutions of most radionuclides, including those no longer issued by NBS or that are currently out of stock, may be submitted to NBS for calibration as described in "Calibration and Test Services of the National Bureau of Standards," NBS Special Publication 250. Requests for these tests should be submitted, with full source information for approval of suitability, to the Radioactivity Section.

The stated accuracies of the older standards are, in general, an estimate of the standard deviation added to an estimate of maximum possible systematic error. The accuracies of more recent standards are based on the 99 percent confidence level of precision, with the same estimate of systematic error.

The International Commission on Radiation Units (ICRU) recommended definition of the activity (A) of a quantity of a radioactive nuclide is the quotient of ΔN by Δt , where ΔN is the number of nuclear transformations that occur in this quantity, in time Δt : ($A = \Delta N / \Delta t$). NBS uses the abbreviation ntps for nuclear transformation per second. In this list both ntps and dps are used; the latter when dps has been used in certificates printed before 1968. The terms: α ps, β ps, β^+ ps, K-x-rays ps, γ ps are used for the emission rates of alpha particles, beta particles, positrons, K-x-rays, and gamma-rays, respectively.

The SRM's listed below, not marked with an asterisk (*), may be ordered singly, without a license, under the general licensing provisions of the Atomic Energy Act of 1954. Those marked by an asterisk are available only under the special licensing provisions of the Atomic Energy Act of 1954.

Alpha-Particle Standards

These SRM's consist of a practically weightless deposit of the nuclide on a thin platinum foil cemented to a monel disk.

SRM	Radionuclide	Approximate Activity at Time of Calibration (Month, Year)	Accuracy (%)
*4900	Polonium-210	100 α ps	
*4901	Polonium-210	200 α ps	
*4902	Polonium-210	500 α ps	
*4906	Plutonium-238	1.4×10^3 to 3.4×10^4 ntps (4/69)	± 0.8
4904-D	Americium-241	2×10^3 to 5×10^4 ntps (2/70)	1.0
4907	Gadolinium-148	IN PREP	

* Available on Request

Beta-Ray and Gamma-Ray Gas Standards

These SRM's contain the Radionuclide in the inactive gas at a pressure of about one atmosphere in a glass break-seal ampoule.

SRM	Radionuclide	Approximate Activity at Time of Calibration (Month, Year)	Accuracy (%)
4935-C	Krypton-85	7×10^7 ntps per mole (7/69)±	3†
*4235	Krypton-85	1 to 2×10^7 ntps (per source) (9/69)	3.8†
4236	Xenon-133	IN PREP	
4306	Xenon-133	IN PREP	
4307	Xenon-133	IN PREP	
4300	Argon-37	4×10^6 ntps per mole (5/72)	4.3
4301	Argon-37	4×10^3 ntps per mole (5/72)	5.8
4302	Argon-39	IN PREP	
4303	Argon-39	IN PREP	
4304	Xenon-131m	IN PREP	
4305	Xenon-131m	IN PREP	

†These are provisional accuracies. When the final accuracies are determined, previous purchasers will be notified.

Beta-Ray, Gamma-Ray, and Electron-Capture Solution Standards

These standard reference materials are contained in flame-sealed ampoules.

SRM	Radionuclide	Approximate Activity or Emission Rate per gram of Solution at Time of Calibration (Month, Year)	Approx. Weight of Solution (gram)	Accuracy (%)
4925	Carbon-14 (benzoic acid in toluene)	2×10^4 dps (7/58)	3	±2.0
4222	Carbon-14 (n-hexadecane)	4×10^4 dps (6/67)	3	3.1
4223	Carbon-14 (n-hexadecane)	4×10^3 dps (6/67)	3	3.1
4224	Carbon-14 (n-hexadecane)	4×10^2 dps (6/67)	3	3.1
4943	Chlorine-36	1×10^4 β ps (1962)	3	2.0
4941-C	Cobalt-57	3×10^5 ntps (3/69)	5.2	1.0
4926	Hydrogen-3 (Water)	9×10^3 dps (9/61)	25	1.0
4927	Hydrogen-3 (Water)	9×10^5 dps (9/61)	3	1.0
4947	Hydrogen-3 (Tritiated toluene) . . .	3×10^5 dps (2/64)	4	1.0
4929-C	Iron-55	2×10^4 K-x-rays ps (4/70)	3.9	2.7
*4226	Nickel-63	1.5×10^6 ntps (5/68)	4.1	1.0
4940-B	Promethium-147	5×10^4 dps (11/67)	3	1.9
*4228	Selenium-75	2.5×10^5 ntps (3/71)	4.6	2.3
4921-C	Sodium-22	1×10^4 β ⁺ ps (8/64)	2.8	1.0
4922-E	Sodium-22	2×10^5 β ⁺ ps (3/67)	5.1	1.4
4229	Aluminum-26	39 ntps (11/71)	4.6	1.1
4230	Chromium-51	IN PREP		
4232	Silver-110m	IN PREP		
4245	Carbon-14 (Na ₂ CO ₃ in H ₂ O)	10 μCi (1972)	5	
4246	Carbon-14 (Na ₂ CO ₃ in H ₂ O)	1 μCi (1972)	5	
4247	Carbon-14 (Na ₂ CO ₃ in H ₂ O)	0.01 μCi (1972)	5	
4949	Iodine-129	IN PREP		
4231	Cobalt-56	IN PREP		
4233	Cesium-137-Barium-137m	IN PREP		
4234	Barium-140-Lanthanum	IN PREP		

Contemporary Standard for Carbon-14 Dating Laboratories

SRM	Description
4990-B	Oxalic acid; no specific activity is given. (One pound of oxalic acid taken from specially prepared material for use as a common contemporary standard against which world-wide measurements can be compared.)

NOTE: These SRM's are shipped parcel post, prepaid to domestic and overseas purchasers.

Environmental Standards

Mixed-radionuclide gamma-ray emission-rate standards have been issued for monitoring radioactive effluents of nuclear-power reactors.

These SRM's contain: Cadmium-109, Cobalt-57, Tin-113-Indium-113m, Cesium-137-Barium-137, Manganese-54, Cobalt-60, and Yttrium-88.

SRM	Type	Approximate Activity (month/year)	Unit (ml)
4242-B	Mixed Radionuclides	0.5 μ Ci total (1/72)	450
4242-B	Mixed Radionuclides	1 μ Ci total (1/72)	50
4244-B	Mixed Radionuclides	10 ³ γ ps total (1/72)	15
4252	Mixed Radionuclides	Test Standard	450
4253	Mixed Radionuclides	Test Standard	50

Gamma-Ray "Point-Source" Standards

This group of Standard Reference Materials is usually prepared by depositing the radioactive material and sealing it between two layers of polyester tape, mounted on an aluminum ring. Exceptions to this procedure are americium, krypton, and thorium SRM's. The americium-241 SRM's, 4211 and 4213, are prepared by electroplating americium onto a 0.010-cm thick platinum foil, which is covered with a 0.005-cm thick aluminum foil. The aluminum-covered source is sandwiched between two layers of 0.036-cm thick polyurethane film tape. The krypton-85 SRM, 4212, is prepared by sealing a krypton-85 impregnated aluminum foil between two glass disks, with an epoxy adhesive. The thorium-228 SRM's, 4205 and 4206, are prepared by depositing and sealing the radionuclide between two layers of gold foil and this sandwich is then sealed between two double layers of polyurethane-film tape.

SRM	Radionuclide	Gamma-Ray Energy (MeV)	Approximate Activity (ntps) at Time of Calibration (month, year)	Accuracy (%)
4211	Americium-241	0.060	4.0 to 18×10^4 (2/70)	± 2.8
*4213	Americium-241	0.060	1.9 to 4.1×10^5 (2/70)	2.8
4202	Cadmium-109	0.088	2×10^6 (12/67)	1.8
*4212	Krypton-85	0.514	6.5 to 37×10^6 (5/71)	2.6
4200-B	Cesium-137	0.662	7×10^4 (12/68)	1.3
4207	Cesium-137	0.662	5×10^5 (12/68)	1.3
4201-B	Niobium-94	0.702, 0.871	4 to 6×10^3 (4/70)	1.5
*4203	Cobalt-60	IN PREP		
*4210	Cobalt-60	1.173, 1.332	2×10^6 (4/69)	1.1
4991-C	Sodium-22	1.274	6×10^4 (4/69)	1.5
4996-B	Sodium-22	1.274	3×10^5 (4/69)	1.5
4205	Thorium-228	2.614	6×10^4 (8/68)	2.2
4206	Thorium-228	2.614	6×10^5 (8/68)	2.2
4240	Bismuth-207	IN PREP		
4214	Cobalt-57	IN PREP		
4215	Mixed Radionuclides	IN PREP		
4216	Mixed Radionuclides	IN PREP		

Radium Gamma-Ray Solution Standards

These samples are contained in flame-sealed glass ampoules.

SRM	Radium Content (in micrograms)	Accuracy (%)
4955	0.1	±3.6
4956	0.2	4.4
4957	0.5	1.8
4958	1.0	1.8
4959	2.0	1.3
4960	5.0	1.3
4961	10	1.1
4962	20	1.1
4963	50	1.1
4964-B	102	0.5

Radium Solution Standards for Radon Analysis

These samples are contained in flame-sealed glass ampoules.

SRM	Approximate Radium Content (gram)	Approx. Wt. Soln. (grams)	Accuracy (%)
4951	10^{-11}	100	±0.3
4950-B	10^{-9}	20	1.0
4953	10^{-8}	20	1.0

Metallurgical Standards

These SRM's are intended for calibrating x-ray diffraction equipment to determine the relative amounts of austenite and iron carbide in steel.

SRM	Type	Form
485	Austenite in Ferrite (4% nominal)	Disk: 20.6 mm Diameter, 2.5 mm thick
493	Spheroidized Iron Carbide (Fe_3C) in Ferrite	Wafer: 29 × 29 × 2.4 mm

Mossbauer Standards

These SRM's are issued for the calibration of the isomer shift of iron compounds and alloys and to provide a uniform basis for presentation of mossbauer isomer shift data.

SRM	Type	Form
725	Sodium Pentacyanonitrosylferrite II (Sodium Nitroprusside) for Isomer Shift of Iron Compounds	Platelet: 1 × 1 × 0.0775 cm
1541	Iron Foil	Foil: 2.5 cm × 2.5 cm × 23 μm

Permittivity Standards

The three solution SRM's (1511, 1512, and 1513) are for calibrating cells and test capacitors used to determine the relative permittivity (dielectric constant) of liquids. The nominal dielectric constants (ϵ) for SRM's 1511, 1512, and 1513 are: 2.0, 10.4, and 35.7, respectively. The four polymer SRM's (1516, 1517, 1518, and 1519) are for calibrating systems used to measure permittivity and related dielectric quantities. These SRM's are disks of a fluorinated ethylene-propylene copolymer and are individually calibrated.

SRM	Type	Unit Size
1511	Cyclohexane	400 ml
1512	1,2-Dichloroethane	400 ml
1513	Nitrobenzene	400 ml
1516	Permittivity	38 mm diameter 2.5 mm thick
1517	Permittivity	38 mm diameter 5 mm thick
1518	Permittivity	51 mm diameter 2.5 mm thick
1519	Permittivity	51 mm diameter 5 mm thick



Engineering Type Standards

These SRM's are intended to relate measurements used for production or quality control data to a central point of reference. The values certified for these materials are in some cases empirical and do not necessarily relate to the National Measurement System.

Standard Rubbers and Rubber-Compounding Material

These SRM's have been prepared to provide the rubber industry with standard materials for rubber compounding. They are useful for the testing of rubber and rubber-compounding materials in connection with quality control of raw materials and for the standardization of rubber testing.

Each material has been statistically evaluated for uniformity by mixing rubber compounds and vulcanizing them in accordance with ASTM Designation D-15 and determining the stress-strain properties of the resulting vulcanizates. Certificates are issued for the rubbers because the properties of different lots are not the same. Replacement lots of rubber-compounding SRM's impart essentially the same characteristics to rubber vulcanizates so that Certificates are not issued for these SRM's.

Standard Rubbers

SRM	Type	Wt/Unit (grams)
385b	Natural	34,000
386g	Styrene-butadiene 1500	34,000
388f	Butyl	34,000
389	Styrene-butadiene 1503	34,000
391	Acrylonitrile-butadiene	25,000

Rubber Compounding Materials

SRM	Type	Wt/Unit (grams)
370d	Zinc Oxide	8,000
371f	Sulfur	6,000
372g	Stearic Acid	3,200
373f	Benzothiazyl disulfide	2,000
374c	Tetramethylthiuram disulfide	2,000
375f	Channel Black	28,000
376a	Light Magnesia	450
377	Phenyl-beta-naphthylamine	600
378a	Oil Furnace Black	28,000
379	Conducting Black	5,500
380	Calcium Carbonate	6,000
381	Calcium Silicate	4,000
382a	Gas Furnace Black	32,000
383	Mercaptobenzothiazole	3,200
384a	N-tertiary-Butyl-2-benzothiazolesulfenamide	3,200

Reference Magnetic Tapes

This SRM is intended for use in evaluating the performance of magnetic computer tapes and maintaining control over their production. Each SRM is individually calibrated and certified.

SRM	Type	Unit of Issue
3200	Secondary standard magnetic tape—computer amplitude reference	Reel/600 ft

Sizing Standards

Glass Spheres for Particle Size

SRM	Type	Size (μm)	Sieve Nos.	Wt/Unit (grams)
1003	Calibrated Glass Spheres	5-30	-----	40-45
1004	Calibrated Glass Beads	34-120	400-140	63
1017a	Calibrated Glass Beads	100-310	140-50	84
1018a	Calibrated Glass Beads	225-780	60-25	74
1019	Glass Spheres	890-2590	18-8	100

Turbidimetric and Fineness Standard (Cement)

This SRM is available to calibrate the Blaine fineness meter according to the latest issue of Federal Test Method Standard 158, Method 2101 or ASTM Designation C204; to calibrate the Wagner turbidimeter according to ASTM Designation C115; and to determine sieve residue according to ASTM Designation C430. Each set consists of twenty sealed vials, each containing approximately 10 grams of cement. This SRM is supplied only in sets of twenty vials or multiples thereof.

SRM	Type	Certification	Unit
114m	Portland Cement	Residue on No. 325 sieve, electroformed wet method Surface area (Wagner turbidimeter) Surface area (Air-permeability) Mean particle diameter (Air-permeability)	Set of 20 vials

Color Standards

The ISCC-NBS Centroid Color Charts

SRM 2106, ISCC-NBS Centroid Color Charts, is available to illustrate a characteristic color for each of the ISCC-NBS color-name blocks in NBS Circular 553. NBS Circular 553, The ISCC-NBS Method of Designating Colors and a Dictionary of Color Names, may be purchased from the Superintendent of Documents, U.S. Government Printing Office, Washington, D.C. 20402, for \$2 (SD Catalog No. C13.4:553). SRM 2106, along with the table containing the history of the color-names project, the centroid number, and the Munsell denotation of each of the 251 color chips included, constitutes a supplement to NBS Circular 553. Each chart set contains 18 constant-hue centroid color charts. These centroid colors represent a systematic sampling of the whole color solid, each color of which has been carefully measured. Each centroid color has its own specification and can be used as a color standard. The centroid color charts can also be used for approximate color specifications wherever the ISCC-NBS color designations are applicable, for statistical studies of trends in industrial color usage, or for planning lines of merchandise intended to have coordinated colors.

SRM	Type	Unit of Issue
2106	Centroid color charts	Set of 18 charts

Paint Pigment Standards for Color and Tinting Strength

SRM	Type	Wt/Unit (grams)
307	Metallic brown	60

Light-Sensitive Papers and Plastic Chips

Light-Sensitive Papers

Standard light-sensitive paper and booklets of standard faded strips of this paper are available for use in standardizing the dosage of radiant energy when testing textiles for color fastness by exposure in commercial carbon-arc fading lamps. The paper is distributed in units of 100 pieces 2 5/8 by 3 1/4 in. The booklets contain six strips of the paper 1 1/4 in wide that have been faded by exposure in the NBS master lamp. A copy of NBS Misc. Publ. 260-41, which describes the preparation and use of the materials, is furnished with each booklet.

SRM	Type	Unit of Issue
700c	Light-sensitive paper	Pkg. of 100 pieces - 2 5/8 in x 3 1/4 in
701c	Standard faded strips	Booklet - 6 strips 1 1/4 in wide

Light-Sensitive Plastic Chips

Standard light-sensitive plastic chips are available for use in calibration and standardization of artificial weathering and fading apparatus. These chips are distributed in two thicknesses (0.060 and 0.124 in) in units of five plates 2 in by 4 1/8 in, and have been standardized by the measurement of the change of transmittance as a function of exposure (in standard fading hours) to the NBS master lamps.

SRM	Type	Unit of Issue
702	Light-sensitive plastic chips . .	Package of 5 chips 0.124 in thick
703	Light-sensitive plastic chips . .	Package of 5 chips 0.060 in thick

Photographic Standards

SRM's 1008 and 1009 are calibrated photographic step tablets of 21 steps that cover the optical density range from 0 to 4 and 0 to 3, respectively.

These step tablets are designated as Type Visual VI-b and are certified for diffuse transmission density in conformance with conditions specified for American National Standard Diffuse Visual Density, Type VI-b "ANSI PH2.19-1959, American National Standard Diffuse Transmission Density." The 21 steps of the 3.5 X 25 cm tablets cover the density range from 0 to 4 (1008) or 0 to 3 (1009) and are individually certified by a method having a precision such that three times the standard deviation of the mean is 1 percent or 0.01, whichever is greater.

SRM 1010a, Microcopy Resolution Test Charts, is used to test the resolving power of cameras or of whole microcopying systems. SRM 1010a consists of five charts printed photographically on paper, and have 26 high-contrast five-line patterns ranging in spatial frequency from one cycle per millimeter to 18 cycles per millimeter. Instructions for the use of the charts are supplied with each order.

SRM	Type	Unit of Issue
1008	Photographic Step Tablet (0-4)	1 tablet, 21 steps
1009	Photographic Step Tablet (0-3)	1 tablet, 21 steps
1010a	Microcopy Resolution Test Charts	Set of 5 charts

Surface Flammability Standard

SRM 1002b, Hardboard Sheet, is issued for checking the operation of radiant panel test equipment in accordance with the procedures outlined in ASTM Standard E162-67. Flame Spread Index, $I_s = 190$; Heat Evolution Factor, $Q = 45.4$.

SRM	Type	Certification	Unit of Issue
1002b	Hardboard Sheet	Flame Spread Index Heat Evolution Factor, Q	190 45.4 Set of 4 6 X 18 X 1/4 inch

Smoke Density Chamber Standards

These SRM's are certified for maximum specific optical density and are issued for performing operational checks of smoke density chambers.

SRM	Type	Maximum Specific Optical Density	Unit of Issue
1006	Non-flaming Exposure Condition (α -cellulose)	Dm(corr) = 170	Set of 3 sheets 11.7 × 9.2 inches
1007	Flaming Exposure Condition (plastic)	Dm(corr) = 455	Set of 3 sheets 9.5 × 9.2 inches

Water Vapor Permeance

This material is intended for use in the measurement of water vapor permeance in accordance with ASTM Method D-96. It may also be useful in other test methods where movement of water vapor across a barrier is involved. These SRM's are made from sheets of poly (ethylene terephthalate) approximately 0.001 inches thick (25.4 μ m). They are certified for water vapor permeance for both dry cup and wet cup procedure.

SRM	Type	Certification	Unit of Issue
707-1	Water Vapor Permeance	Dry Cup - 0.66 perm	12 sheets, 6 in diameter
707-2	Water Vapor Permeance	Wet Cup - 0.72 perm	6 sheets, 10 × 12 inches

Internal Tearing Resistance Standard Paper

This SRM is available for calibration of instruments used for the determination of the internal tearing resistance of paper according to methods ASTM Designation D689 and TAPPI Standard T414. Sufficient material is furnished in each unit to provide 40 or more measurements. Initial distribution is in a set of twelve packages, one package shipped at approximately monthly intervals. Packages are also available on a four month cycle. The tearing strength value of the material is approximately 40 g. The exact value will be given in the certificate accompanying the standard.

This SRM is sold only on a subscription basis in sets of four packages or multiples thereof.

SRM	Type	Unit of Issue
704a	Internal tearing resistance paper	Sets of 4 packages

Linerboard Standard for Tape Adhesion Testing

This material is intended as a uniform source of linerboard for use under ASTM Designation D2860, Procedure A: Adhesion of Pressure Sensitive Tape to Fiberboard at 90 Degree Angle and Constant Stress.

SRM	Type	Unit
1810	Linerboard for Tape Adhesion Testing	Package of 50 sheets

Research Materials

Research Materials (RM's) are in addition to and distinct from the Standard Reference Materials (SRM's) issued by NBS. The distinctions between Research Materials and Standard Reference Materials are in the information supplied with them and purpose for which they are used. Unlike SRM's, the RM's are not issued with Certificates of Analysis; rather they are accompanied by a "Report of Investigation," the sole authority of which is the author of the report. A Research Material is intended primarily to further scientific or technical research on that particular material. One of the principal considerations in issuing an RM is to provide homogeneous material so that an investigator in one laboratory can be assured that the material he has is the same as that being investigated in a different laboratory.

High Purity Materials

RM-1C Ultra-purity aluminum single crystal cubes (1 cm on a side) are intended for use in studies of a variety of solid state phenomena for which both extreme purity and knowledge of crystallographic orientation are required; e.g., in studies of electron spin resonance, De Haas-Van Alphen effect, cyclotron resonance, and in a variety of studies relating to the Fermi surface and the transport properties of aluminum.

RM-1R Ultra-purity aluminum polycrystalline rods (4.2 mm in diameter and 25.4 mm long) are intended for use in research on the mechanical and physical properties of extremely pure aluminum: e.g., in the determination of resistivity as a function of strain at cryogenic temperatures to facilitate the design of cryogenic magnets or superconductor stabilizing elements.

Phosphors


These materials are issued without Certification. NBS Technical Note 417, Spectral Emission Properties of NBS Standard Phosphor Samples under Photo-Excitation, is issued with these materials, and is equivalent to the "Report of Investigation" issued with Research Materials. They are issued so that those interested in developing methods of measurement for phosphor materials can work on a common source of materials. NBS Technical Note 417 may be purchased from Superintendent of Documents, U.S. Government Printing Office, Washington, D.C. 20402, for 25 cents, by SD Catalog No. C13.46:47.

SRM	Type	Wt/Unit (grams)
1020	Zinc sulfide phosphor	14
1021	Zinc silicate phosphor	28
1022	Zinc sulfide phosphor	14
1023	Zinc-cadmium sulfide phosphor (Ag activator)	14
1024	Zinc-cadmium sulfide phosphor (Cu activator)	14
1025	Zinc phosphate phosphor	28
1026	Calcium tungstate phosphor	28
1027	Magnesium tungstate phosphor	28
1028	Zinc silicate phosphor	28
1029	Calcium silicate phosphor	14
1030	Magnesium arsenate phosphor	28
1031	Calcium halophosphate phosphor	28
1032	Barium silicate phosphor	28
1033	Calcium phosphate phosphor	28

General Materials

General Materials (GM's) are being distributed by NBS to meet industry needs. These materials have been standardized either by some Government agency other than NBS, or by some standards-making body such as the American Society for Testing and Materials (ASTM), the American National Standards Institute (ANSI), and the Organization for International Standardization (ISO). For this class of materials, NBS acts only as a distribution point and does not participate in the standardization of these materials.

- GM-1 Hydrogen in Steel Standards were produced and certified by The Welding Institute in Cambridge, England, and are distributed in the United States by NBS. GM-1 is a set of 15 cylinders, 5 each of H1, H2, and H3, containing nominally 0.05, 0.10, and 0.20 ml hydrogen, respectively. The cylinders are 6.35 mm in diameter and about 30 mm long, weighing approximately 6 grams.
- GM-2 Hydrogen in Steel Standards were produced and certified by the Welding Institute in Cambridge, England, and are distributed in the United States by NBS. GM-2 is a set of 15 cylinders, 5 each of H4, H5, and H6, containing nominally 0.20, 0.60, and 1.10 ml hydrogen, respectively. The cylinders are 12.7 mm in diameter and about 30 mm long, weighing approximately 22 grams.
- GM-5 Nickel and Vanadium in Residual Oil was produced and analyzed under the sponsorship of the Western Oil and Gas Association and the American Petroleum Institute, and is distributed by NBS. The assigned values for nickel and vanadium are 93 and 79 ppm, respectively. GM-5 is issued in 475 ml units.
- GM-2007 Attapulugus clay is distributed by NBS on request of the ASTM Committee D-2007. It is an adsorbant type clay, 30 to 60 mesh, having adsorptive characteristics as specified by ASTM D-2007.



ALPHABETICAL INDEX

PAGE

ACIDIMETRIC	42	AMERICIUM	
Benzoic acid		Radioactivity	68
Potassium phthalate, acid		AMMONIUM DIHYDROGEN PHOSPHATE	
ALLOYS CHEMICAL COMPOSITION		Fertilizer	48
(See also individual metals)		ANALYZED GASES	45
Aluminum		Carbon dioxide in nitrogen	
Chip form	31	Hydrocarbons in air	
Wire (neutron monitor density)	53	Oxygen in nitrogen	
Cobalt		ANALYZED LIQUIDS	
Chip form	31	Permeation tubes	45
Wire (neutron monitor density)	53	Nickel and vanadium in residual oil	78
Copper (see also copper)		Sulfur in distillate fuel oil	45
Chip form	31	Sulfur in residual fuel oil	45
Microprobe	41	ANALYZED SOLIDS	45
Solid form	34	Sulfur in coal	
Ferro (steelmaking)	30	Trace mercury in coal	
Gold		ARSENIC	
High-purity	40	Trioxide (oxidimetric value)	42
Microprobe wire	41	ASSAY STANDARDS	
High temperature		(See specific constituent)	
Solid form	26	BARIUM	
Iron		Metallo-organic compound	44
Chip and granular form	28	BASIMETRIC	
Microprobe	41	Tris(hydroxymethyl)aminomethane	42
Solid form	28	BAUXITE (aluminum ore)	47
Lead		BENZOIC ACID	
Chip form	32	Acidimetric	42
Solid form	36	Calorimetric	61
Magnesium		Microchemical	42
Molybdenum		BIOLOGICAL STANDARDS	43
Microprobe	41	Liver, bovine	
Nickel		Orchard leaves	
Chip form	32	Tomato leaves	
Oxides	32	Tuna, albacore	
Solid form	36	BISMUTH	
NONFERROUS ALLOYS		Radioactivity	68
Platinum	40	BORON	
Doped, wire		Boric acid	42,53
High purity, wire		Boric acid (B-10 enriched)	53
Selenium	33	Metallo-organic compound	44
Granular form		CADMIUM	
Silver	41	Metallo-organic compound	44
Microprobe		Radioactivity	68
Solder		Vapor pressure standard	62
Chip form	32	CALCIUM	
Solid form	36	Carbonate (clinical)	43
Tin		Carbonate (rubber compounding)	73
Chip form	33	Metallo-organic compound	44
Solid form	36	Molybdate (steelmaking)	30
Titanium		Silicate (rubber compounding)	73
Chip form	33	CALORIMETRIC	61
Solid form	36	Combustion calorimetric	
Gas-in	37	Heat source	
Tungsten	41	CANE SUGAR (Sucrose)	42
Microprobe		CARBIDES	
Spelter	37	Silicon (refractory)	50
Zinc		CARBON	
Chip form	33	Channel black (rubber compounding)	73
High-purity	40	Conducting black (rubber compounding)	73
Solid form	37	Gas furnace black (rubber compounding)	73
Zirconium		Oil furnace black (rubber compounding)	73
Chip form	33	Carbon dioxide in nitrogen	45
Solid form	38	CARBON STEELS	16
ALUMINUM		CAST IRONS	28
Al-Co alloy (neutron density monitor wire)	53	Blast furnace	
Al-Si alloy	31	Car wheel	
Casting alloy	31	Chip form	
Freezing point standard	60	Copper-nickel-chromium	
Metallo-organic compound	44	Ductile	
Ore (bauxite)	47	Gray	
Radioactivity	69	High phosphorus	
Wrought alloy	31	Nickel-chromium	
		Nickel-chromium-molybdenum*	

Nodular		ISCC-NBS centroid color charts	74
Solid form		Light-sensitive paper	75
White		Light-sensitive plastic chips	75
CAST STEEL	28	Tristimulus integrator systems - spectrophotometer	65
CEMENTS		COPPER	
Portland (chemical composition)	48	Beryllium copper alloys	34
Turbidimetric and fineness standard	74	Brass	
CERAMIC MATERIALS		Aluminum	34
Carbides	50	Cartridge	34
Glasses	50	Free Cutting	34
Minerals	48	Naval	34
Refractories	50	Red	34
CERIUM		Sheet	31
Radioactivity	68	Bronze	
CESIUM		Cast	31
Radioactivity	68	Commercial	34
CHEMICAL STANDARDS		Leaded tin	31
Primary	42	Silicon	31
Acid potassium phthalate		Coating thickness	55
Arsenic trioxide		Concentrate	47
Benzoic acid		Freezing point standard	60
Boric acid		Gilding metal	34
Potassium chloride		Isotopic reference standard	53
Potassium dichromate		Linear thermal expansion	63
Sodium oxalate		Metallo-organic compound	44
Sucrose		Microprobe (gold-copper wires)	41
Tri(hydroxymethyl)aminomethane		Mill heads	47
Uranium oxide		Mill tails	47
Intermediate purity		Nickel silver alloy (Cu58-Ni12-Zn29)	31
Selenium	33	DENSITY STANDARDS	
Zinc	40	Neutron monitor wire	53
CHLORINE		2,2,4-trimethylpentane	58
Isotopic reference standard	53	Smoke density, flaming	76
Radioactivity	68	Smoke density, non-flaming	76
CHROMIUM		DEXTROSE (Glucose)	42
Ferro (steelmaking)	30	DIELECTRIC	72
Metallo-organic compound	44	Cyclohexane	
Isotopic reference standard	53	1,2 dichloroethane	
Radioactivity	68	Nitrobenzene	
Refractories	50	DIFFERENTIAL THERMAL ANALYSIS	62
CLAYS		DUCTILE IRON	28
Atapulgu (ASTM D-2007)	78	ENVIRONMENTAL STANDARDS	45
Flint	48	Analyzed gases	
Plastic	48	Analyzed liquids	
CLINICAL LABORATORY STANDARDS	43	Analyzed solids	
Bilirubin		ELASTICITY	59
Calcium carbonate		EMITTANCE	66
Cholesterol		FELDSPARS	48
Creatinine		Soda	
Cortisol		Potash	
Glass filters for spectrophotometry		FERTILIZER	
Liquid filters for spectrophotometry		Ammonium dihydrogen phosphate	48
D-Glucose		Phosphate rock	47
D-Mannitol		Potassium nitrate	48
Lithium carbonate		FILTERS	65
Potassium chloride		Glass for spectrophotometry	
Sodium chloride		Liquid for spectrophotometry	
Tris(hydroxymethyl)aminomethane		Spectrophotometer - tristimulus integrator systems	
Tris(hydroxymethyl)aminomethane HCl		FLAMMABILITY STANDARDS	
Urea		Flame spread index (hardboard sheet)	75
Uric acid		Smoke density, flaming	76
VMA		Smoke density, non-flaming	76
COATING THICKNESS	55	FLUORSPAR	47
Copper on steel		Assay	
Gold on copper		Geological	
Gold on copper-clad glass-epoxy laminate		FREEZING POINT STANDARDS	60
Gold on Fe-Ni-Co glass sealing alloy		Defining fixed points	
Gold on nickel		Tin	
Nickel on brass		Zinc	
Nickel on steel		Determined reference points	
Tin on steel		Aluminum	
COBALT		Copper	
Metallo-organic compound	44	Lead	
Radioactivity	68	Tin	
See Alloys		Zinc	
COLOR STANDARDS		GASES, ANALYZED	45
Paint pigment	74	Carbon dioxide in nitrogen	
		Hydrocarbons in air	


Oxygen in nitrogen	
Sulfur dioxide permeation tubes	
GASES, IN METALS	38
Iron alloys	
Steels	
Titanium alloys	
Zirconium alloys	
GENERAL MATERIALS	78
Clay, attapulugus (ASTM D-2007)	
Hydrogen in steel	
Nickel and vanadium in residual oil	
GEOLOGICAL	51
Fluorspar	47
GLASS	
Boron, high	50
Boron, low	50
Filters for spectrophotometry	65
Chemical composition	50
Lead-barium	50
Opal	50
Spheres for sieve calibration (sizing standards)	74
Trace elements in	51
Viscosity standards	58
GLUCOSE (Dextrose)	
Clinical	43
Primary chemical	42
GOLD	
Coating thickness	57
on epoxy	
on Fe-Ni-Co glass sealing alloy	
on nickel	
High-purity	40
Microprobe standards (gold-copper wire)	41
Microprobe standards (gold-silver wire)	41
Vapor pressure standard	62
HEAT SOURCE	61
Zirconium-barium chromate	
HIGH-PURITY METALS	40
Gold, wire and rod	
Platinum, doped	
Platinum, wire	
Zinc	
HIGH TEMPERATURE ALLOYS	26
HYDROCARBON BLENDS	46
HYDROGEN	
in steel (GM-1, GM-2)	78
in titanium	37
Radioactivity standards	68
INTEGRATOR SYSTEMS (COLOR FILTERS)	65
ION-ACTIVITY	54
Ion-selective electrodes	
pD standards	
pH standards	
IRON	
Chip form	28
Solid form	28
Metallo-organic compound	44
Ore	47
Radioactivity	68
Thermal conductivity	63
IODINE	
Radioactivity	68
ISOTOPIC REFERENCE STANDARDS	53
Boron, natural	
Boron, enriched	
Bromine	
Chlorine	
Chromium	
Copper	
Lead, equal atom, (206/208)	
Lead, natural	
Lead, radiogenic, 92% lead-206	
Magnesium	
Rubidium	
Silver	
Strontium	

KRYPTON	
Radioactivity	68
LEAD	
Alloys	32, 36
Bearing metal	32, 36
Freezing point	60
Isotopic reference standards	53
Metallo-organic compound	44
Solder	32, 36
LIGHT SENSITIVE STANDARDS	75
Light sensitive papers	
Paper, standard faded strips	
Light sensitive plastic chips	
LIMESTONE	48
Argillaceous	
Dolomitic	
LINERBOARD STANDARD FOR TAPE ADHESION	
TESTING	76
LIQUIDS, ANALYZED	
Nickel and vanadium in residual oil	78
Sulfur in residual fuel oil	45
Sulfur in distillate fuel oil	45
LITHIUM	
Metallo-organic compound	44
Ores	47
LOW ALLOY STEELS	
Chip form	16
Solid form	21
MAGNESIUM	
Alloy	32
Isotopic reference standard	53
Magnesia, light (MgO) (rubber compounding)	73
Metallo-organic compound	44
Refractory (burned magnesite)	50
MAGNETIC	
Susceptibility standard	64
Tape, computer amplitude reference - secondary standard	73
MANGANESE	
Metallo-organic compound	44
Ore	47
MARAGING STEELS	26
MERCURY	
Metallo-organic compound	44
Trace element in coal	46
MELTING POINTS	
Alumina	61
METALLO-ORGANIC COMPOUNDS	44
Aluminum	
Barium	
Boron	
Cadmium	
Calcium	
Chromium	
Cobalt	
Copper	
Iron	
Lead	
Lithium	
Magnesium	
Manganese	
Mercury	
Nickel	
Phosphorus	
Potassium	
Silicon	
Silver	
Sodium	
Strontium	
Tin	
Vanadium	
Zinc	
METALS	
Alloys - see index entry	
Elements - see index entry	
High-purity	40
Gas in	38

Microprobe	41	ORES	47
Vapor pressure	62	Bauxite (aluminum)	
Thermal conductivity	63	Iron, sibley	
Freezing points	60	Lithium	
METALLURGICAL STANDARDS	71	Lepidolite	
Austenite in ferrite		Petalite	
Iron carbide in ferrite		Spodumene	
MICROCHEMICAL STANDARDS	42	Manganese	
Acetanilide		Phosphate rock	
Anisic acid		Zinc concentrate	
Benzoic acid		OXALIC ACID	
o-Bromobenzoic acid		Radioactivity	68
m-chlorobenzoic acid		OXIDIMETRIC	42
Cystine		Arsenic trioxide	
p-fluorobenzoic acid		Potassium dichromate	
Nicotinic acid		Sodium oxalate	
Triphenyl phosphate		OXIDES	
Phosphate		Magnesium (rubber compounding)	73
Urea		Nickel	32
MICROCOPY RESOLUTION TEST CHARTS	75	Titanium	48
MICROPROBE STANDARDS	41	Uranium	42
Gold-copper		OXYGEN IN	38
Gold-silver		Ingot iron	
Iron - 3% silicon		Nitrogen	
Tungsten-20% molybdenum		Steel	
MINERALS	48	Titanium	
Clays		Zirconium	
Flint		PAINT	74
Plastic		PAPER	
Feldspars		Internal tearing resistance standard	76
Potash		Light sensitive	75
Soda		Faded strips	75
Limestone		PERMEANCE, WATER VAPOR	76
Argillaceous		PERMEATION TUBE	
Dolomitic		Sulfur dioxide	45
Ores	47	PERMITTIVITY	72
Bauxite (aluminum)		pD STANDARDS	54
Fluorspar		Disodium hydrogen phosphate	
Iron, sibley		Potassium dihydrogen phosphate	
Lithium		Sodium bicarbonate	
Manganese		Sodium carbonate	
Phosphate rock		pH STANDARDS	54
Zinc		Acid potassium phthalate	
MOLECULAR WEIGHT STANDARDS	59	Borax	
Polyethylene, branched		Disodium hydrogen phosphate	
Polyethylene, linear		Potassium dihydrogen phosphate	
Polystyrene, narrow molecular weight distribution		Potassium hydrogen tartrate	
Polystyrene, broad molecular weight distribution		Potassium tetroxalate	
MOLYBDENUM	41	Sodium bicarbonate	
Microprobe	41	Sodium carbonate	
MOSSBAUER STANDARDS	71	Tris(hydroxymethyl)aminomethane	
NEUTRON DENSITY MONITOR WIRE		Tris(hydroxymethyl)aminomethane HCl	
Cobalt in aluminum	53	PHOSPHATE	
NICKEL		Ammonium dihydrogen	48
Alloys		Rock	47
Chip form	32	PHOSPHORS	77
Solid form	36	Barium silicate	
Coating thickness	55	Calcium halophosphate	
on brass		Calcium phosphuate	
on steel		Calcium silicate	
Metallo-organic compound	44	Calcium tungstate	
Oxides	32	Magnesium arsenate	
Radioactivity	68	Magnesium tungstate	
NIOBIUM		Metallo-organic compound	
Radioactivity	68	Zinc-cadmium sulfide (Ag activator)	
NITROGEN IN	38	Zinc-cadmium sulfide (Cu activator)	
Cast iron		Zinc phosphate	
Ingot iron		Zinc silicate	
Steel		Zinc sulfide	
Titanium		PHOSPHORUS	
Zirconium		Ferro (steelmaking alloy)	30
NONFERROUS ALLOYS	31	Metallo-organic compound	44
NUCLEAR MATERIALS	52	PHOTOGRAPHIC	75
Neutron density monitor wire		Microcopy resolution test charts	
Plutonium assay		Step tablets	
Plutonium isotopic			
Uranium assay			
Uranium isotopic			

PLATINUM	40	SILVER	
Doped, wire		Alloys (microprobe)	41
High-purity, wire		Isotopic reference standard	53
PLUTONIUM	52	Metallo-organic compound	44
Metal assay		Radioactivity	68
Standard matrix		SIZING STANDARDS	74
Sulfate tetrahydrate		Cement	
Radioactivity		Calibrated glass spheres	
POLONIUM		Turbidimetric and fineness (cement)	
Radioactivity standard	68	SMOKE DENSITY CHAMBER STANDARDS	76
POLYMER		SODIUM	
Molecular weight	59	Bicarbonate (pD and pH standard)	54
Permittivity	72	Bromide (isotopic reference standard)	53
POLYETHYLENE,	59	Carbonate (pD and pH standard)	54
Branched		Dihydrogen phosphate (pD and pH standard)	54
Linear		Isotopic reference standard	53
POLYSTYRENE, MOLECULAR WEIGHT	59	Metallo-organic compound	44
Broad distribution		Oxalate (oxidimetric value)	42
Narrow distribution		Radioactivity	68
POTASSIUM		Tetaborate (borax) (pH standard)	54
Acid phthalate (pH standard)	54	SOLDER	
Dichromate (oxidimetric value)	42	Chip form	32
Dihydrogen thosphate (pD and pH standards)	54	Solid form	36
Hydrogen tartrate (pH standard)	54	SOLIDS, ANALYZED	45
Metallo-organic compound	44	SPECTROPHOTOMETRIC STANDARDS	65
Nitrate (fertilizer)	48	Color filters	
Tetroxalate (pH standard)	54	Glass filters	
		Liquid filters	
PRIMARY CHEMICALS	42	SPECULAR SPECTRAL REFLECTANCE	67
Arsenic trioxide		Aluminum on glass	
Benzoic acid		Gold on glass	
Boron		STAINLESS STEELS	18, 24
Dextrose		STEELMAKING ALLOYS	30
Potassium chloride		Calcium molybdate	
Potassium dichromate		Ferrochromium (high carbon)	
Potassium acid phthalate		Ferrochromium (low carbon)	
Sodium oxalate		Ferroniobium	
Sucrose		Ferrophosphorus	
Uranium oxide		Ferrosilicon (50%)	
PROMETHIUM		Silicon, refined	
Radioactivity	68	STEEL, THERMAL CONDUCTIVITY	63
RADIOACTIVITY STANDARDS	68	STEELS, CAST	28
Alpha-particle standards		STEELS, CHIP FORM	
Beta-particle and gamma-ray gas standards		High alloy	18
Beta-particle, gamma-ray, and electron-capture		Low alloy	16
solution standards		Low alloy, special	16
Contemporary standard for Carbon-14 dating		Plain carbon	16
laboratories		Stainless	18
Environmental standards		Tool	18
Gamma-ray "point-source" standard		STEELS, GAS-IN	38
Radium gamma-ray solution standards		STEELS, GRANULAR FORM	18
Radium solution standards for radon analysis		STEELS, SOLID FORM	
RADIUM		Cast	28
Radioactivity	68	High temperature alloy	26
RADON		Low alloy	21
Radioactivity	68	Low alloy, special	22
REFLECTANCE, SPECULAR SPECTRAL	67	Maraging	26
REFRACTIVE INDEX		Oxygen	27
2,2,4-trimethylpentane	66	Stainless	24
REFRACTORY MATERIALS	50	Tool	26
RESEARCH MATERIALS	70	STRONTIUM	
Aluminum ultra-purity		Isotopic reference	53
Phosphors		Metallo-organic compound	44
RUBBER	73	SUCROSE	42
RUBBER COMPOUNDING STANDARDS	73	Primary chemical	
RUBIDIUM		SULFUR	
Isotopic reference standard	53	Elemental (rubber compounding)	73
SELENIUM		in residual fuel oils	44
Metal	33	in coal	45
Radioactivity	68	SUPERCONDUCTIVE THERMOMETRIC FIXED POINT	
SILICON		DEVICE	60
Ferro (steelmaking alloy)	30	SURFACE FLAMMABILITY STANDARD	75
		TAPE, MAGNETIC - SECONDARY REFERENCE	73
Metallo-organic compound	44	TEMPERATURE	
Refractories (SiO ₂)	50	Calorimetry	61
		Differential thermal analysis	62
		Heat source	61

THERMAL CONDUCTIVITY	63	TUNGSTEN	
THERMAL EMITTANCE STANDARDS	66	Microprobe	41
THERMAL EXPANSION	63	TURBIDIMETRIC AND FINENESS STANDARD ...	74
THERMOCOUPLE MATERIALS	63	URANIUM	52
THORIUM		Metal	
Radioactivity	68	Oxide	
TIN		Depleted	
Bearing metal	33	Enriched	
Chip form	33	Primary standard	42
Coating thickness	55	VANADIUM	
Freezing point	60	Metallo-organic compound	44
		VAPOR PERMEANCE, WATER	76
Metallo-organic compound	44	VAPOR PRESSURE STANDARD	62
Solid form	36	VISCOSITY STANDARDS	
TITANIUM		Glass	58
Chip form	33	WATER VAPOR PERMEANCE	76
Dioxide	48		
Gases in	38	Radioactivity	
Metallo-organic compound	44	ZINC	
		Freezing point	60
Solid form	36	High purity	40
Unalloyed for oxygen	37	Intermediate purity	40
TRACE ELEMENTS		Metallo-organic compound	44
Biological	43	Oxide (rubber compounding)	73
in glass	51	Spelter modified	37
Mercury in coal	46	ZIRCONIUM	
		Metal	38
		Zircaloy-2	33



NUMERICAL INDEX TO SRM CERTIFICATES

The Certificate Date listed is the current version of the Certificate. Those dates followed by the letter "P" indicate that it is a Provisional Certificate.

In general, Provisional Certificates are issued for Standard Reference Materials before all of the values have been certified, but after a sufficient number of values are certified so that the material is a valuable standard for the intended purpose. As additional values are certified, the Provisional Certificate may be revised and when and if all of the values are certified, the final Certificate is issued.

When new or revised Certificates are issued for SRM's they are announced in the SRM Price and Availability List which supplements this catalog. SRM purchasers whose Certificate shows an earlier date than listed below may obtain the current version of the Certificate from the Office of Standard Reference Materials, Room B314, Chemistry Building, National Bureau of Standards, Washington, D.C. 20234.

One or more of the following letters apply to materials where no date is listed.

- A. Individually Certified
- B. The material is issued with "Instructions for Use" in lieu of a certificate
- C. This material is not certified, refer to page reference for details
- D. Material is in preparation
- E. Research Material: Issued with a "Report of Investigation"
- F. General Material: Information provided, but not certified by NBS
- G. Set of SRM's: Issued with Certificates for the individual SRM's
- H. Prepared on Request: Refer to page reference for details

SRM	Certificate Date	Page	SRM	Certificate Date	Page	SRM	Certificate Date	Page
1b	11-8-66P	48	54d	9-20-57	33	103a	9-28-62	50
3b	3-11-69	28	55e	2-2-61	28	104	7-10-33	50
4j	1-25-71	28	57	D	30	105	4-22-37P	16
5L	11-9-70	28	58a	D	30	106b	3-24-61	16
6g	11-9-70	28	59a	11-6-69	30	107b	11-1-62	28
7g	10-5-59	28	64b	1-8-62	30	112	5-26-37	50
8j	4-10-72	16	65d	10-30-58	16	113a	D	47
10g	10-14-65	16	69a	12-28-51P	47	114m	D	74
11h	9-1-66P	16	70a	3-26-65P	48	115a	4-19-62	28
12h	3-7-66	16	71	2-1-29	30	120b	9-1-72	47
13g	11-18-65P	16	72f	7-20-65	16	121d	7-7-71P	18
14e	11-9-65	16	73c	7-13-66	18	122e	11-6-69	28
15g	10-1-65P	16	76a	D	50	123c	7-7-71P	18
16e	3-23-66P	16	77a	D	50	124d	11-6-70	31
17	4-63	42	78a	D	50	125b	4-10-70	16
19g	9-30-64	16	79a	12-6-71	47	126c	12-6-72	18
20g	10-23-70	16	81a	D	48	127b	1-25-68P	32
25c	4-25-66P	47	82b	4-12-66	28	131b	7-16-69	18
27e	6-21-66	47	83c	4-16-70P	42	132b	11-7-72P	18
30f	12-23-66P	16	84h	7-9-69	42	133a	6-8-56	18
32e	4-5-57	16	85b	5-9-57	31	134a	5-6-57	18
33d	6-10-55	16	87a	11-17-69P	31	136c	3-24-70P	42
36b	7-18-69	16	88a	1-31-67P	48	139a	12-17-59	16
37e	8-28-58	31	89	8-27-32	50	140b	4-16-70	42
39i	7-15-68	61	90	10-1-28	30	141b	11-23-65	42
40h	4-24-69P	42	91	6-15-31	50	142	7-9-69	42
41a	12-26-67	42	92	4-30-71P	50	143b	7-9-69	42
42g	7-18-72	60	93a	9-29-72	50	147	7-9-69	42
43h	D	60	94b	7-28-55	33	148	12-22-70	42
44f	1-5-73	60	97a	10-8-69	48	152a	10-11-65	16
45d	12-6-71	60	98a	10-8-69	48	153a	1-5-60	18
49e	12-6-71	60	99a	3-26-65P	48	154b	11-15-72	48
50c	6-25-57	18	100b	8-18-59	16	155	10-1-46	16
51b	10-20-57	16	101f	5-19-70	18	157a	10-16-58	31
53e	1-20-70	32				158a	8-8-61	31

SRM	Certificate Date	Page	SRM	Certificate Date	Page	SRM	Certificate Date	Page
160b	8-4-69P	18	371f	C	73	493	8-28-70	71
162a	12-18-58	32	372g	C	73	592	4-1-61	46
163	1-11-68	18	373f	C	73	593	4-1-61	46
165a	D	48	374c	C	73	594	4-1-61	46
166c	3-17-70	18	375f	C	73	595	4-1-61	46
168	1-12-56P	31	376a	C	73	596	4-1-61	46
171	11-27-51	32	377	C	73	597	4-1-61	46
173a	4-25-66P	33	378a	C	73	598	4-1-61	46
174	8-20-70P	33	379	C	73	599	4-1-61	46
176	8-20-70P	33	380	C	73	607	D	51
178	7-14-69	16	381	C	73	608	G	51
180	3-31-71	47	382a	C	73	609	G	51
181	8-20-70P	47	383	C	73	610	8-8-72P	51
182	8-20-70P	47	384a	C	73	611	8-8-72P	51
183	8-20-70P	47	385b	8-15-67	73	612	8-8-72P	51
184	4-10-57	31	386g	10-10-69	73	613	8-8-72P	51
185e	12-5-72	54	388f	4-7-72	73	614	8-8-72P	51
186lc	9-1-70	54	389	6-2-66	73	615	8-8-72P	51
186llc	9-1-70	54	391	2-14-69	73	616	8-8-72P	51
187b	9-9-70	54	404a	5-5-65	21	617	8-8-72P	51
188	1-10-64	54	405a	5-5-65	21	618	G	51
189	1-10-64	54	407a	5-5-65	21	619	G	51
191	7-15-68	54	408a	5-5-65	21	620	5-17-72P	50
192	7-15-68	54	409b	5-5-65	21	625	4-24-64	37
193	10-16-72	48	413	5-5-65	21	626	4-24-64	37
194	10-16-72	48	414	5-5-65	21	627	4-24-64	37
195	D	30	417a	5-5-65	21	628	4-24-64	37
196	11-9-70	30	418	5-5-65	21	629	4-24-64	37
198	1-6-60	50	418a	5-5-65	21	630	4-24-64	37
199	1-6-60	50	420a	5-5-65	21	631	4-16-70P	37
217b-5	9-28-62	58,61,66	427	5-5-65	21	633	D	48
217b-8S	9-28-62	58,61,66	431	6-23-50P	36	634	D	48
217b-25	9-28-62	58,61,66	432	6-23-50P	36	635	D	48
217b-50	9-28-62	58,61,66	433	6-23-50P	36	636	D	48
307	11-3-44	74	434	6-23-50P	36	637	D	48
330	D	47	435	6-23-50P	36	638	D	48
331	D	47	436	9-22-59	26	639	D	48
332	D	47	437	9-22-59	26	641	1-26-60P	36
333	D	47	438	9-22-59	26	642	1-26-60P	36
335	4-7-66	16	439	9-22-59	26	643	1-26-60P	36
337	8-16-66	16	440	9-22-59	26	644	1-26-60P	36
339	7-21-65	18	441	9-22-59	26	645	1-26-60P	36
340	11-9-70	30	442	10-25-65	24	646	1-26-60P	36
341	3-26-62	28	443	10-25-65	24	654a	3-22-71P	36
342	1-5-62	28	444	10-25-65	24	661	8-15-72P	22
342a	4-27-70	28	445	1-19-66	24	662	8-15-72P	22
343	5-28-62	18	446	1-19-66	24	663	8-15-72P	22
344	10-2-63	18	447	1-19-66	24	664	8-15-72P	22
345	1-23-64	18	448	1-19-66	24	665	8-15-72	22
346	7-14-65	18	449	1-19-66	24	666	G	22
348	1-17-66P	18	450	1-19-66	24	667	G	22
349	5-15-59P	33	461	5-5-65P	22	668	G	22
350	4-15-58	42	462	5-5-65P	22	671	9-12-60	32
352	8-7-61	37	463	5-5-65P	22	672	9-12-60	32
353	8-7-61	37	464	5-5-65P	22	673	9-12-60	32
354	8-7-61	37	465	5-5-65P	22	680L-1	12-28-67P	40
355	9-8-66	37	466	5-5-65P	22	680L-2	12-28-67P	40
356	9-8-66	37	467	5-5-65P	22	681L-1	12-28-67P	40
360a	1-30-67	33	468	5-5-65P	22	681L-2	12-28-67P	40
361	7-26-70P	16	479	9-29-72	41	682	2-1-71P	40
362	7-27-70P	16	480	11-22-68P	41	683	7-9-68P	40
363	9-30-70P	16	481	2-14-69	41	685-R	9-26-68P	40
364	6-16-71P	16	482	6-6-69	41	685-W	9-26-68P	40
365	7-28-70P	28	483	6-11-71	41	700c	B	75
366	G	16	485	7-30-70P	71	701c	10-1-71	75
370d	C	73						

SRM	Certificate Date	Page	SRM	Certificate Date	Page	SRM	Certificate Date	Page
702	11-1-66	75	767	D	60	953-L3	3-12-69	53
703	11-1-66	75	803a	5-5-65	21	960	5-12-72	52
704a	Undated	76	D803a	5-5-65	21	975	3-11-65	53
705	1-17-72	59	804a	5-5-65	21	976	3-24-65	53
706	Undated(1963)	59	805a	5-5-65	21	977	3-24-65	53
707	1-17-72	76	D805a	5-5-65	21	978	3-24-65	53
709	D	58	807a	5-5-65	21	979	5-3-66	53
710	6-29-62	58	D807a	5-5-65	21	980	1-31-67	53
711	7-1-64	58	808a	5-5-65	21	981-3	6-1-68	53
712	10-4-65	58	809b	5-5-65	21	984	7-27-70	53
713	10-4-65	58	D809b	5-5-65	21	987	3-6-72P	53
714	10-4-65	58	810a	5-5-65	21	988	10-4-72	53
715	9-7-66	58	817	5-5-65	21	999	9-6-72	42
716	9-7-66	58	820a	5-5-65	21	1002b	12-71	75
717	11-18-69	58	D820a	5-5-65	21	1003	7-1-65P	74
718	4-28-72	54	821	5-5-65	21	1004	4-3-72	74
720	8-26-70	62	827	5-5-65	21	1006	10-4-71	76
723	7-14-70	42	D836	9-22-59	26	1007	10-4-71	76
724	4-14-72	61	837	9-22-59	26	1008	12-5-72	75
725	3-1-71P	71	D837	9-22-59	26	1009	12-5-72	75
726	1-31-67	33	838	9-22-59	26	1010a	B	75
728	7-9-68	33	D838	9-22-59	26	1011	6-18-71P	48
731-L1	7-31-72	63	839	9-22-59	26	1013	3-16-70P	48
731-L2	7-31-72	63	D839	9-22-59	26	1014	6-18-71P	48
731-L3	7-31-72	63	840	9-22-59	26	1015	6-18-71P	48
733	12-30-71	63	D840	9-22-59	26	1016	1-17-69P	48
734-S	6-14-71	63	841	9-22-59	26	1017a	9-24-71	74
734-L1	6-14-71	63	D841	9-22-59	26	1018a	D	74
734-L2	6-14-71	63	845	1-19-66	24	1019	1-17-66	74
735-S	3-24-72	63	D845	1-19-66	24	1020	E	77
735-M1	3-24-72	63	D846	1-19-66	24	1021	E	77
735-M2	3-24-72	63	D847	1-19-66	24	1022	E	77
735-L1	3-24-72	63	849	1-19-66	24	1023	E	77
735-L2	3-24-72	63	D849	1-19-66	24	1024	E	77
736-L1	11-19-66	63	850	1-19-66	24	1025	E	77
736-L2	11-19-66	63	D850	1-19-66	24	1026	E	77
736-L3	11-19-66	63	911	10-20-67P	43	1027	E	77
737-L1	D	63	912	9-24-68	43	1028	E	77
737-L2	D	63	913	9-24-68	43	1029	E	77
737-L3	D	63	914	9-24-68	43	1030	E	77
739-L1	5-12-71	63	915	3-4-69	43	1031	E	77
739-L2	5-12-71	63	916	3-10-71P	43	1032	E	77
739-L3	5-12-71	63	917	11-18-70	43	1033	E	77
740	2-14-70	60	918	1-22-71	43	1051b	7-15-68	44
741	7-18-72	60	919	9-27-72	43	1052b	3-1-68	44
742	6-5-70	61	920	1-17-72	43	1053a	1-23-70	44
745	5-14-69	62	921	D	43	1055b	7-23-68	44
746	8-11-70	62	922	6-23-71P	43	1057b	8-5-68	44
747	D	62	923	6-23-71P	43	1059b	7-15-67	44
748	8-10-70	62	924	2-23-72P	43	1060a	4-24-64	44
749	D	62	925	D	43	1061c	9-9-70P	44
755	11-6-69P	62	930a	12-4-72	43, 54	1062a	4-24-64	44
756	11-6-69P	62	931	1-17-72	43, 54	1063a	4-24-64	44
758	6-11-71	62	944	2-14-69	52	1064	4-24-64	44
759	6-11-71	62	945	9-6-72P	52	1065b	11-1-67	44
760	6-11-71	62	946	12-3-71P	52	1066a	4-23-69	44
763-1	D	64	947	12-3-71P	52	1069b	2-13-69	44
763-2	D	64	948	9-1-72P	52	1070a	4-24-64	44
763-3	D	64	949d	9-1-72P	52	1071a	4-24-64	44
764-1	D	64	950a	10-27-69	52	1073b	7-15-67	44
764-2	D	64	951	2-28-69	53	1074a	5-13-66	44
765-1	D	64	952	2-28-69	53	1075a	10-25-67	44
765-2	D	64	953	3-12-69	53	1076	4-24-64	44
765-3	D	64	953-L1	3-12-69	53	1077a	2-10-68	44
766-1	D	64	953-L2	3-12-69	53			

SRM	Certificate Date	Page	SRM	Certificate Date	Page	SRM	Certificate Date	Page
1078a	7-25-72	44	1139	7-13-70	28	1343	A	56
1079b	2-26-69	44	1140	2-24-70	28	1344	A	56
1080	2-26-69	44	1141	2-24-70	28	1345	A	56
1089	G	27	1142	2-24-70	28	1346	A	56
1090	4-2-69	27	1143	12-22-69	28	1351	A	56
1091	4-2-69	27	1144	12-22-69	28	1352	A	56
1092	4-2-69	27	1147	12-22-69	28	1353	A	56
1093	3-4-69P	27	1148	12-22-69	28	1361	A	56
1094	6-12-69P	27	1149	12-22-69	28	1362	A	56
1095	10-23-70	27	1152	2-9-66	24	1363	A	56
1096	6-13-72	27	1154	2-9-66	24	1364	A	56
1097	5-26-72P	27	1155	8-4-69P	24	1365	A	56
1098	D	27	1156	4-5-66P	26	1366	A	56
1099	7-28-70P	27	1159	2-18-69P	36	1367	A	56
1101	7-31-64P	34	1160	2-18-69P	36	1368	A	56
C1101	7-31-64P	34	1165	5-5-65P	22	1369	A	56
1102	7-31-64P	34	1166	5-5-65P	22	1370	A	56
C1102	7-31-64P	34	1167	5-5-65P	22	1371	A	57
1103	8-13-65	34	1171	7-7-71P	24	1372	A	57
C1103	8-13-65	34	1172	7-7-71P	24	1373	A	57
1104	8-13-65	34	1185	2-4-59P	24	1374	A	57
C1104	8-13-65	34	1206-2	11-6-70P	26	1375	A	57
1105	7-14-65P	34	1207-1	11-6-70P	26	1376	A	57
C1105	8-13-65	34	1207-2	11-6-70P	26	1377	A	57
1106	11-17-69P	34	1208-1	11-6-70P	26	1378	A	57
C1106	11-17-69P	34	1208-2	11-6-70P	26	1381	A	57
1107	11-17-69P	34	1209	G	26	1382	A	57
C1107	11-17-69P	34	1210	11-28-62P	38	1383	A	57
1108	11-17-69P	34	1261	8-16-72P	22	1384	A	57
C1108	11-17-69P	34	1262	8-16-72P	22	1385	A	57
1109	7-14-65P	34	1263	8-16-72P	22	1386	A	57
C1109	8-18-70P	34	1264	8-16-72P	22	1398	A	57
1110	8-18-70P	34	1265	8-16-72P	22	1399	A	57
C1110	8-18-70P	34	1266	G	22	1402	4-11-66	66
1111	8-18-70P	34	1301	A	56	1403	4-11-66	66
C1111	8-18-70P	34	1302	A	56	1404	4-11-66	66
1112	1-17-69P	34	1303	A	56	1405	4-11-66	66
C1112	1-17-69P	34	1304	A	56	1406	4-11-66	66
1113	1-17-69P	34	1305	A	56	1407	4-11-66	66
C1113	1-17-69P	34	1306	A	56	1408	4-11-66	66
1114	1-17-69P	34	1307	A	56	1409	4-11-66	66
C1114	1-17-69P	34	1308	A	56	1420	12-16-65	66
1115	2-19-70P	34	1309	A	56	1421	12-16-65	66
C1115	2-19-70P	34	1310	A	56	1422	12-16-65	66
1116	2-19-70P	34	1311	A	56	1423	12-16-65	66
C1116	2-19-70P	34	1312	A	56	1424	12-16-65	66
1117	2-19-70P	34	1313	A	56	1425	12-16-65	66
C1117	2-19-70P	34	1314	A	56	1427	12-16-65	66
1118	7-14-65P	34	1315	A	56	1428	12-16-65	66
C1118	2-26-70P	34	1316	A	56	1440	12-16-65	66
1119	7-14-65P	34	1317	A	56	1441	12-16-65	66
C1119	2-26-70P	34	1318	A	56	1442	12-16-65	66
1120	7-14-65P	34	1319	A	56	1443	12-16-65	66
C1120	2-26-70P	34	1320	A	56	1444	12-16-65	66
1121	12-19-66P	34	1331	A	56	1445	12-16-65	66
C1121	12-19-66P	34	1332	A	56	1475	12-2-71	59
1122	12-19-66P	34	1333	A	56	1476	11-6-69	59
C1122	12-19-66P	34	1334	A	56	1511	4-24-69	72
1123	12-19-66P	34	1335	A	56	1512	7-2-69	72
C1123	12-19-66P	34	1336	A	56	1513	9-17-69	72
1131	1-25-68P	36	1337	A	56	1516	7-9-69	72
1132	1-30-70	36	1338	A	56	1517	7-9-69	72
1134	4-30-70	22	1339	A	56	1518	7-9-69	72
1135	7-27-72	22	1341	A	56	1519	7-9-69	72
1138	7-13-70	28	1342	A	56	1541	3-1-71P	71

SRM	Certificate Date	Page	SRM	Certificate Date	Page	SRM	Certificate Date	Page
1571	10-1-71P	43	2318	A	57	4935-C	7-69P	69
1573	D	43	2331	A	57	4940-B	1-4-68	69
1577	4-15-72P	43	2332	A	57	4941-C	4-69	69
1579	9-29-72	46	2333	A	57	4943	1962	69
1601	6-20-67	45	2334	A	57	4947	2-10-64	69
1602	6-20-67	45	2335	A	57	4949	D	69
1603	6-20-67	45	2336	A	57	4950-B	5-68	70
1604a	1-25-71	45	2338	A	57	4951	11-56	70
1605	3-1-68	45	2339	A	57	4953	5-68	70
1606	3-1-68	45	2340	A	57	4955	3-1-68	70
1607	3-1-68	45	3200	A(5-5-71)	73	4956	3-1-68	70
1608	3-1-68	45	4200-B	A(3-69)	70	4957	A(3-1-68)	70
1609	3-1-68	45	4201-B	A(6-70)	70	4958	A(3-1-68)	70
1610	10-31-69	45	4202	A(3-1-68)	70	4959	A(3-1-68)	70
1611	10-31-69	45	4203	D	70	4060	A(3-1-68)	70
1613	10-31-69	45	4205	A(9-68)	70	4961	A(3-1-68)	70
1621	12-11-67	45	4206	A(9-68)	70	4962	A(3-1-68)	70
1622	12-11-67	45	4207	A(3-69)	70	4963	A(3-1-68)	70
1623	4-7-71	45	4210	A(5-69)	70	4964-B	A(12-31-65)	70
1624	4-7-71	45	4211	A(5-70)	70	4990-B	B	69
1625	A(12-1-70)	45	4212	A(6-71)	70	4991-C	A(5-69)	70
1626	A(8-12-71)	45	4213	A(5-70)	70	4996-B	A(5-69)	70
1627	A(8-12-71)	45	4214	D	70	U-0002	7-30-70P	52
1630	11-3-71P	46	4215	D	70	U-005	4-21-69P	52
1631	D	46	4216	D	70	U-010	4-21-69P	52
1651	11-12-68	61	4222	8-67	69	U-015	4-21-69P	52
1652	11-12-68	61	4223	8-67	69	U-020	4-21-69P	52
1653	11-12-68	61	4224	8-67	69	U-030	4-21-69P	52
1654	4-7-71	61	4226	8-69	69	U-050	4-21-69P	52
1810	5-5-72	76	4228	4-71	69	U-100	6-23-66P	52
2001	A(4-16-71)	67	4229	A(5-72)	69	U-150	8-5-66P	52
2002	A(4-16-71)	67	4230	D	69	U-200	6-1-66P	52
2003	A(4-16-71)	67	4321	D	69	U-350	5-23-66P	52
2005	A(4-16-71)	67	4232	D	69	U-500	5-24-66P	52
2006	A(4-16-71)	67	4233	D	69	U-750	2-11-66P	52
2007	A(4-16-71)	67	4234	D	69	U-800	2-11-66P	52
2008	A(4-16-71)	67	4235	A(10-69P)	69	U-850	2-11-66P	52
2101-5	A	65	4236	D	69	U-900	2-11-66P	52
2106	C	74	4240	D	70	U-930	2-11-66P	52
2141	8-28-70	42	4242-B	A(5-72)	70	U-970	7-9-70P	52
2142	9-1-70	42	4243-B	A(5-72)	70			
2143	D	42	4244-B	A(5-72)	70			
2144	D	42	4245	D	69			
2186-I	5-28-68	54	4246	D	69			
2186-II	5-28-68	54	4247	D	69			
2191	5-28-68	54	4300	5-72	69	RM-1C	E	77
2192	5-28-68	54	4301	5-72	69	RM-1R	E	77
2201	2-22-71	54	4302	D	69			
2202	2-22-71	54	4303	D	69			
2301	A	57	4304	D	69			
2302	A	57	4305	D	69	GM-1	F(6-14-71)	78
2303	A	57	4306	D	69	GM-2	F(6-14-71)	78
2304	A	57	4307	D	69	GM-5	F	78
2305	A	57	4900	H	68	GM-2007	F	78
2306	A	57	4901	H	68			
2307	A	57	4902	H	68			
2308	A	57	4904-D	A(4-70)	68			
2311	A	57	4906	A(5-69)	68			
2312	A	57	4907	D	68			
2313	A	57	4921-C	8-6-64	69			
2314	A	57	4922-E	A(5-67)	69			
2315	A	57	4925	7-58	69			
2316	A	57	4926	9-3-61	69			
2317	A	57	4927	9-3-61	69			
			4929-C	6-70	69			

OTHER NBS PUBLICATIONS IN THIS SERIES

- NBS Spec. Publ. 260, Catalog of Standard Reference Materials, July 1970. 75 cents.* (Supersedes NBS Misc. Publ. 260, January 1968 and NBS Misc. Publ. 241, March 1962.)
- NBS Misc. Publ. 260-1, Standard Reference Materials: Preparation of NBS White Cast Iron Spectrochemical Standards, June 1964. 30 cents.*
- NBS Misc. Publ. 260-2, Standard Reference Materials: Preparation of NBS Copper-Base Spectrochemical Standards, October 1964. 35 cents.*
- NBS Misc. Publ. 260-3, Standard Reference Materials: Metallographic Characterization of an NBS Spectrometric Low-Alloy Steel Standard, October 1964. 20 cents.* (Out of print).
- NBS Misc. Publ. 260-4, Standard Reference Materials: Sources of Information on Standard Reference Materials, February 1965. 20 cents.* (Out of print).
- NBS Misc. Publ. 260-5, Standard Reference Materials: Accuracy of Solution X-Ray Spectrometric Analysis of Copper-Base Alloys, March 1965. 25 cents.* (Out of print).
- NBS Misc. Publ. 260-6, Standard Reference Materials: Methods for the Chemical Analysis of White Cast Iron Standards, July 1965. 45 cents.*
- NBS Misc. Publ. 260-7, Standard Reference Materials: Methods for the Chemical Analysis of NBS Copper-Base Spectrochemical Standards, October 1965. 60 cents.*
- NBS Misc. Publ. 260-8, Standard Reference Materials: Analysis of Uranium Concentrates at the National Bureau of Standards, December 1965. 60 cents.* (Out of print).
- NBS Misc. Publ. 260-9, Standard Reference Materials: Half Lives of Materials Used in the Preparation of Standard Reference Materials of Nineteen Radioactive Nuclides Issued by the National Bureau of Standards, November 1965. 15 cents.*
- NBS Misc. Publ. 260-10, Standard Reference Materials: Homogeneity Characterization on NBS Spectrometric Standards II: Cartridge Brass and Low-Alloy Steel, December 1965. 30 cents.*
- NBS Misc. Publ. 260-11, Standard Reference Materials: Viscosity of a Standard Lead-Silica Glass, November 1966. 25 cents.*
- NBS Misc. Publ. 260-12, Standard Reference Materials: Homogeneity Characterization of NBS Spectrometric Standards III: White Cast Iron and Stainless Steel Powder Compact, September 1966. 20 cents.*
- NBS Misc. Publ. 260-13, Standard Reference Materials: Mossbauer Spectroscopy Standard for the Chemical Shift of Iron Compounds, July 1967. 40 cents.*
- NBS Misc. Publ. 260-14, Standard Reference Materials: Determination of Oxygen in Ferrous Materials — SRM 1090, 1091, and 1092, September 1966. 30 cents.*
- NBS Misc. Publ. 260-15, Standard Reference Materials: Recommended Method of Use of Standard Light-Sensitive Paper for Calibrating Carbon Arcs Used in Testing Textiles for Colorfastness to Light, June 1967. 20 cents.*
- NBS Spec. Publ. 260-16, Standard Reference Materials: Homogeneity Characterization of NBS Spectrometric Standards IV: Preparation and Microprobe Characterization of W-20% Mo Alloy Fabricated by Powder Metallurgical Methods, January 1969. 35 cents.*
- NBS Spec. Publ. 260-17, Standard Reference Materials: Boric Acid; Isotopic and Assay Standard Reference Materials, February 1970. 65 cents.*
- NBS Spec. Publ. 260-18, Standard Reference Materials: Calibration of NBS Secondary Standard Magnetic Tape (Computer Amplitude Reference) Using the Reference Tape Amplitude Measurement "Process A", November 1969. 50 cents.*
- NBS Spec. Publ. 260-19, Standard Reference Materials: Analysis of Interlaboratory Measurements on the Vapor Pressure of Gold (Certification of Standard Reference Material 745), January 1970. 30 cents.*
- NBS Spec. Publ. 260-20, Standard Reference Materials: Preparation and Analysis of Trace Element Glass Standards. (In preparation)
- NBS Spec. Publ. 260-21, Standard Reference Materials: Analysis of Interlaboratory Measurements on the Vapor Pressures of Cadmium and Silver, January 1971. 35 cents.*
- NBS Spec. Publ. 260-22, Standard Reference Materials: Homogeneity Characterization of Fe-3Si Alloy, February 1971. 35 cents.*

- NBS Spec. Publ. 260-23, Standard Reference Materials: Viscosity of a Standard Borosilicate Glass, December 1970. 25 cents.*
- NBS Spec. Publ. 260-24, Standard Reference Materials: Comparison of Redox Standards, January 1972. \$1.*
- NBS Spec. Publ. 260-25, Standard Reference Materials: A Standard Reference Material Containing Nominally Four Percent Austenite, February 1971. 30 cents.*
- NBS Spec. Publ. 260-26, Standard Reference Materials: National Bureau of Standards—U.S. Steel Corporation Joint Program for Determining Oxygen and Nitrogen in Steel, February 1971. 50 cents.*
- NBS Spec. Publ. 260-27, Standard Reference Materials: Uranium Isotopic Standard Reference Materials, April 1971. \$1.25.*
- NBS Spec. Publ. 260-28, Standard Reference Materials: Preparation and Evaluation of SRM's 481 and 482 Gold-Silver and Gold-Copper Alloys for Microanalysis, August 1971. \$1.*
- NBS Spec. Publ. 260-29, Standard Reference Materials: Calibration of NBS Secondary Standard Magnetic Tape (Computer Amplitude Reference) Using the Reference Tape Amplitude Measurement "Process A-Model 2", June 1971. 60 cents.*
- NBS Spec. Publ. 260-30, Standard Reference Materials: Standard Samples Issued in the USSR (A Translation from the Russian), June 1971. \$1.*
- NBS Spec. Publ. 260-31, Standard Reference Materials: Thermal Conductivity of Electrolytic Iron SRM 734 from 4 to 300 K, November 1971. 35 cents.*
- NBS Spec. Publ. 260-32, Standard Reference Materials: The Cooperative Study of Temperature Scale Standards for DTA by ICTA and NBS. (In preparation)
- NBS Spec. Publ. 260-33, Standard Reference Materials: Comparison of Original and Supplemental SRM 705, Narrow Molecular Weight Distribution Polystyrene, H. L. Wagner, May 1972. 35 cents.*
- NBS Spec. Publ. 260-34, Standard Reference Materials: Thermoelectric Voltage, April 1972. 40 cents.*
- NBS Spec. Publ. 260-35, Standard Reference Materials: Thermal Conductivity of Austenitic Stainless Steel, SRM 735 from 5 to 280 K, April 1972. 35 cents.*
- NBS Spec. Publ. 260-36, Standard Reference Materials: A Referee Method for the Determination of Calcium in Serum. SRM 915, May 1972. \$1.25.*
- NBS Spec. Publ. 260-37, Standard Reference Materials: Methods of Analysis of NBS Clay Standards, June 1972. 75 cents.*
- NBS Spec. Publ. 260-38, Standard Reference Materials: Preparation and Calibration of Standards of Spectral Specular Reflectance, May 1972. 60 cents.*
- NBS Spec. Publ. 260-39, Standard Reference Materials: The Eddy Current Decay Method for Resistivity Characterization of High-Purity Metals, May 1972. 55 cents.*
- NBS Spec. Publ. 260-40, Standard Reference Materials: Selection of Thermal Analysis Temperature Standards Through a Cooperative Study (SRM 758, 759, 760), August 1972. 65 cents.*
- NBS Spec. Publ. 260-41, Standard Reference Materials: Use of Standard Light-Sensitive Paper for Calibrating Carbon Arcs used in Testing Textiles for Colorfastness to Light, August 1972. 30 cents.*
- NBS Spec. Publ. 260-42, Standard Reference Materials: The Characterization of Linear Polyethylene, SRM 1475, September 1972. 45 cents.*
- NBS Spec. Publ. 260-43, Standard Reference Materials: Preparation and Homogeneity Characterization of an Austenitic Iron-Chromium-Nickel Alloy, November 1972. 45 cents.*
- NBS Spec. Publ. 260-44, Standard Reference Materials: Preparation and Use of Superconductive Fixed Point Devices, SRM 767, December 1972. 75 cents.*
- NBS Spec. Publ. 260-45, Standard Reference Materials: Powdered Lead-Based Paint, SRM 1579. (In preparation)

*Send order with remittance to: Superintendent of Documents, U.S. Government Printing Office, Washington, D.C. 20402. Remittance from foreign countries should include an additional one-fourth of the purchase price for postage.

U.S. DEPT. OF COMM. BIBLIOGRAPHIC DATA SHEET	1. PUBLICATION OR REPORT NO. NBS SP260, 1973 Catalog	2. Gov't Accession No.	3. Recipient's Accession No.
4. TITLE AND SUBTITLE Standard Reference Materials		5. Publication Date April 1973	
		6. Performing Organization Code	
7. AUTHOR(S) R. W. Seward, Technical Representative		8. Performing Organization	
9. PERFORMING ORGANIZATION NAME AND ADDRESS NATIONAL BUREAU OF STANDARDS DEPARTMENT OF COMMERCE WASHINGTON, D.C. 20234		10. Project/Task/Work Unit No. 3020585	
		11. Contract/Grant No.	
12. Sponsoring Organization Name and Address Same as No. 9		13. Type of Report & Period Covered Interim	
		14. Sponsoring Agency Code	
15. SUPPLEMENTARY NOTES Supersedes NBS Special Publication 260, July 1970 Edition			
16. ABSTRACT (A 200-word or less factual summary of most significant information. If document includes a significant bibliography or literature survey, mention it here.) This Catalog lists and describes the Standard Reference Materials (SRM's), Research Materials (RM's), and General Materials (GM's) currently distributed by the National Bureau of Standards, as well as many of the materials currently in preparation. SRM's are used to calibrate measurement systems and to provide a central basis for uniformity and accuracy of measurement. The unit and quantity, the type, and the certified characterization are listed for each SRM, as well as directions for ordering. The RM's are not certified, but are issued to meet the needs of scientists engaged in materials research. RM's are issued with a "Report of Investigation", the sole authority of which is the author of the report. The GM's are standardized by some agency other than NBS. NBS acts only as a distribution point and does not participate in the standardization of these materials. Announcements of new and renewal SRM's, RM's, and GM's are made in the semi-annual supplements to this Catalog, SRM Price and Availability List, the NBS Technical News Bulletin, and in scientific and trade journals.			
17. KEY WORDS (Alphabetical order, separated by semicolons) Analysis; characterization, composition; General Materials; properties; Research Materials; Standard Reference Materials			
18. AVAILABILITY STATEMENT <input checked="" type="checkbox"/> UNLIMITED. <input type="checkbox"/> FOR OFFICIAL DISTRIBUTION. DO NOT RELEASE TO NTIS.		19. SECURITY CLASS (THIS REPORT) UNCLASSIFIED	21. NO. OF PAGES 96
		20. SECURITY CLASS (THIS PAGE) UNCLASSIFIED	22. Price \$1.25, domestic postpaid \$1.00, GPO Bookstore

NBS TECHNICAL PUBLICATIONS

PERIODICALS

JOURNAL OF RESEARCH reports National Bureau of Standards research and development in physics, mathematics, and chemistry. Comprehensive scientific papers give complete details of the work, including laboratory data, experimental procedures, and theoretical and mathematical analyses. Illustrated with photographs, drawings, and charts. Includes listings of other NBS papers as issued.

Published in two sections, available separately:

• Physics and Chemistry (Section A)

Papers of interest primarily to scientists working in these fields. This section covers a broad range of physical and chemical research, with major emphasis on standards of physical measurement, fundamental constants, and properties of matter. Issued six times a year. Annual subscription: Domestic, \$17.00; Foreign, \$21.25.

• Mathematical Sciences (Section B)

Studies and compilations designed mainly for the mathematician and theoretical physicist. Topics in mathematical statistics, theory of experiment design, numerical analysis, theoretical physics and chemistry, logical design and programming of computers and computer systems. Short numerical tables. Issued quarterly. Annual subscription: Domestic, \$9.00; Foreign, \$11.25.

TECHNICAL NEWS BULLETIN

The best single source of information concerning the Bureau's measurement, research, developmental, cooperative, and publication activities, this monthly publication is designed for the industry-oriented individual whose daily work involves intimate contact with science and technology—for engineers, chemists, physicists, research managers, product-development managers, and company executives. Includes listing of all NBS papers as issued. Annual subscription: Domestic, \$6.50; Foreign, \$8.25.

NONPERIODICALS

Applied Mathematics Series. Mathematical tables, manuals, and studies.

Building Science Series. Research results, test methods, and performance criteria of building materials, components, systems, and structures.

Handbooks. Recommended codes of engineering and industrial practice (including safety codes) developed in cooperation with interested industries, professional organizations, and regulatory bodies.

Special Publications. Proceedings of NBS conferences, bibliographies, annual reports, wall charts, pamphlets, etc.

Monographs. Major contributions to the technical literature on various subjects related to the Bureau's scientific and technical activities.

National Standard Reference Data Series. NSRDS provides quantitative data on the physical and chemical properties of materials, compiled from the world's literature and critically evaluated.

Product Standards. Provide requirements for sizes, types, quality, and methods for testing various industrial products. These standards are developed cooperatively with interested Government and industry groups and provide the basis for common understanding of product characteristics for both buyers and sellers. Their use is voluntary.

Technical Notes. This series consists of communications and reports (covering both other-agency and NBS-sponsored work) of limited or transitory interest.

Federal Information Processing Standards Publications. This series is the official publication within the Federal Government for information on standards adopted and promulgated under the Public Law 89-306, and Bureau of the Budget Circular A-86 entitled, Standardization of Data Elements and Codes in Data Systems.

Consumer Information Series. Practical information, based on NBS research and experience, covering areas of interest to the consumer. Easily understandable language and illustrations provide useful background knowledge for shopping in today's technological marketplace.

BIBLIOGRAPHIC SUBSCRIPTION SERVICES

The following current-awareness and literature-survey bibliographies are issued periodically by the Bureau:

Cryogenic Data Center Current Awareness Service (Publications and Reports of Interest in Cryogenics).

A literature survey issued weekly. Annual subscription: Domestic, \$20.00; foreign, \$25.00.

Liquefied Natural Gas. A literature survey issued quarterly. Annual subscription: \$20.00.

Superconducting Devices and Materials. A literature survey issued quarterly. Annual subscription: \$20.00.

Send subscription orders and remittances for the preceding bibliographic services to the U.S. Department of Commerce, National Technical Information Service, Springfield, Va. 22151.

Electromagnetic Metrology Current Awareness Service (Abstracts of Selected Articles on Measurement Techniques and Standards of Electromagnetic Quantities from D-C to Millimeter-Wave Frequencies). Issued monthly. Annual subscription: \$100.00 (Special rates for multi-subscriptions). Send subscription order and remittance to the Electromagnetic Metrology Information Center, Electromagnetics Division, National Bureau of Standards, Boulder, Colo. 80302.

Order NBS publications (except Bibliographic Subscription Services) from: Superintendent of Documents, Government Printing Office, Washington, D.C. 20402.

U.S. DEPARTMENT OF COMMERCE
National Bureau of Standards
Washington, D.C. 20234

OFFICIAL BUSINESS

Penalty for Private Use, \$300

POSTAGE AND FEES PAID
U.S. DEPARTMENT OF COMMERCE
COM 215

